



**<23 km**  
Range per Charge



**147 kg**  
Weight Capacity



**122.5 cm**  
Turning Radius



**GO GO®**  
**ENDURANCE Li**

## For People Who Are Always On the GO

The Go Go® Endurance Li is a highly convenient travel companion that provides an effortless mode of transportation for people who are always on the go. This 4-wheel scooter can easily be disassembled into 5 portable pieces, making it easy to transport or store in the trunk of a car. With the lightweight lithium-ion battery pack, you can travel up to 23 km on a single charge, giving you plenty of freedom to explore your surroundings and reach more destinations.

Visit [pridemobility.com.au](http://pridemobility.com.au) for more information







Long-Lasting, Lithium-ion Battery



Feather-Touch Disassembly



CTS Suspension on Front and Rear

|  |  |
|--|--|
| RANGE PER CHARGE (MAX) <sup>1</sup>      | Up to 23 km  |
| MAXIMUM WEIGHT CAPACITY                  | 147 kg   |
| TURNING RADIUS                           | 122.5 cm   |
| MAXIMUM SPEED <sup>1</sup>               | 6.7 km/h   |
| LENGTH & WIDTH <sup>2</sup>              | L: 107 cm W: 52 cm   |
| SEAT TYPE<br>WIDTH & DEPTH<br>(STANDARD) | Compact, foldable seat<br>W: 43 cm (17") D: 43 cm (17")<br>W: 50.8 cm (20") D: 45.7 cm (18") |
| SEAT-TO-DECK HEIGHT                      | 41.2 cm - 45 cm  |
| SEAT-TO-GROUND HEIGHT                    | 53.9 cm - 59 cm  |
| GROUND CLEARANCE <sup>5</sup>            | 3.6 cm at motor  |

|  |                                    |
|--|------------------------------------|
| OVERALL WEIGHT (NOT INC. BATTERY) <sup>2</sup> | 42.4 kg                            |
| WEIGHT OF HEAVIEST PIECE                       | 18.45 kg (front section)           |
| SUSPENSION                                     | front & rear CTS                   |
| 3W/4W OPTION                                   | 4W                                 |
| FRONT TYRES                                    | 8" solid                           |
| REAR TYRES                                     | 8" solid                           |
| BATTERY <sup>4</sup> REQUIREMENTS              | 16Ah Lithium-ion (2 x 8Ah) (385Wh) |
| BATTERY CHARGER                                | Off-board, 3.5A Lithium-ion        |
| BATTERY WEIGHT                                 | 5.4 kg each                        |
| COLOUR OPTIONS <sup>6</sup>                    | ■ Garnet Red ■ Euro Grey           |

1. Varies with user weight, terrain type, battery amp hour (Ah), battery charge, battery condition and tyre condition.
2. Includes standard seat weight.
3. Width of base (not including seat).
4. Standard 16Ah battery does not meet airline requirements. To make this scooter airline compliant please remove one of the 8Ah batteries before travelling on an aircraft.
5. Does not account for anti-tip clearance.
6. All shroud panels included

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of 3%.



Transport of this item on commercial aircraft is subject to federal regulations and airline policy. Please consult your airline prior to your flight to ensure your Go Go<sup>®</sup> Endurance Li and its 1x 8Ah lithium-ion battery are permitted.