Model	BA140
Weight Capacity	400 lb (181 kg)
Maximum Speed ²	Up to 8 mph (12.8 km/h)
Ground Clearance ³	1.57" at motor skid plate 2.68" at front bumper 4.45" at mid-deck
Turning Radius ³	53"
Dimensions ³	Length: 49.6" Width: 27.83"
Seat-to-Deck Height ³	18.11"
Front & Rear Tires ³	9" solid
Suspension	Full suspension
Range Per Charge ^{2,5}	20Ah Battery Pack: Up to 14.5 mi (23.3 km) at 220 lb (100 kg) Up to 11.25 mi (18.1 km) at 400 lb (181 kg) 40Ah Battery Pack: Up to 29 mi (46.7 km) at 220 lb (100 kg) Up to 22.5 mi (36.2 km) at 400 lb (181 kg)
Total Weight ⁶	183.4 lb (83.2 kg) with 20 Ah battery 197.8 lb (89.7 kg) with 40 Ah battery
Standard Seating ³	Type: Stadium seat with slider Material: Black Vinyl Width: 18.31" Depth: 17.32" Weight: 34.8 lb (15.8 kg)
Drive System	Rear-wheel drive, sealed transaxle, 36-volt DC motor
Dual Braking System	Regenerative and electromechanical
Battery Require- ments ⁴	Type: Lithium-ion battery pack Size: 20Ah or 40Ah
Battery Pack Weight ⁷	20 Ah: 20.8 lb (9.4 kg) 40 Ah: 35.2 lb (16 kg)
Battery Charger	4 amp, off-board charger

- 1 Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
- 2 Varies with user weight, terrain type, battery amp-hour rating (Ah), battery charge, battery condition, tire type and tire condition. This specification can be subject to a variance of +10%, -5%.
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 4 Lithium-ion, AGM, Gel-Cell, or Hybrid Gel (Hybrid Gel is defined as Nano Silica Gel Electrolyte with High Quality Glass Mat Separators) type required.
- 5 Tested in accordance with ANSI/RESNA, WC Vol 2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculation based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 6 Includes standard seat.
- 7 Battery weight may vary based on manufacturer.








TOP SPEED

SUSPENSION

CAPACITY

CA Bandit Trifold 052224




BUILT FOR THE BOLD.

It's time to discover what's really out there. With impressive ground clearance, a top speed of 8 mph (12.8 km/h), and a lithium-ion battery pack which powers up to a 29-mi (46.7 km) range, the Baja® Bandit is built for exploring the unknown and confidently tackling rugged terrain away from the beaten path.

FEATURES

- 4.5" ground clearance
- 400 lb (181 kg) weight capacity
- Competitive 53" turning radius
- Infinite tiller adjustment lever
- USB port for charging your devices
- Full LED lighting package for maximum visibility
- Full bumper package (standard)
- Supplemental hand brake for additional stopping power
- Feather-touch disassembly
- Seat post with full suspension
- The upgradable dual lithium-ion battery configuration will deliver up to a 29-mi (46.7 km) range on a single charge (14.5 mi (23.3 km) with the standard single lithium-ion battery)
- 3-year warranty on lithium-ion battery



ACCESSORIES

- 50/60/70" lap belt
- Cup holder
- Cell phone holder
- Rearview mirror
- Saddlebag
- Safety flag
- Weather cover

Only one rear accessory can be used at a time.

- Walker holder
- 3" single or double cane or crutch holder
- Forearm crutch holder
- Wishbone crutch holder
- Oxygen holder
- Quad cane holder
- Rear basket

Cup holder, rear basket and saddlebag can be ordered as a bundle.