

Baja™ Wrangler® 2

Model Number	BA400
Weight Capacity	450 lb / 204.1 kg
Maximum Speed ¹	11.5 mph/18.5 km/h at 200 lb/90.7 kg 10.6 mph/16.9 km/h at 450 lb/204 kg
Ground Clearance ³	8.25" at mid-frame 5.5" at motor
Turning Radius ³	94.75"
Overall Length ³	61.25"
Overall Width ³	32.5"
Seat-to-Ground Height ³	29-30.5"
Seat-to-Deck Height ³	20.5-22"
Front Tires ³	14.5" pneumatic
Rear Tires ³	14.5" pneumatic
Suspension	Front and rear (CTS Suspension)
Range Per Charge ^{1,2}	75Ah: 16.5 mi/26.5 km at 200 lb/90.7 kg 13.3 mi/21.3 km at 450 lb/204.1 kg 100 Ah: 21.4 mi/34.5 km at 200 lb/90.7 kg 17.3 mi/27.8 km at 450 lb/204.1 kg
Total weight ⁵	274 lb / 124.3 kg
Standard Seating ³	Black vinyl high-back with limited recline memory foam seat with sliders Weight: 48.6 lb/22.1 kg Size: 18"x18-20"
Drive System	24-volt DC dual in-line motors, rear-wheel drive
Braking System	Regenerative and electromechanical
Battery Requirements ^{4,6}	(2) 75 AH (included) or 100 AH ⁷ Weight: 48 lb/21.8 kg (75 AH) 60.5 lb/27.4 kg (100 AH)
Battery Charger	8-amp, off-board
Warranty	Frame: lifetime limited Drivetrain: 2-year limited Electronics: 2-year limited Battery: 13-months limited

¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.

² Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

³ Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.

⁴ Battery weight may vary based on battery manufacturer.

⁵ Includes standard seat and batteries.

⁶ AGM or gel cell type required.

⁷ Optional.



Beamsville, Ontario

pridemobility.ca |    

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

CA_BajaWrangler2_010224



bajamobility.ca |    

THE SPIRIT OF ADVENTURE

18" x 18-20" black vinyl, limited recline memory foam seat with manual adjustment, heavy-duty seat post, sliders and height, width, depth, and angle adjustable armrests and 60" lap belt

Optional seat sizes:
16" x 16-18"
20" x 18-20"
20" x 20-22"



Rugged and built to last, the Baja™ Wrangler® 2 is everything you want in an outdoor scooter. Equipped with front and rear suspension, this aggressive scooter is designed for cruising on walking trails or driving around your property. Make your adventures legendary with the Baja™ Wrangler® 2!

Features

- CTS Suspension includes industry-leading front and rear suspension for a smoother ride
- Steering dampening system absorbs uncontrolled movement over rugged terrain, limits fatigue and reduces energy needed to control handles
- Black vinyl, limited recline memory foam seat with manual adjustment, heavy-duty seat post, sliders and height, width, depth, and angle adjustable armrests and 60" lap belt
- User friendly LED console with variable speed control displays time, temperature, miles driven, and trips taken plus horn and light buttons
- Adjustable delta tiller with ergonomic, wraparound handles and XLR charging port
- USB charger for charging portable devices
- Easily accessible tie-down points (for transport of unoccupied scooter)
- Full LED lighting package includes headlights, hazard lights, and directional signals
- Rear backup sensor LED lights indicate when you are too close to an object when in reverse
- Hydraulic hand brake with standard regenerative and electro-mechanical braking system
- Easy battery installation

Colors



Desert Sand

True Blue

Optional Seating

- 16" x 16-18", 20" x 18-20" and 22" x 20-22" black vinyl, limited recline, memory foam seat, with sliders

Accessories

- 3" single or double cane/crutch holder
- 50/70" Lap belt
- Cup holder
- Cell phone holder
- Forearm crutch holder
- Front basket
- Oxygen holder
- Quad cane holder
- Rear basket
- Additional rearview mirror
- Saddlebag
- Walker holder
- Wishbone crutch holder



Backup sensor with LED lights



User friendly console



Front headlights