

|   |   |
|---|---|
| Model   | S710DXW   |
| Weight Capacity                               | 350 lb / 159 kg   |
| Maximum Speed <sup>1</sup>                    | Up to 7.43 mph / 11.96 km/h   |
| Ground Clearance <sup>3</sup>                 | 2.13" at motor  |
| Turning Radius <sup>3</sup>                   | 51.7"   |
| Overall Length <sup>3</sup>                   | 47.05"  |
| Overall Width <sup>3</sup>                    | 25.94"  |
| Seat-to-Ground Height <sup>3</sup>            | 23.175"-25.175"   |
| Seat-to-Deck Height <sup>3</sup>              | 17.25"-19.25"   |
| Front Tires <sup>3</sup>                      | 10" pneumatic   |
| Rear Tires <sup>3</sup>                       | 10" pneumatic   |
| Suspension                                    | Front and rear (CTS Suspension)   |
| Range Per Charge <sup>1,2</sup>               | Up to 19.8 mi / 31.8 km @ 220 lb / 100 kg<br>Up to 16.2 mi / 26 km @ 350 lb / 159 kg                                  |
| Total Weight <sup>3,4,5</sup>                 | Without batteries: 165.3 lb / 75 kg<br>With batteries: 227.9 lb / 103.4 kg  |
| Heaviest Piece When Disassembled <sup>3</sup> | 68.6 lb / 31.1 kg (front section)   |
| Standard Seating <sup>3</sup>                 | Type: High-Back Reclining<br>Weight: 43.8 lb / 19.9 kg<br>Material: black vinyl<br>Dimensions: Width: 18", Depth: 18" |
| Drive System                                  | Rear-wheel drive, sealed transaxle,<br>24-volt DC motor, 120 Amp controller   |
| Dual Braking System                           | Regenerative and<br>electromechanical   |
| Battery Requirements <sup>4,6</sup>           | Size: (2) 50 Ah included<br>Weight: 31.3 lb / 14.2 kg each  |
| Battery Charger                               | Off-board, 5 amp  |
| Warranty                                      | Frame: lifetime limited<br>Drivetrain: 2-year limited<br>Electronics: 2-year limited<br>Battery: 13-months            |

1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.

2 Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.

4 Battery weight may vary based on battery manufacturer.

5 Includes standard seat weight.

6 AGM or gel cell type required.

**PRIDE  
MOBILITY®**

pridemobility.ca



NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.



**AFFORDABLE  
COMFORT**



**LIVE YOUR BEST®**

pridemobility.ca



# BUILT FOR STYLE, COMFORT, AND RANGE



Rear LED Lights



Front and Rear CTS Suspension



LED Accent Lighting on Upper Tiller

Built for style and comfort, the Colt Sport can get you where you need to be smoothly and quickly, whether near or far. The ample storage options make it easy to take along whatever you need for your outing with room to spare for your purchases.

## Standard Features:

- Delta Tiller with push-pull throttle control
- Large front basket
- Dual rear-view mirrors
- Full LED light package with turn-signals
- Automatic brake lighting
- Tiller with puddle light, built-in storage, and cup-holder
- Comfort Track Suspension (CTS)
- Standard pneumatic tires
- 120 Amp controller
- High-capacity batteries
- High-back Captain's seat

## Available Colours\*\*



Stone  
Grey

Black  
Pearl

## Accessories:

- Cell phone holder
- Cup holder
- Saddle bag
- Rear basket
- Safety flag
- 50", 60", 70" lap belt
- Stadium seat\*\*\*
- Oxygen tank holder
- Cane/crutch holder (single, double)
- Cane holders (Quad, Single, Double, and Forearm)
- Wishbone crutch holder
- Walker holder
- Weather cover

\*Only one rear accessory can be used at a time.

\*\*Matte finish.

\*\*\* Seat slide may be limited by rear accessories.