

COMFORT AND STYLE

Pride[®] Elegance Collection Lift Chairs are designed to be the ultimate in lift chair comfort, style and performance. Available in several beautiful fabrics and styles to fit any décor while utilizing a state-ofthe-art lift system.

> Easy-to-use hand control with USB charging port

ERGY

BATTERY BACKUP



Scan here for more information on the **Elegance Collection**







WORRY-FREE OPERATION

FEATURES & BENEFITS

- Industry-first Total Comfort Seating features a unique spring design and one-piece soft foam seat for outstanding relaxation (not available on LC-108)
- Premium furniture grade laminate/hardwood frames for added strength and durability
- Easy-to-use hand control with USB charger port (USB charging port not available when selecting the
 optional heat and massage)
- Standard head and arm covers
- Large, dual pockets for extra storage
- Back has access for adjusting the fill for comfort and fit (Back cushion of the LC-450 and LC-108 is fully removable and adjustable)
- Quiet and smooth lift system









	2	-		
MODEL	LC-108	LC-570S	LC-570M	LC-570L
POSITION	Infinite-position Trendelenburg Chaise Lounger 5'4" - 6'0"	3-Position Chaise Lounger 5′3″ and below	3-Position Chaise Lounger 5'4" - 5'9"	3-Position Chaise Lounger 5'10″ an above
SIZE	Medium	Small	Medium	Large
WEIGHT CAPACITY	375 lbs.	375 lbs.	375 lbs.	375 lbs.
SEAT-TO-FLOOR	19″	18.5″	19.5″	20.5″
SEAT DEPTH	20″	19.5″	19.5″	21.5″
SEAT WIDTH	19.5″	17″	18.5″	20.5″
TOP OF BACK TO SEAT	30″	23.5″	26″	27″
FULLY PADDED CHAISE	Yes (Removable and adjustable)	Yes	Yes	Yes
SEATING TYPE	Sinuous Spring	Total Comfort	Total Comfort	Total Comfort
DISTANCE FROM WALL	26.5″	16.5″	19″	19.5″
OVERALL WIDTH	30.5″	33.5″	34.5″	35.5″
CHAIR WEIGHT	136 lbs.	110 lbs.	113 lbs.	114 lbs.
WARRANTY	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated
PACIFIC STONE WALNUT BLACK CHERRY CHESTNUT				

CLOUD 9 FABRICS Standard

Beamsville, ON / 888-570-1113 / pridemobility.ca

Note: Measurements may vary +/- 1". The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.



BEST

в

YOUR

LIVE

CA_Elegance_SS_030724

OPTIONS*

- Cleaning kit
- Head pillow
- Footrest extension

for some options.

* Additional charges apply for all

options. Extra lead time may apply