

Part Number	SC443E
Weight Capacity	300 lbs.
Maximum Speed ²	Up to 6.8 km/h (4.23mph)
Ground Clearance ⁴	1.29" at motor, 3.5" under deck
Turning Radius ⁴	39.2"
Overall Length ⁴	41.5"
Overall Width ^{4,6}	19.45"
Seat-to-Deck Height	16.38 – 17.83"
Seat-to-Ground Height	20.75" - 22.25"
Tires	Front: 2"x 7" Rear: 2.5"x 8"
Range ^{2,3}	18 Ah: Up to 18.1 km at 300 lbs. Up to 22.6 km at 220 lbs.
Product Weight	83.6 lbs. without batteries 103 lbs. with 12 Ah batteries 111 lbs. with 18 Ah batteries
Heaviest Piece ⁴	36.2 lbs. (front section)
Standard Seating	Type: Foldable, mesh back Material: Black vinyl Dimensions: 17"x17" Weight: 23 lbs.
Motor Output	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Brakes	Regenerative & electromechanical
Batteries	(Two) 12-volt, sealed lead-acid, non-spillable deep-cycle
Battery Weight ⁵	12 Ah: 20 lbs. (standard) 18 Ah: 27.4 lbs.
Battery Charger	2 amp off-board charger
Warranty	Lifetime limited on frame; 2-year limited on Electronics; 2-year limited on Drivetrain Batteries: 6 months

- 1 Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
- 2 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- 3 Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 4 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- 5 Battery weight may vary based on manufacturer.
- 6 Dimensions listed are for scooter only and do not include seats or arms. Width may vary depending on seat options and accessories.

Featuring a new
deck design with
more leg room!



Shown in
Standard
Cardinal Red

**PRIDE[®]
MOBILITY**

pridemobility.ca | 888-570-1113

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.



pridemobility.ca

GO GO[®]

ELITE TRAVELLER GEN 3

Featuring a new deck design with **more leg room!**



FEATURING EXCLUSIVE PATENTED



Convenient charging port in the tiller.



Intuitive control panel framed by a soft-touch delta tiller handle.



Tiller angle adjustment knob.

Go Go Elite Traveller[®] Gen 3 features EZ Turn technology and an ergonomically designed deck for improved leg room and pass-through space. A breathable mesh-back seat, wraparound delta tiller, large front basket and feather-touch disassembly also enhance comfort and convenience.

FEATURES

- NEW sloped deck design for improved comfort and leg room
- Front and rear pathway lights
- EZ Turn technology provides unrivaled tight turning without sacrificing the stability of a 4-wheel scooter
- Feather-touch disassembly allows for easy separation into 5 lightweight pieces for transport or storage
- Large front basket included
- Interchangeable shroud panels can be replaced with accessory colors shown below

INCLUDED COLOR PANELS (MATTE)



Cardinal Red

OPTIONAL ACCESSORY COLOR PANELS (MATTE)



Black



Iceberg Blue



Robin's Egg Blue



Sugar Plum



Chell-o

ACCESSORIES*

- XLR
- USB Charger
- Cell phone holder
- Cup holder
- Lap belt
- Oxygen holder
- Rear basket
- Rear view mirror
- Safety flag
- Saddle bag
- Single/double rear mount 3" cane/crutch holder
- Under seat storage
- Walker holder

* Only one rear accessory can be used at a time.