

The Essential Collection from Pride® Power Lift Recliners provides exceptional value and outstanding comfort. Several fabric options will match any décor while delivering rest and worry-free operation.



LC -104


ENERGY  
EFFICIENT  
WITH BUILT-IN  
BATTERY BACKUP  
(BATTERIES REQUIRED)



LC -104



Easy Touch  
Controllers



Scan here for more  
information on the  
**Essential Collection®**



FOLLOW US ON SOCIAL MEDIA

**PRIDE®**  
MOBILITY  
LIVE YOUR BEST®

## FEATURES & BENEFITS

- Easy-to-use hand control adjusts your chair to find your perfect position
- Large dual pockets offer ample storage space that is easily accessible
- Smooth and durable lift and scissor mechanism for long-lasting and quiet performance
- Back has access for adjusting the fill for comfort and fit.
- Integrated battery backup ensures chair is operational when the power goes out (Batteries not included)
- Energy-efficient power box equipped with self-testing electronics and 9' cord
- Lifetime warranty on steel lift frame
- Footrest extension standard (LC-107 only)

## OPTIONS\*

- Cleaning kit
- Non-skid leg levelers

\*Additional charges apply for all options. Extra lead time may apply for some options.



MODEL	LC-102	LC-104	LC-106	LC-107
POSITION	3-Position Chaise Lounger 5'2" and below	3-Position Chaise Lounger 5'4" - 5'9"	3-Position Chaise Lounger 5'4" - 5'9"	3-Position Chaise Lounger 5'4" - 5'9"
SIZE	Super Petite	Medium	Medium	Medium
WEIGHT CAPACITY	250 lbs.	325 lbs.	325 lbs.	325 lbs.
SEAT-TO-FLOOR	16.5"	17.5"	18"	19"
SEAT WIDTH	17"	19.5"	20"	20.5"
SEAT DEPTH	16.5"	18"	19"	19"
TOP OF BACK TO SEAT	24.5"	23"	25"	25"
FULLY PADDED CHAISE	Yes	Yes	Yes	Yes
SEATING TYPE	Sinuous Spring	Sinuous Spring	Sinuous Spring	Sinuous Spring
DISTANCE FROM WALL	16"	16.5"	18"	22"
OVERALL WIDTH	26"	29.5"	34"	34.5"
CONTROLLER	Easy Touch Controller (No USB)	Easy Touch Controller with USB	Easy Touch Controller with USB	Easy Touch Controller with USB port and Infinite
CHAIR WEIGHT	106 lbs.	106 lbs.	106 lbs.	116 lbs.
WARRANTY	1-Year Limited	1-Year Limited	1-Year Limited	1-Year Limited

JAVA

MIDNIGHT

PEWTER

## DARBY FABRICS

Beamsville, ON / 888-570-1113 / pridemobility.ca

Note: Measurements may vary +/- 1". The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.



**PRIDE**<sup>®</sup>  
MOBILITY  
LIVE YOUR BEST<sup>®</sup>