

Model	Go Chair [®] MED
Weight Capacity	300 lb / 136 kg
Maximum Speed ²	Up to 3.8 mph / 6 km/h
Ground Clearance ⁴	1.75" at motor
Turning Radius ⁴	25.5"
Overall Length ⁴	33" with footplate 28.5" without front rigging
Overall Width ⁴	23"
Seat-to-Ground Height	22"-23"
Front Wheels	6" solid casters
Drive Wheels	9" solid
Rear Wheels	2" solid anti-tips
Suspension	N/A
Range Per Charge ^{2,3}	Up to 12.1 mi / 19.5 km at 200 lb / 91 kg
Component Weights ^{4,5,6}	Front Section: 30.75 lb / 13.9 kg Rear Section: 36.5 lb / 16.5 kg Seat: 37 lb / 16.8 kg or 44 lb / 20 kg Batteries: 31.4 lb / 14.2 kg (battery box)
Standard Seating	Type: Memory Foam, Medium Back Material: Black vinyl Dimensions: 18"W x 18"D, or 20"W x 18"D
Drivetrain	Two-motor, rear-wheel drive
Braking System	Regenerative and electromechanical
Standard Electronics	40 amp, PG GC3
Battery Requirements	20 Ah battery pack
Battery Charger	Off-board, 2 amp
Warranty	Frame: 5-year limited Drivetrain: 24-month limited Electronics: 24-month limited Batteries: 13-month limited

- 1 Pride FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
- 2 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
- 3 Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 4 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 5 Battery weight may vary based on battery manufacturer.
- 6 Seat weight includes the joystick.

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

(Canada) 888-570-1113

pridemobility.ca |      

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

CA_GOCHAIRMED-Trifold_110424

GOCHAIR[®]
MED

TRAVEL
MADE SIMPLE



FDA Class II
Medical Device¹

LIVE YOUR BEST[®]

pridemobility.ca |      

POWERFUL, PORTABLE INDEPENDENCE.



FDA Class II Medical Device¹



Disassembles quickly for transport or storage



Optional swing-away joystick swings out at an angle for convenient access



Available with convenient cup holder

The Go Chair® MED is #1 for portability! With convenient feather-touch disassembly, the Go Chair MED allows you to experience lightweight travel and independence on the go. Compact, portable and perfect for travel, the Go Chair® MED is a feature-rich power chair that disassembles easily in just a few steps!

FEATURES

- Memory Foam, Medium Back 18x18 or 20x18 5-point swivel seat, with 60" lap belt
- Compact turning radius
- Large perforated footplate
- One-hand feather-touch disassembly

OPTIONS

- Swing-away joystick for convenient access

AVAILABLE COLOURS



ACCESSORIES*

- Cup holder
- Oxygen tank holder
- Lap belts in 50" or 70"
- Single- or Double-3" cane / crutch holder
- Forearm crutch holder
- Saddle bag
- Safety flag
- XLR USB charger
- Weather cover

* **Only one rear accessory can be used at a time.**