



Model	Go Go Elite Traveller® 2 (4-wheel)
Part Number	SC442E
Weight Capacity	300 lb (136 kg)
Maximum Speed ²	Up to 6.4 km/h (4 mph)
Ground Clearance ⁴	1.38" at motor
Turning Radius ⁴	37"
Overall Length ⁴	39.4"
Overall Width ⁴	19.5"
Seat-to-Deck Height ⁴	16.4" - 17.9"
Tires ⁴	Front: 7" Rear: 8"
Range Per Charge ^{2,3}	22.1 km (13.8 mi) at 220 lb (99.79 kg)
Total Weight Without Batteries ⁴	83 lb (37.7 kg)
Heaviest Piece When Disassembled ⁴	35 lb (15.9 kg) (front section)
Standard Seating ⁴	Type: Folded Weight: 22.5 lb (10.2 kg) Material: Black vinyl Dimensions: 17"x17"
Motor Output	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Brakes	Regenerative and electromechanical
Battery Requirements	(Two) 12-volt, sealed lead-acid, non-spillable deep-cycle
Battery Weight ⁵	18 Ah: 29.6 lb (13.4 kg)
Battery Charger	2 amp off-board charger
Warranty	Lifetime limited on frame 2-year limited on electronics 2-year limited on drivetrain 13-month on batteries

- 1 Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
- 2 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- 3 Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 4 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- 5 Battery weight may vary based on manufacturer.



PRIDE
MOBILITY
LIVE YOUR BEST®

pridemobility.ca |    

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

CA_GoGoET2_Trifold_062724

GO GO
ELITE TRAVELLER®²

4-wheel model

KEEP IT TIGHT

37"
TURNING RADIUS

FDA Class II Medical Device¹

LIVE YOUR BEST®

pridemobility.ca |    

KEEP IT TIGHT



37"
TURNING
RADIUS

New compact, foldable seat with breathable mesh back and armrests

Feather-touch disassembly



Front basket (standard)

Removable footboard basket

Easy-to-grip tiller adjustment knob

Non-scutting flat-free tires

8 color choices
(Black comes standard and choice of one additional colour at no-charge)

FDA Class II
Medical Device¹

KEEP IT TIGHT

Thanks to Pride's patented, EZ Turn Technology, the Go Go Elite Traveller® 2 4-Wheel boasts the tightest turning radius of any 4-wheel scooter on the market. This feature properly equips you to ride with confidence, no matter the destination. You'll enjoy convenient standard features such as feather-touch disassembly, a built-in USB charger and multiple storage baskets. Add it all up, and this scooter is ideal for all your travels, all day long.

FEATURES

- 300 lb (136 kg) weight capacity
- Up to 6.4 km/h (4 mph) max speed
- Improved turning radius of 37"
- Non-scutting, flat-free tires
- Delta tiller with adjustment knob and charging port
- Feather-touch disassembly allows the unit to be easily separated into 5 lightweight pieces
- Front frame-mounted seat post for optimal stability
- Dual voltage charger can charge battery pack on-board or off-board
- Removable footboard basket for additional storage
- Front basket (standard)

MATTE SHROUD COLOURS



ACCESSORIES*

- 12 Ah battery box
- 2.1A XLR USB charger
- Cup holder
- Saddle bag
- Safety flag
- Under seat storage
- Weather cover
- Rear basket
- Rear mount 3" double cane/crutch holder
- Single rear mount 3" cane/crutch holder
- Forearm crutch holder
- Oxygen holder
- Walker holder
- Wishbone crutch holder

* Only one rear accessory can be used at a time.



New breathable mesh back



Improved turning radius of 37" for tight indoor maneuvering



Feather-touch disassembly for quick transport or storage