| Model | Go-Go® Ultra X |
|--|---|
| Part Number | 3-wheel: S39 4-wheel: S49 |
| Weight Capacity | 300 lbs. |
| Maximum Speed ² | 3-wheel: Up to 4.3 mph 4-wheel: Up to 4.1 mph |
| Ground Clearance 4 | 1.25" |
| Turning Radius ⁴ | 3-wheel: 30.75" 4-wheel: 40.25" |
| Overall Length ⁴ | 40.25" |
| Overall Width ⁴ | 19.25" |
| Seat-to-Deck Height | 16.5 - 18" |
| Seat-to-Ground Height | 21.5-23" |
| Tires | Front: 8" Rear: 8" |
| Range Per Charge ^{2,3} | 3-wheel: Up to 14.4 miles (200 lbs.) 3-wheel: Up to 10.7 miles (300 lbs.) 4-wheel: Up to 11.5 miles (200 lbs.) 4-wheel: Up to 8.9 miles (300 lbs.) |
| Total Weight Without Batteries | 3-wheel: 76.25 lbs. 4-wheel: 81.5 lbs. |
| Heaviest Piece When Disassembled ⁴ | 3-wheel: 29.5 lbs. (front section) 4-wheel: 34.5 lbs. (front section) |
| Standard Seating | Type: Compact, foldable Weight: 23 lbs. Material: Black vinyl Dimensions: 17"x17" |
| Motor Output | 24-volt DC motor, sealed transaxle |
| Brakes | Regenerative and electromechanical |
| Battery Requirements ¹ | (Two) 12 AH Weight: 9 lbs. each |
| Battery Charger | Off-board, 2 amp |
| Warranty | Lifetime limited on frame 2-year limited on electronics 2-year limited on drivetrain |

- 1 Battery weight may vary based on manufacturer.
- 2 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- 3 Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 4 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or 3%).



NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or -3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.





The Go-Go® Ultra X is loaded with features such as one-hand disassembly and a convenient dropin battery box for worry-free travel. With a 300 lb. weight capacity and impressive speed and range, the Go-Go Ultra X is an exceptional value.

Additional Features

- Easy-to-grip tiller adjustment knob
- Standard front basket
- · Feather-touch disassembly easily disassembles into 5 lightweight pieces for transport and storage
- Charger port conveniently located on tiller to charge your scooter
- Dual voltage charger permits charging the battery pack on-board or off-board for added convenience
- Pride's exclusive black, non-scuffing tires

Shroud Color

Black

Accessories*

- Weather cover
- 12 Ah battery box
- Saddle bag
- Cup holder
- Safety flag
- Under seat storage
- 2.1A XLR USB charger
- Cell phone holder Rear basket

- Walker holder
- Oxygen holder
- Wishbone crutch holder
- Single rear mount 3" cane/crutch holder
- Rear mount 3" double cane/crutch holder
- Forearm crutch holder



* Only one rear accessory can be used at a time.