

Specifications

Weight Capacity	300 lbs.
Maximum Speed ¹	Up to 4.0 mph (elevated and non-elevated position)
Ground Clearance ²	1" at motor bracket
High-Back Seat to Floor Range	21.25" – 33"
Turning Radius ²	20.75"
Overall Size ^{2, 6}	Length: 41.93" Width: 23.19"
Drive Controller	PG VR2
Drive Wheels	10" flat free
Caster Wheels	6" flat free
Suspension	Active-Trac® with Mid-Wheel 6® Technology
Range ^{1, 4}	Up to 16.8 miles (U1) Up to 19.2 miles (40 Ah)
Component Weights ⁵	Base: 230 lbs. ⁷ w/ U1 batteries 245 lbs. ⁷ w/ 40 Ah batteries U1 Batteries: 24.5 lbs. 40 Ah Batteries: 32.5 lbs.
Drivetrain	Two motor, Mid-Wheel 6® drive
Brakes	"Intelligent Braking" (electronic regenerative, disc park brakes)
Batteries ³	(2) 12 V, U1 (2) 40 Ah
Battery Charger	5 amp, off-board
Warranty	5 year limited (frame) 13 month limited (batteries) 2-year limited (electronics) 2-year limited (drive motors)

- 1) Varies with base model, user weight, terrain type, battery amp hour rating (Ah), battery charge, battery condition, motors, controller type, tire type, and tire condition. This specification is subject to a variance of +10%, -5%.
- 2) Due to manufacturing tolerances and continued product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 3) AGM or Gel-Cell type required.
- 4) Tested in accordance with ANSI/RESNA, WC Vol. 2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 5) Battery weight may vary based on manufacturer.
- 6) Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.
- 7) Weight includes base, seat and batteries.

NOTE: This product conforms to all applicable ANSI/RESNA, ISO 7176 Series, and EN12184 standards. All specifications subject to change without notice.

Features

- High-back or solid seat pan, reclining sport seat. Seat features 2" of depth adjustment. Choice of 16"x16"-18", 18"x18"-20", 20"x20"-22"
- Angle, height and depth adjustable armrests
- Height adjustable headrest
- Newly re-designed lift features:
 - · 2" lower seat-to-floor height enables user to sit comfortably at standard table or desk
 - Raises and lowers in 11 seconds by using the Air switch or VR2 joystick
- Additional features:
 - Drives safely at 4 mph in the elevated position
 - Integrated foot platform with rubber insert and 2"of adjustment for comfort in both the elevated and non-elevated positions
 - Bright LED marker lights offers visibility in all environments
 - Active-Trac® suspension with Mid-Wheel 6® Technology allows for maximum stability in the elevated and non-elevated positions
 - Single-sided caster forks for enhanced maneuverability
 - Under seat front storage compartment keeps your valuables accessible

Accessories

- Alternative size lap belt: 50", 70", 80", 90"
- Cane holder
- Cell phone holder
- Cup holder
- Flag kit
- Oxygen holder
- Rear basket
- Saddle bag
- Walker holder
- Weather cover
- · XLR USB Charger

Options

- 40 amp AGM batteries
- Swing-away joystick
- · Stealth® cushions

Colors









Ruby

Red

Silver Blue

New Matte Colors



White





Black Pearl

Tanzanite

Garnet Red

Authorized Pride® Provider