

<b>Part Number</b>	Jazzy® Carbon 27
<b>Weight Capacity</b>	300 lb (136 kg)
<b>Maximum Speed</b> <sup>1</sup>	Up to 3.7 mph (6 km/h)
<b>Ground Clearance</b>	1" at motor
<b>Turning Radius</b> <sup>3</sup>	33.4"
<b>Overall Length</b> <sup>3</sup>	35.04"
<b>Overall Width</b> <sup>3</sup>	23.43"
<b>Suspension</b>	Full front suspension only
<b>Tires</b>	Front: 6.89" Rear: 8.27"
<b>Range Per Charge</b> <sup>1,2</sup>	Up to 9.3 mi (15 km) per charge
<b>Product Weight</b> <sup>5</sup>	27.8 lb (12.6 kg) without batteries
<b>Standard Seating</b> <sup>3</sup>	Dimensions: 17" x 17" Solid seat pan with removable cushion and mesh sling back.
<b>Motor Output</b>	24V 250W Brushless
<b>Brakes</b>	Regenerative and electromechanical
<b>Battery Requirements</b>	Type: Lithium battery box Size: 12 Ah
<b>Battery Weight</b> <sup>4</sup>	4.63 lb (2.1 kg)
<b>Battery Charger</b>	2 Amp
<b>Warranty</b>	5 years on frame; 1-year limited on electronics 1-year limited on drivetrain 1-year limited on battery

#### POWER MOBILITY DEVICE, NOT CODED BY DME PDAC

- 1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
- 2 Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 4 Battery weight may vary based on battery manufacturer.
- 5 Includes standard seat weight. Product is in transport configuration with all accessories removed.
- 6 Transport of this item on commercial aircraft is subject to federal regulations and airline policy. Please consult your airline prior to your flight to ensure your Jazzy® and its lithium-ion battery are permitted.

\* Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.

**PRIDE  
MOBILITY®**  
pridemobility.ca



**Jazzy®**  
carbon 27

Only  
**27 lb**  
without  
the battery



FDA Class II  
Medical Device\*

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.



Folds and  
unfolds for easy  
transport



Quick disconnect  
joystick with  
LCD display



Full Front  
suspension



**Jazzy**  
carbon 27

Only  
**27<sup>lb</sup>**  
without  
the battery

Shown with Optional headrest  
and attendant control bracket

17" x 17" seat  
with mesh back

Adjustable lap belt  
(Standard)

Durable, lightweight  
carbon fiber frame

Mesh cup holder  
(Standard)

Compact base  
folds & unfolds  
for easy transport

FDA Class II  
Medical Device\*

Full front  
suspension



Folds & unfolds  
for easy transport



Quick disconnect joystick  
with LCD display



Full front  
suspension



Accessible  
release lever

## PUSHING THE BOUNDARIES OF LIGHTWEIGHT MOBILITY

World-class Jazzy® quality meets the next level of lightweight mobility with the Jazzy® Carbon 27 power wheelchair. Weighing 27.8 lb without the battery, this unit folds in a flash, is easily stored for convenient transport, and has a durable carbon fiber frame. Travel worry-free and experience the Jazzy® difference.

### FEATURES

- Unit weighs 27.8 lb without the battery
- 300-lb weight capacity
- Durable carbon fiber frame
- Front suspension for a smooth ride
- Top speed of up to 3.7 mph
- Folded dimensions: 27.5" L x 23.4" W x 13.38 H
- Easy fold design secures with a strap
- Flip-up lightweight carbon fiber arms
- 12 Ah, airline compliant<sup>6</sup> lithium battery (standard)
- Quick disconnect LCD joystick with digital display, USB-C port & puddle light
- Mesh cup holder, cushion and adjustable lap belt (standard)

### COLOR OPTIONS

Night Sky

Black Cherry

### OPTIONAL ACCESSORIES

- Additional 12 Ah airline-compliant lithium battery<sup>6</sup>
- Additional 20 Ah lithium battery (COMING SOON)
- Battery travel bag
- Headrest and attendant control bracket
- Jazzy® travel kit
- Simplicity cushion – 17" x 17"



Jazzy® travel kit shown