

Part Number	Jazzy® Carbon 27X
Weight Capacity	275 lb (125 kg)
Maximum Speed ¹	Up to 3.7 mph (6 km/h)
Ground Clearance ³	3"
Turning Radius ³	35.4"
Overall Length ³	41"
Overall Width ³	23"
Seat-to-Ground Height	19.7"
Tires	Front: 7" Rear: 10"
Range Per Charge ^{1,2}	Up to 8 mi (13 km) per charge
Product Weight ^{3,5}	27.5 lb (13 kg) without batteries
Standard Seating ³	Padded sling back with sling seat and removable seat cushion
Folded Dimensions	34" L x 12" W x 30" H
Motor Output	24V 150W
Brakes	Regenerative and electromechanical
Battery Requirements ⁶	24V 10Ah
Battery Weight ⁴	3.5 lb (1.6 kg)
Battery Charger	3.5A
Warranty	5-year limited on frame; 2-year limited on electronics 2-year limited on drivetrain 13-month on battery

POWER MOBILITY DEVICE, NOT CODED BY DME PDAC

1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.

2 Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.

4 Battery weight may vary based on battery manufacturer.

5 Includes standard seat weight. Product is in transport configuration with all accessories removed.

6 Transport of this item on commercial aircraft is subject to federal regulations and airline policy. Please consult your airline prior to your flight to ensure your Jazzy® and its lithium-ion battery are permitted.

* Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.

**PRIDE
MOBILITY®**

pridemobility.ca



NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

Jazzy®
carbon 27X

Weighs just
27lb
without the
battery



FDA Class II
Medical Device*



FDA Class II
Medical Device*

CA_JazzyCarbo27X_Trifold_011926

Jazzy[®] carbon 27X

Weighs just
27lb
without the
battery



Programmable quick
disconnect joystick

Airline compliant⁶
lithium battery
(10Ah only)

Durable, lightweight
carbon fiber frame

Adjustable
foot rests

FDA Class II
Medical Device*

Storage pocket
in seat back

Adjustable lap belt
(Standard)

Removable cushion
(Standard)

Carbon fiber
cup holder
(Standard)

X-frame design
folds vertically
for easy transport

YOUR MOBILITY X-FACTOR.

The Jazzy[®] Carbon 27X folds vertically in just a few seconds thanks to its X-frame design, and it weighs a mere 27.5 lb without the battery. Comprised of super-strong carbon fiber, this power wheelchair is durable enough to withstand all your adventures and light enough to keep you moving.

FEATURES

- Unit weighs 27.5 lb without the battery
- 275-lb weight capacity
- Top speed up to 3.7 mph
- X-frame design allows the unit to fold nearly flat in vertical position
- Advanced carbon fiber frame offers durability and light weight
- Folded dimensions: 34"L x 12"W x 30"H
- Flip up lightweight carbon fiber arms
- Adjustable Foot Rests
- Airline-compliant lithium battery** (10 Ah only)
- Storage pocket in seat back
- Cushion and adjustable lap belt (standard)
- Programmable quick disconnect joystick (standard)
- Carbon fiber Storage with Cup Holder (standard)
- Collapsible handles for easy storage

COLOR OPTIONS

Black



Folds & unfolds
for easy transport



Advanced durable
carbon fiber frame



X-frame design for
vertical folding



Latch for locking unit
open or closed

OPTIONAL ACCESSORIES

Additional 10 Ah
airline-compliant battery**
Battery Travel Bag
Jazzy Travel Kit (Coming Soon)
Military Patches
Simplicity Cushion 18" x 18"
Mesh Cup Holder