

Model	Jazzy Elite HD
Weight Capacity	450 lbs. / 204 kg.
Turning Radius ²	24"
Overall Width ²	24.25"
Overall Length ²	35.5"
Maximum Speed ¹	4 mph / 6.4 km/h
Ground Clearance ²	3"
Front Wheels ²	3" solid anti-tips
Drive Wheels ²	14" foam filled
Rear Wheels ²	6" solid casters
Drivetrain	Two motor, front-wheel drive
Braking System	Regenerative and electro-mechanical
Suspension Type	Limited
High-back Seat Sizes ²	Up to 24" x 24"
High-Back Seat-to-Floor Range ²	20" - 22"
Synergy Seating ²	Up to 24" x 24"
Synergy Seat-to-Floor Range ²	15.75" - 17.75"
Per Charge Range ^{1,4}	Up to 14.4 ml / 23 km
Standard Electronics	50A, Dynamic LiNX
Component Weights ^{2,5}	Base 99 lbs. / 45 kg. Seats: 18" wide - 51.8 lbs. / 23.6 kg. 20" wide - 54 lbs. / 24.5 kg. 22" wide - 57.6 lbs. / 26.2 kg. Batteries 37 lbs. / 16.8 kg. each
Battery Requirements ³	(2) 12 volt, NF-22 batteries
Battery Charger	5 amp, off-board
Warranty	5-year limited on frame; 24-months limited on electronics; 24-months limited on drive motors; 13-months on batteries.
Warranty	Frame: Lifetime limited Drivetrain: 13-month Electronics: 13-month Battery: 13-month

- 1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- 2 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- 3 AGM or gel cell type required.
- 4 Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 5 Battery weight may vary based on manufacturer.



Beamsville, Ontario • 888-570-1113

pridemobility.ca |    

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

The information contained herein is correct at time of publication; we reserve the right to alter specifications without prior notice.

0004874-113022



Jazzy
ELITE HD

**HEAVY DUTY
FRONT-WHEEL
DRIVE**

FDA Class II
Medical Device

**NORTH AMERICA'S
FAVORITE POWER CHAIR**

pridemobility.ca |    

ENHANCED OUTDOOR PERFORMANCE



The **Jazzy® Elite HD** features two-motor, front-wheel drive technology and heavy-duty construction. The Jazzy Elite HD features 14" drive wheels for enhanced outdoor performance combined with a front-wheel drive design to provide tight turns for optimal indoor maneuverability in small spaces. Large front wheels add extra absorption for superb climbing capabilities and combine with front anti-tip wheels for handling various terrains with ease.

- Universal, fixed, in-line joystick mount
- Black, flat-free, non-marking tires
- Lap belt

Colors



Optional Seating

Black vinyl, high-back seat, available sizes:

- Up to 24"W x 24"D

Synergy Seating, available sizes:

- Up to 24"W x 24"D

Options

- Bell shaped foot platform
- Elevating leg rests
- Heel loops
- Height adjustable foot platform bracket
- Residual limb support
- Swing-away legrests
- Swing-away joystick

Accessories

- Cane/crutch holder
- Cell phone holder
- Cup holder
- Flag holder
- Lap belt (50", 70", 80", 90")
- Military patches
- Oxygen tank holder
- Rear basket
- Saddle bag
- Walker holder
- Weather cover
- XLR USB charger



Front-wheel drive design



Easy-access free wheel levers



50A, Dynamic LiNX controller