	Jazzy® EVO 613	
Models	Jazzy® EVO 613 Jazzy® EVO 613 Li	
Weight Capacity	300 lbs.	
Turning Radius ³	22"	
Overall Width ³	22″	
Overall Length ³	41" with footrest (35" without)	
Maximum Speed ²	Up to 5.1 mph	
Ground Clearance ³	2.5" at motor	
Front Casters	5" solid	
Drive Wheels	13" solid	
Rear Casters	5″ solid	
Drivetrain	Two motor, in-line, mid-wheel drive	
Braking System	Regenerative and electromechanical	
Suspension Type	Active-Trac [®] (ATX)	
Seat Sizes	High-back: 16" × 16-18" 18" × 18-20" 20" × 18-20"	Synergy®: 16-20″ x 16-20″
Seat-to-Floor Range	High-back: 19.75" - 21.25" Synergy®: 17.5" - 18.5"	
Per Charge Range ^{2,5}	EVO 613: 14 miles @ 300 lbs. 17.1 miles @ 200 lbs. EVO 613 Li: 17.6 miles @ 300 lbs. 21 miles @ 200 lbs.	
Standard Electronics	40A, PG GC3	
Component Weights	Base: 94 lbs. High-back seat: 48 lbs. (18" seat) Synergy [®] seat: 62.8 lbs. ⁷ (2) U1 Batteries: 23.5 lbs. each ⁶ (1) Lithium battery: 18.5 lbs. ⁶	
Total Weight ¹	EVO 613: 189 lbs. EVO 613 Li: 160.5 lbs.	
Battery Requirements ⁴	EVO 613: (2) 12 volt, U-1 EVO 613 Li: (1) Lithium iron phosphate, 25.6V, 32.4Ah	
Battery Charger	3.5 amp, off-board	
Warranty	Frame: 5-year limited Drivetrain: 2-year limited Electronics: 2-year limited SLA Battery: 13-month limited Lithium Battery: 3-year limited	

1 Total weight includes base, high-back seat and batteries.

- 2 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- 4 AGM or gel cell type are required for U-1 batteries.
- 5 Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 6 Battery weight may vary based on manufacturer.
- 7 Weight listed includes TRU-Comfort seating, armrests and controller.



Jazzy[®] EVO 613 shown in Sugar Plum with Synergy[®] seating system



Beamsville, Ontario • 888-570-1113 pridemobility.ca | 🖸 🗗 🕑 🚱

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or -3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.



OUR NARROWEST **JAZZY[®] YET!**



ALWAYS INSIST ON LITHIUM! pridemobility.ca

OUR NARROWEST JAZZY®



Always insist on lithium!

The all new Jazzy[®] EVO 613 series is here! The Jazzy[®] EVO 613 series provides the freedom to navigate small spaces and tight corners indoors with a narrow base width of only 22". The Jazzy[®] EVO 613 Li uses a lithium battery that provides longer life and more convenient charging.

Features

- Maximum speed up to 5.1 mph
- 13" mid-wheel drive tires
- Narrow base width of only 22"
- Removable matte-finish color-through shrouds
- Access to the controller and batteries without the user needing to exit the chair
- Depth-adjustable, limited recline, memory foam single seat post allows the user to swivel
- Pride's patented Active-Trac[®] (ATX) Suspension
- Mid-Wheel 6[®] Technology

Available Colors



Accessories

- Cane holder
- Cell phone holder
- Cup holder
- Flag kit
- Oxygen tank holder
- Rear basket
- Saddlebag
- Walker holder
- Weather cover
- 2.1A XLR USB charger

Options

- Elevated leg rests
- Swing-away leg rests
- Unoccupied transit kit
- Solid seat pan
- Synergy[®] seating