

Model	Jazzy® EVO 614HD ¹
Weight Capacity	450 lbs. / 204 kg.
Turning Radius ³	22"
Overall Width ³	27"
Overall Length ³	40.5"
Maximum Speed ²	Up to 4.88 mph / 7.85 km/h
Ground Clearance ³	3"
Front Casters ³	6"
Drive Wheels ³	13.75" perforated
Rear Casters ³	6"
Drivetrain	Two in-line HD motors
Braking System	Regenerative and electromechanical
Suspension Type	Full independent front and rear Active-Trac® ATX Suspension
High-back Seat Sizes ³	18" x 18-20" 20" x 18-20" 22" x 20-22"
Seat-to-Floor Range ³	22-24"
Per Charge Range ^{2,5}	Up to 13 miles / 20.9 km at 450 lbs. / 204 kg Up to 21.5 miles / 34.6 km at 200 lbs. / 90.7 kg
Standard Electronics	60A Dynamic Linx
Component Weights ^{3,6}	Base: 125.7 lbs. / 57 kg Seat: 57.8 lbs. / 26.2 kg Batteries: (2) 36 lbs. / 16.3 kg each
Total Weight ^{3,7}	255.5 lbs. / 115.9 kg
Battery Requirements	(2) NF-22
Battery Charger	5 amp, off-board
Warranty	Frame: 5-year limited Drivetrain: 24-month limited Electronics: 24-month limited Battery: 13-month limited

¹ Jazzy® Power Chairs are FDA Class II Medical Devices designed to aid individuals with mobility impairments.

² Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).

³ Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

⁴ AGM or gel cell type required.

⁵ Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

⁶ Battery weight may vary based on manufacturer.

⁷ Total weight includes base, seat and batteries.

PRIDE
MOBILITY
LIVE YOUR BEST®

(Canada) 888-570-1113

pridemobility.ca |    

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

0004778-011623

Jazzy®
EVO 614HD

**OUR MOST ADVANCED
HEAVY DUTY POWER
CHAIR YET!**



FDA Class II
Medical Device¹

CANADA'S FAVORITE POWER CHAIR

pridemobility.ca

OUR MOST ADVANCED HEAVY DUTY POWER CHAIR YET!

Height, width, depth
and angle adjustable
armrests

Optional
swing-away
joystick

Depth-adjustable,
limited recline,
memory foam
high-back seat
with swivel

Removable
color-through
shrouds

Heavy duty
seat post

13.75" perforated
drive wheels

Access to controller
and batteries without
the user needing to
exit the chair



Front and rear independent
Active-Trac® ATX Suspension

FDA Class II
Medical Device¹

The Future of Heavy Duty Mobility

The Jazzy® EVO 614HD has tall 13.75-inch drive wheels, the centerpiece for Pride's beloved Mid-Wheel 6® platform. Active-Trac® ATX Suspension adapts to the terrain, so grassy fields and gravel paths feel smoother than they ever have before. And now, with a new memory foam seat, your Jazzy EVO 614HD is ready to provide all-day comfort right out of the box.

- 3" ground clearance
- 13.75" perforated mid-wheel drive tires
- Front and rear independent Active-Trac® ATX Suspension
- Removable matte-finish color-through shrouds
- Depth-adjustable, limited recline, memory foam heavy duty single seat post allows the user to swivel
- Access to controller and batteries without the user needing to exit the chair
- Standard lap belt (60" - 18" or 20" W, 70" - 20" W)

Available Colors



Accessories

- Cane/crutch holder
- Cell phone holder
- Cup holder
- Flag kit
- Oxygen tank holder
- Rear basket
- Saddlebag
- Stealth® cushions
- Walker holder
- Weather cover
- 2.1A XLR USB charger
- 50-90" lap belts

Options

- Elevated leg rests
- Solid seat pan, limited recline, high-back with swivel
- Heel loops
- Swing-away joystick
- Swing-away leg rests
- Unoccupied transit
- Residual limb support



Front and rear independent
Active-Trac® ATX Suspension



Large footplate with removable
rubber cover for easy cleaning



High-back memory foam seat



Scan QR code
for more info.