

VivaLift!®
Atlas Plus
Collection
PLR-2985M

MEDIUM



VivaLift!®
Elegance
Collection
PLR-975L

LARGE



Heritage
Collection
LC-358XXL

XXL



Sizing Guide

The complete reference guide for
choosing the perfect size
power lift recliner



PRIDE
MOBILITY
LIVE YOUR BEST®



pridemobility.ca

VivaLIFT!
POWER RECLINER COLLECTION

PRIDE
POWER RECLINER
COLLECTION

A Comfy Fit

Our power lift recliners are as varied in size and shape as the people who enjoy them. Taking the following detailed measurements prior to a purchase will ensure a proper fit and deliver maximum satisfaction.

USER HEIGHT

USER WEIGHT

SEAT TO TOP OF BACK

Measure the vertical distance from the base of a seat to the top of the user's head while seated.

SEAT DEPTH

Measure the horizontal distance from the back of a seat to the back of the user's knee while seated.

FLOOR TO TOP OF SEAT

Measure the vertical distance from the floor to the top of the back of the user's knee while seated.

WIDTH BETWEEN ARMS

Measure the horizontal distance between the hips while seated (add a couple inches for wiggle room).



VIVALIFT!® SIZES*

CHAIR SIZE

USER HEIGHT

Small/Petite Wide 5'3" and below

Medium 5'4"-5'9"

Large/Large Tall 5'10" and above

PRIDE® POWER LIFT RECLINER SIZES*

CHAIR SIZE

USER HEIGHT

Super Petite 5'2" and below

Small/Petite Wide 5'3" and below

Medium 5'4"-5'9"

Large/Large Tall 5'10" and above

XL 5'6" - 6'1"

XXL 5'6" - 6'1"

* All dimensions and sizing estimates are recommendations and do not necessarily represent the requirements of every individual user.

VivaLIFT!
POWER RECLINER COLLECTION



TECHNOLOGY

Our lift chairs are equipped with innovative, state-of-the-art technology that meets your comfort needs. VivaLift!® Power Recliners feature a power headrest and power lumbar for a truly relaxing experience.

COMFORT

Experience customizable comfort at the touch of a button with the Elegance Collection! Stretch out comfortably with the footrest extension and enjoy peace of mind with the lithium battery backup.

DURABILITY

Engineered to last, each lift chair features a durable, all-steel lift frame backed by a lifetime warranty. A heavy-duty actuator and scissor mechanism deliver quiet and smooth performance. Plus, a premium furniture-grade hardwood frame offers a lifetime of strength and durability.