

PERFORMANCE FABRICS



Uf ULTRAFABRICS® ULTRALEATHER™

Ultrafabrics® Ultraleather™ on Pride® Power Lift Recliners offer a richly detailed grain pattern expressive of premium grade leather. Ultrafabrics' polyurethane products are manufactured utilizing Tukumi™ technology, a 4-layer laboratory process to create a luxurious, durable and comfortable polyurethane fabric.

FEATURES AND BENEFITS

- Rich feel and natural grain pattern expressive of soft luxury-grade leather
- Clinically proven more resistant to temperature changes than vinyl and genuine leather for comfort
- Antimicrobial surface cleans effortlessly with soap and water



CRYPTON® | super fabric

Crypton® Super Fabrics from Pride® Power Lift Recliners are not only luxurious, durable and contemporary, but also resist spills, stains, odors and bacteria. Nothing gets past Crypton's unique patented integrated fiber process, protecting your upholstery inside and out.

FEATURES AND BENEFITS

- 100% water resistant – spills just wipe away
- Only disinfectable fabrics in the U.S. approved by the EPA
- Permanent technology that will not come off the fabrics
- 5-year warranty on moisture barrier



staKleen® PERFORMANCE FABRIC

Lexis Sta-Kleen® Performance Fabric is a revolution in faux leather upholstery – the first stain resistant urethane fabric. Sta-Kleen offers permanent protection, created with a proprietary bonding process inherent in the manufacturing process of the product. It is not a topical application; it will last the lifetime of the upholstery.

FEATURES AND BENEFITS

- Look and feel of real leather, but with care-free durability
- Exceptional strength, abrasion resistance, and UV stability
- "Stain shield" allows spills to wipe off with a cloth, and it's anti-bacterial



QUALITY CONSTRUCTION



Pride's exclusive high-quality power lift recliner designs offer:

- Premium furniture-grade laminate/hardwood frames for added strength and durability
- Heavy-duty lift actuator and scissor mechanism for long chair life
- Quiet and smooth lift system
- High resiliency fiber-fill for maximum comfort
- Fabric options to fit any need and décor

In addition, Pride's industry leading power lift recliner designs are easy to use, available in a wide variety of stylish, high-quality fabric selections and provide the highest level of durability and longevity available.

HAND CONTROLS						
EXAMPLES OF APPLICABLE HAND CONTROLS AS INDICATED IN THE GRID			OPTIONAL HAND CONTROLS			
EASY TOUCH	EASY TOUCH WITH USB PORT	EASY TOUCH INFINITE-POSITION WITH USB PORT	TOGGLE	TOGGLE WITH HEAT AND MASSAGE	DUAL-MOTOR DELUXE HEAT AND MASSAGE	STANDARD HEAT AND MASSAGE

PRIDE® POWER RECLINER COLLECTION



North America's #1 Power Lift Recliner

PRIDE® POWER RECLINER COLLECTION

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

Beamsville, Ontario | 888-570-1113 | pridemobility.ca |

North America's #1 Power Lift Recliner

PRIDE
POWER
RECLINER
COLLECTION

ELEGANCE
COLLECTION


HERITAGE
COLLECTION

essential
COLLECTION



MODEL	LC-108	LC-450C	LC-570S	LC-570M	LC-570L	LC-358PW	LC-358S	LC-358M	LC-358L	LC-358XL	LC-358XXL	LC-102	LC-106	LC-107
Size	Medium	Medium	Small	Medium	Large	Petite Wide	Small	Medium	Large	Heavy-Duty	Very Heavy-Duty	Super Petite	Medium	Medium
Position	Infinite-position, Trendelenburg Chaise Lounger 5'4" - 6'	3-position Chaise Lounger 5'4"-5'9"	3-position Chaise Lounger 5'3" and below	3-position Chaise Lounger 5'4"-5'9"	3-position Chaise Lounger 5'10" and above	3-position Chaise Lounger 5'3" and below	3-position Chaise Lounger 5'3" and below	3-position Chaise Lounger 5'4"-5'9"	3-position Chaise Lounger 5'10" and above	3-position Chaise Lounger 5'6"-6'1"	2-position Chaise Lounger 5'6"-6'1"	3-position Chaise Lounger 5'2" and below	3-position Chaise Lounger 5'4"-5'9"	Infinite-position, Trendelenburg Chaise Lounger 5'4"-5'9"
Weight Capacity	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.	500 lbs.	600 lbs.	250 lbs.	325 lbs.	325 lbs.
Seat-to-Floor	19"	20"	19"	20"	21"	19"	19"	20"	21"	20"	20"	16"	18"	20"
Seat Depth	19"	20"	19"	20"	22"	19"	19"	20"	22"	22"	22"	18"	19"	20"
Seat Width	20"	20"	19"	20"	22"	24"	18"	20"	22"	26"	30"	16"	20"	20"
Top of Back to Seat	29"	27"	24.5"	26.5"	27"	23.5"	24"	24.5"	25.5"	29"	29"	23.5"	27.5"	27.5"
Back Style	Sewn pillow back (back upholstery is fully removable and adjustable) (lumbar pillow standard)	Blown fiber split back (back upholstery is fully removable and adjustable)	Blown fiber pillow back	Blown fiber waterfall back	Blown fiber waterfall back	Blown fiber button back	Blown fiber button back	Blown fiber button back	Blown fiber button back	Blown fiber button back	Blown fiber button back	Dual pillow back	Sewn pillow back	Sewn pillow back (lumbar pillow standard)
Fully Padded Chaise	YES (seat pad is fully removable and adjustable)	YES (seat pad is fully removable and adjustable)	YES	YES	YES	YES	YES	YES	YES	YES (removable, fully padded chaise)	YES (removable, fully padded chaise)	YES	YES	YES
Footrest Extension	YES (standard)	YES (optional)	YES (optional)	YES (optional)	YES (standard)	YES (standard)	YES (standard)	YES (standard)	YES (standard)	N/A	N/A	N/A	N/A	YES (standard)
Heat and Massage	N/A	YES (optional)	YES (optional)	YES (optional)	YES (optional)	YES (optional)	YES (optional)	YES (optional)	YES (optional)	YES (optional, toggle hand control needed)	N/A	N/A	N/A	N/A
Seat Type	Sinuuous spring	Total Comfort ¹	Total Comfort ¹	Total Comfort ¹	Total Comfort ¹	Total Comfort ¹	Total Comfort ¹	Total Comfort ¹	Total Comfort ¹	Sinuuous spring	Sinuuous spring	Sinuuous spring	Sinuuous spring	Sinuuous spring
Distance from Wall	24"	18"	16.5"	19"	24"	18"	18"	17"	19"	19"	14"	16"	18"	26"
Overall Width	31.5"	37"	35"	36"	38"	40"	35"	35"	37"	39.5"	43"	25.5"	32.75"	35"
Chair Weight	136 lbs.	119 lbs.	110 lbs.	113 lbs.	114 lbs.	130 lbs.	118 lbs.	117 lbs.	120 lbs.	159 lbs.	168 lbs.	96 lbs.	96 lbs.	116 lbs.
Hand Control	Easy touch Infinite-position with USB port	Easy touch with USB port	Easy touch with USB port	Easy touch with USB port	Easy touch with USB port	Easy touch with USB port	Easy touch with USB port	Easy touch with USB port	Easy touch with USB port	Easy touch with USB port	Easy touch with USB port	Easy touch	Easy touch	Easy touch Infinite-position with USB port
Warranty	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	1-Yr. Limited	1-Yr. Limited	1-Yr. Limited

FABRIC SELECTION GUIDE



ADDITIONAL NOTES:

- All chairs have head and arm covers standard
- All chairs have an engineered furniture-grade laminate frame
- Lifetime Limited on lift, recline and scissor mechanism and wood frame on all power lift recliners (Exceptions: Essential Collection (1-year), LC-358XL and LC-358XXL (7-year))
- Transformer features self-diagnostic electronics with 9' cord
- All chairs have integrated emergency battery backup (batteries not included)
- All chairs can have fiber fill added or removed for custom comfort adjustment
- All chairs are shipped in an impact-resistant, reinforced, corrugated box

ICON KEY

Indicates shortest available lead times

PARTNER TRADEMARKS

POSITION GUIDE (The lifted, seated and reading positions are achievable on all Pride® Power Lift Recliners.)

PRODUCT VIEW	LIFTED	SEATED	READING	NAPPING	RELAXING	INFINITE POSITION	TRENDELENBERG
2-Position	The 2-position TV recline is achievable on all Pride® Power Lift Recliners. This is a comfortable position for performing many relaxing activities including reading, watching TV, working on your tablet or for conversation with family and friends.						
3-Position	The 3-position full recline is available on most Pride® Power Lift Recliner models and puts you in a position ideal for napping or relaxing.						
Infinite Position	The infinite-position lay flat recline is perfect for taking a nap. Infinite-position Pride® Power Lift Recliners enable the backrest and footrest to move independently from one another so that a variety of positions can be reached. Infinite positioning allows for easier social engagement, eating and crafting activities while having the footrest fully elevated.						
Trendelenburg	The Trendelenburg position, elevates one's feet higher than his or her head for increased comfort and relaxation.						

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.