

PERFORMANCE FABRICS



Uf ULTRAFABRICS® ULTRALEATHER™

Ultrafabrics' UltraLeather™ on Pride® Power Lift Recliners offer a richly detailed grain pattern expressive of premium grade leather. Ultrafabrics' polyurethane products are manufactured utilizing Tukumi™ technology, a 4-layer laboratory process to create a luxurious, durable and comfortable polyurethane fabric.

FEATURES AND BENEFITS

- Rich feel and natural grain pattern expressive of soft luxury-grade leather
- Clinically proven more resistant to temperature changes than vinyl and genuine leather for comfort
- Antimicrobial surface cleans effortlessly with soap and water



CRYPTON® super fabric

Crypton® Super Fabrics from Pride® Power Lift Recliners are not only luxurious, durable and contemporary, but also resist spills, stains, odors and bacteria. Nothing gets past Crypton's unique patented integrated fiber process, protecting your upholstery inside and out.

FEATURES AND BENEFITS

- 100% water resistant – spills just wipe away
- Only disinfectable fabrics in the U.S. approved by the EPA
- Permanent technology that will not come off the fabrics
- 5-year warranty on moisture barrier



staKleen® PERFORMANCE FABRIC

Lexis Sta-Kleen® Performance Fabric is a revolution in faux leather upholstery – the first stain resistant urethane fabric. Sta-Kleen offers permanent protection, created with a proprietary bonding process inherent in the manufacturing process of the product. It is not a topical application; it will last the lifetime of the upholstery.

FEATURES AND BENEFITS

- Look and feel of real leather, but with care-free durability
- Exceptional strength, abrasion resistance, and UV stability
- "Stain shield" allows spills to wipe off with a cloth, and it's anti-bacterial



QUALITY CONSTRUCTION



Pride's exclusive high-quality power lift recliner designs offer:

- Premium furniture-grade laminate/hardwood frames for added strength and durability
- Heavy-duty lift actuator and scissor mechanism for long chair life
- Quiet and smooth lift system
- High resiliency fiber-fill for maximum comfort
- Fabric options to fit any need and décor

In addition, Pride's industry leading power lift recliner designs are easy to use, available in a wide variety of stylish, high-quality fabric selections and provide the highest level of durability and longevity available.

HAND CONTROLS						
EXAMPLES OF APPLICABLE HAND CONTROLS AS INDICATED IN THE GRID			OPTIONAL HAND CONTROLS			
EASY TOUCH	EASY TOUCH WITH USB PORT	EASY TOUCH INFINITE-POSITION WITH USB PORT	TOGGLE	TOGGLE WITH HEAT AND MASSAGE	DUAL-MOTOR DELUXE HEAT AND MASSAGE	STANDARD HEAT AND MASSAGE

PRIDE® POWER RECLINER COLLECTION



North America's #1 Power Lift Recliner

PRIDE® POWER RECLINER COLLECTION

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

Beamsville, Ontario | 888-570-1113 | pridemobility.ca |

North America's #1 Power Lift Recliner

PRIDE
POWER
RECLINER
COLLECTION

ELEGANCE
COLLECTION

HERITAGE
COLLECTION

essential
COLLECTION



COLLECTION	ELEGANCE COLLECTION				HERITAGE COLLECTION						ESSENTIAL COLLECTION						
SIZE	Medium	Small	Medium	Large	Petite/Wide	Small	Medium	Large	Heavy Duty	Very Heavy Duty	Super Petite	Medium	Medium	Medium			
MODEL NUMBER	LC-108	LC-570S	LC-570M	LC-570L	LC-358PW	LC-358S	LC-358M	LC-358L	LC-358XL	LC-358XXL	LC-102	LC-104	LC-106	LC-107			
POSITION	Infinite Position	3-position	3-position	3-position	3-position	3-position	3-position	3-position	3-position	2-Position	3-position	3-position	3-position	Infinite Position			
STANDARD HAND CONTROL	Infinite-Potion Easy Touch with USB	Easy Touch with USB	Easy Touch with USB	Easy Touch with USB	Easy Touch with USB	Easy Touch with USB	Easy Touch with USB	Easy Touch with USB	Easy Touch with USB	Easy Touch with USB	Easy Touch	Easy Touch with USB	Easy Touch with USB	Infinite-Potion Easy Touch with USB			
BACK TYPE	Sewn pillow back (back upholstery is fully removable and adjustable)	Blown fiber pillow back (Pillow are individually accessible to adjust fill)	Blown fiber pillow back (Pillow are individually accessible to adjust fill)	Blown fiber pillow back (Pillow are individually accessible to adjust fill)	Blown fiber button back	Blown fiber button back	Blown fiber button back	Blown fiber button back	Blown fiber button back	Blown fiber button back	Dual pillow back	Sewn pillow back	Sewn pillow back	Sewn pillow back (lumbar pillow standard)			
SEAT TYPE	Sinuus spring (seat pad is fully removable and adjustable)	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Sinuus spring (seat pad is fully removable and adjustable)	Sinuus spring (seat pad is fully removable and adjustable)	Sinuus spring	Sinuus spring	Sinuus spring			
FOOTREST EXTENSION	Standard	Optional	Optional	Standard	Standard	N/A	Standard	Standard	N/A	N/A	N/A	N/A	N/A	Standard			
ADDITIONAL FEATURES	Standard Head and Arm Covers, and Lumbar Roll	Standard Head and Arm Covers	Standard Head and Arm Covers	Standard Head and Arm Covers	Standard Head and Arm Covers	Standard Head and Arm Covers	Standard Head and Arm Covers	Standard Head and Arm Covers	Standard Head and Arm Covers	Standard Head and Arm Covers	Standard Head and Arm Covers	N/A	Standard Head and Arm Covers	Standard Head and Arm Covers, and Lumbar Roll			
USER HEIGHT RANGES**	5'4" - 6"	5'3" and below	5'4" - 5'9"	5'10" and above	5'3" and below	5'3" and below	5'4" - 5'9"	5'10" and above	5'6" - 6'1"	5'6" - 6'1"	5'2" and below	5'4" - 5'9"	5'4" - 5'9"	5'4" - 5'9"			
SEAT-TO-FLOOR	19"	18.5"	19.5"	20.5"	19"	19"	20"	21"	20"	20"	16.5"	17.5"	18"	19"			
SEAT WIDTH	19.5"	17"	18.5"	20.5"	24"	18.5"	20"	22.5"	26"	30"	17"	19.5"	20"	20.5"			
SEAT DEPTH	20"	19.5"	19.5"	21.5"	18.5"	19"	20"	21.5"	22"	22"	16.5"	18"	19"	19"			
SEAT DEPTH "WITH FOOTREST ELEVATED"	38.5"	36.5"	36.5"	42"	38.5"	38.5"	40"	41"	38.5"	39"	32.5"	33.5"	38"	37"			
TOP OF BACK TO SEAT	30"	23.5"	26"	27"	23.5"	24.5"	26"	25"	28"	29"	24.5"	23"	25"	25"			
CHAIR WIDTH	30.5"	33.5"	34.5"	35.5"	39.5"	33.5"	33.5"	36"	39.5"	43"	26"	29.5"	34"	34.5"			
CHAIR DEPTH	36.5"	37"	37.5"	40"	36.5"	37"	37"	37"	41"	43"	33.5"	36"	36"	35.5"			
CHAIR HEIGHT	48"	41"	44.5"	45.5"	41.5"	41.5"	43.5"	44"	46"	47"	40"	39.5"	43"	43"			
DISTANCE FROM WALL	26.5"	16.5"	19"	19.5"	17"	17"	17"	18"	18.5"	10"	16"	16.5"	18"	22"			
CHAIR WEIGHT	136 lb	110 lb	113 lb	114 lb	130 lb	118 lb	117 lb	120 lb	159 lb	168 lb	106 lb	106 lb	106 lb	116 lb			
WEIGHT CAPACITY	375 lb	375 lb	375 lb	375 lb	375 lb	375 lb	375 lb	375 lb	500 lb	600 lb	250 lb	325 lb	325 lb	325 lb			
WARRANTY	7-yr Pro-Rated	7-yr Pro-Rated	7-yr Pro-Rated	7-yr Pro-Rated	7-yr Pro-Rated	7-yr Pro-Rated	7-yr Pro-Rated	7-yr Pro-Rated	7-yr Pro-Rated	7-yr Pro-Rated	1-yr Limited	1-yr Limited	1-yr Limited	1-yr Limited			
FABRIC SELECTION GUIDE	STANDARD	Cloud 9 BLACK CHERRY PACIFIC STONE WALNUT				Cloud 9 BLACK CHERRY PACIFIC STONE WALNUT						Cloud 9 STONE WALNUT		Darby JAZA MIDNIGHT PEWTER			
	OPTIONAL					Lexis Sta-Kleen™ MUSHROOM BLACK BURGUNDY CHESTNUT			Crypton Aria™ ESPRESSO SAND COOL GRAY RED		Ultraleather™ BUFF FUDGE GARNET PECAN CHARCOAL		Lexis Sta-Kleen™ CHESTNUT				

ADDITIONAL NOTES:

- All chairs have an engineered furniture-grade laminate frame
- Lifetime Limited on lift, recline and scissor mechanism and wood frame on all power lift recliners (Exceptions: Essential Collection, LC-358XL and LC-358XXL. See warranty insert for additional details.)
- Transformer features self-diagnostic electronics with 9' cord
- All chairs have integrated emergency battery backup (batteries not included)
- All chairs can have fiber fill added or removed for custom comfort adjustment
- All chairs are shipped in an impact-resistant, reinforced, corrugated box

PARTNER TRADEMARKS



POSITION GUIDE (The lifted, seated and reading positions are achievable on all Pride® Power Lift Recliners.)

PRODUCT VIEW	LIFTED	SEATED	READING	NAPPING	RELAXING	INFINITE POSITION
2-Position	These chairs will achieve the Lifted, Seated, Reading and Relaxing positions					
3-Position	The 3-position full recline chair achieves Lifted, Seated, Reading, Relaxing, and Napping positions					
Infinite Position	The infinite-position chairs offer independent adjustment of the seat and back allowing for adjustment anywhere within the achievable seat and back angles of that particular model. The infinite-position chair allows for the deepest recline available.					

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.