

Model	SC134
Weight Capacity	500 lb (226 kg)
Maximum Speed²	Up to 9.3 mph (15 km/h)
Ground Clearance^{3,4}	3.5" to motor, 4.4" at mid-deck
Turning Radius³	82.28"
Overall Length³	55.9"
Overall Width³	27.05"
Seat-to-Ground Height³	24.2" - 25.6"
Seat-to-Deck Height³	16.5" - 18.4"
Front Tires³	13" pneumatic
Rear Tires³	13" pneumatic
Suspension	Full suspension
Range Per Charge^{2,4}	Up to 23.4 mi (37.7 km) at 220 lb (100 kg)
Total Weight^{3,5}	201.2 lb (91.2 kg) heaviest piece 253.8 lb (115.1 kg) with seat 298.8 lb (135.5 kg) with batteries 351.4 lb (159.4 kg) with seat and batteries
Heaviest Piece When Disassembled⁴	N/A (Unit does not disassemble)
Standard Seating³	Type: High back Weight: 52.6 lb (23.9 kg) Material: Black vinyl Dimensions: 20" W x 18"-20" D
Drive System	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Dual Braking System	Regenerative and electromechanical
Battery Requirements⁵	Type: (2) 12V, sealed lead-acid, non-spillable deep-cycle Size: 75Ah Weight: 48.8 lb (22.1 kg) each
Battery Charger	8 amp, off-board
Warranty	Frame: Lifetime limited Drivetrain: 2-year limited Electronics: 2-year limited Battery: 13 months

1. Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
2. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of +10%, -5%.
3. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
4. Tested in accordance with ANSI/RESNA, WC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
5. Battery weight may vary based on battery manufacturer.

NOW AVAILABLE IN
4 CAMO OPTIONS



**PRIDE
MOBILITY®**

pridemobility.ca



NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

CA_PX4_Trifold_030526



MORE ROOM TO ROAM



FDA Class II Medical Device¹




**WEIGHT
CAPACITY**



**CRUISING
SPEED**



SUSPENSION



**23.4 MILE
RANGE**

LIVE YOUR BEST®

pridemobility.ca



MORE ROOM TO BREATHE

Memory foam seat with 6" slider for adjustable support and premium, all-day comfort

Delta steering wheel with water-resistant LCD console

NFC-enabled key fob for enhanced security



USB charger port

Standard front basket and rear view mirrors

Enhanced LED headlight

Optional storage pod can be mounted to the rear of the seat

Pneumatic tires mounted on split rims for off-road thrills

FDA Class II Medical Device¹


BIG PLANS REQUIRE THE BEST RIDE.

Durable and reliable, this mobility scooter offers supreme comfort thanks to full suspension and memory foam seating in several sizes. The innovative PX4 provides room to breathe – no matter your destination. Don't be afraid to think BIG!

FEATURES

- NFC key fob for enhanced security. Simply place or tap near the receiver to start
- 500 lb. weight capacity
- Impressive top speed of up to 9.3 mph
- Cruises up to 23.4 miles on a single battery charge
- Depth adjustable 18"x18-20" memory foam seat with 6" slider with additional seat size options that include 20"x18-20", 22"x20-22", and 32"x18-20".
- Delta steering wheel with LCD console
- Dual rear view mirrors and front basket (standard)
- Water-resistant user interface
- Full LED lighting package for maximum visibility
- Pneumatic tires mounted on split rims for excellent handling and traction on every terrain
- Supplemental hand brake for additional stopping power
- Choice of nine unique vinyl accent color options, including 4 camo options, with a 7-year outdoor durability rating

ACCENT COLORS



ACCESSORIES**

- Wishbone Crutch Holder
- 3" Single/Double Cane or Crutch Holder
- Cane Crutch Holder
- Cell Phone Holder
- Cup Holder*
- Forearm Crutch Holder
- Lap Belt
- Military Patches
- Oxygen Tank Holder
- Rear Basket*
- Saddlebag*
- Safety Flag
- Storage Pod
- Walker Holder
- Weather Cover

* Cup holder, rear basket and saddlebag can be ordered as a bundle.

**Only one rear accessory can be used at a time.

LCD digital display

Enhanced LED headlight

Full suspension