

# VivaLIFT!®

POWER RECLINER COLLECTION¹

THE MOST  
COMFORTABLE  
RECLINER  
IN THE WORLD!

## THE 4-MOTOR RECLINER



New! Ergonomically designed toggle remote with two programmable memory positions



Shown in Canyon Walnut Fabric

Three different heat pads located in shoulders, lumbar and seat



## RADIANCE¹

Whether you cozy up with a beverage and your favorite book or kick back for a relaxing evening of television, this power recliner is right for you! Experience the warmth of the **Radiance Collection** from VivaLift!® Power Recliners. With three different heating pads in the shoulders, lumbar and seat, you can adjust your heat level accordingly and soothe your tired muscles. Keep your smart phone close at hand while you charge it with the wireless charging station. And if the power goes out, know you have peace of mind with the lithium battery backup.



Scan here for  
more information  
on **VivaLift!®**



**PRIDE**®  
MOBILITY  
LIVE YOUR BEST®

# RADIANCE<sup>1</sup>

**VivaLIFT!**<sup>®</sup>  
POWER RECLINER COLLECTION<sup>1</sup>

COLLECTION	RADIANCE <sup>1</sup>	RADIANCE <sup>1</sup>	RADIANCE <sup>1</sup>
SIZE	Small	Medium	Large/Tall
MODEL	PLR-3955S	PLR-3955M	PLR-3955LT
POSITION	Infinite-Position	Infinite-Position	Infinite-Position
HAND CONTROL	Toggle Remote	Toggle Remote	Toggle Remote
BACK TYPE	Deluxe Biscuit	Deluxe Biscuit	Deluxe Biscuit
SEAT TYPE	Total Comfort	Total Comfort	Total Comfort
STANDARD FEATURES	Yes	Yes	Yes
ADDITIONAL FEATURES			
USER HEIGHT RANGE <sup>2</sup>	5'4" and below	5'4"-5'10"	5'10" and above
SEAT-TO-FLOOR	18"	19.5"	22.5"
SEAT WIDTH	20"	22.5"	23.5"
SEAT DEPTH	18"	21"	22"
SEAT DEPTH WITH FOOTREST ELEVATED	40.5"	38"	39"
TOP OF BACK TO SEAT	25"	28"	31"
CHAIR WIDTH	37"	39"	40"
CHAIR DEPTH	41"	42"	40"
CHAIR HEIGHT	41"	43"	49"
DISTANCE FROM WALL	15"	15"	21.5"
CHAIR WEIGHT	156LB/70.8KG	158LB/71.7KG	178LB/80.7KG
WEIGHT CAPACITY	400LB/158.1KG	400LB/158.1KG	400LB/158.1KG
WARRANTY	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated

## STANDARD FEATURES



- The power headrest lifts your head, neck and shoulders for maximum support.



- The upward lift offers convenient stand assist for getting in and out of the recliner.



- The power lumbar enables you to personalize your comfort.



- The footrest extension allows you to relax and stretch out completely.



- The backrest moves up and down independently to help you find the perfect position.



- Infinite positioning enables the backrest and footrest to move independently so that a variety of positions can be reached including raising your feet above your heart.



- The ergonomically designed Toggle Remote has two programmable memory position buttons, performs the up/down, recline functions, and has a convenient USB charger port that allows you to charge your smart device.



- Lithium battery backup ensures your recliner is operational when the power goes out.

## ADDITIONAL FEATURES

- Three heating pads located in the shoulder, lumbar and seat allow you to adjust your heat levels for your comfort
- Hideaway wireless charging<sup>3</sup> station conveniently charges your phone
- Hideaway cup holder keeps your beverage close at hand


- VivaLift!® Power Recliners are FDA Class II Medical Devices designed to aid individuals with mobility impairments.
- User height ranges should be used as a guide only due to varying factors of each individual user.
- The wireless charger is only compatible with smartphones which have wireless charging capabilities. Additionally, no other smart devices, (i.e., smartwatches) should be used with the wireless charger. The charging pad may not be suitable for all phone types. If the charging pad is not maintaining a suitable charge level for your device, the standard phone charger should be used.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.



= Indicates shortest available lead time

## STANDARD FABRIC: CANYON



OCEAN

SILT

STEEL

WALNUT



**PRIDE**<sup>®</sup>  
MOBILITY  
LIVE YOUR BEST<sup>®</sup>

pridemobility.ca