

Genesis Collection¹




Evolution Collection¹




MODEL	LC-150S¹	LC-150M¹	LC-435S¹	LC-435M¹	LC-435L¹	LC-435XL¹	LC-435XXL¹
Position	Traditional 3-Position	Traditional 3-Position	Enhanced 3-Position	Enhanced 3-Position	Enhanced 3-Position	Traditional 3-Position	Traditional 3-Position
Size	Small	Medium	Small	Medium	Large	Heavy-Duty	Very Heavy-Duty
Suggested Height Range³	5’3” & Below	5’4” - 5’9”	5’3” & Below	5’4” - 5’9”	5’9” - 6’	5’8” - 6’1”	5’8” - 6’1”
Weight Capacity	325 lbs.	325 lbs.	400 lbs.	400 lbs.	400 lbs.	500 lbs.	600 lbs.
Back Style	Split Back	Split Back	Split Back	Split Back	Split Back	Split Back	Split Back
Seating Type	Sinuous Spring	Sinuous Spring	Coil Pocketed Spring	Coil Pocketed Spring	Coil Pocketed Spring	Total Comfort²	Total Comfort²
Footrest Extension	N/A	N/A	Standard	Standard	Standard	N/A	N/A
Fully Padded Chaise	Yes	Yes	Yes	Yes	Yes	Split Ottoman	Split Ottoman
Additional Features	No	No	Cup Holder, Storage Bin, USB C&A Charger Connections			Lithium Battery Backup	Lithium Battery Backup
Seat to Floor Height	18”	20”	18”	20”	22”	20”	20”
Seat Width	20”	22”	20”	22”	24”	25.5”	29.5”
Seat Depth	18”	20”	18”	20”	22”	22”	22”
SeatDepth(withfootrestextended)	37”	38”	39”	41”	43.5”	38”	38”
Top of Back to Seat	24”	26”	24”	26”	28”	27”	27”
Overall Width	32”	33”	38.5”	40.5”	42”	41.5”	46.5”
Overall Depth	35.5”	36”	36”	39”	43.5”	40”	40”
Overall Height	41.5”	44.5”	41”	46”	48.5”	46.5”	46”
Distance from Wall	17”	18.5”	16”	17.5”	19”	19”	19”
Chair Weight	111 lbs.	116 lbs.	136 lbs.	141 lbs.	146 lbs.	187 lbs.	196 lbs.
Warranty	2-Yr. Limited	2-Yr. Limited	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated
FABRIC SELECTION GUIDE	Standard	Sophie   LATTE  SAPPHIRE  SHALE  CABERNET		Cosi   MOCHA  GRAVEL  BLUE DIAMOND  SYRAH		Cosi   MOCHA  GRAVEL	
	Optional	Lexis Vinyl®   CHESTNUT		Sta-Kleen Vinyl®   CHESTNUT		Sta-Kleen Vinyl®   CHESTNUT	

PARTNER TRADEMARKS



ICON KEY

 = Indicates shortest available lead times

NOTE: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

CA_SittingPretty_E2E_072325

FOOTNOTES:

1. Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.

2. Unique spring design and a one-piece soft foam seat for outstanding relaxation.

3. Height ranges should be used only as a guide due to varying factors of each individual user.

Like all furniture, Pride® Power Lift Recliners are individually built and there may be differences in actual measurements due to variations in fabric, upholstery, and filling, in addition to the method of measurement.

ADDITIONAL NOTES:

• All chairs feature premium Furniture-Grade Laminate/Hardwood frames.

• All chairs have a Lifetime Warranty on lift mechanism, shell and recline mechanism (except LC-435XXL).






• All chairs have integrated emergency battery backup (batteries not included).

• All chairs can have fiber fill added or removed for custom comfort adjustment.

• All chairs are shipped in an impact-resistant, reinforced, corrugated box.

POSITION GUIDE

PRODUCTVIEW

 LIFTED  SEATED  READING  NAPPING  RELAXING

Traditional 3-Position

The Traditional 3-Position full recline is available on most Pride® Power Lift Recliner models and puts you in a position ideal for napping or relaxing.

Enhanced 3-Position

The Enhanced 3-Position experience allows the leg rest to extend before the back reclines, offering additional resting positions.

The lifted, seated and reading positions are achievable on all Pride® Power Lift Recliners.



New Sitting Pretty power lift recliners from Pride Mobility® exude style and sophistication. These three-position recliners offer superior craftsmanship and functionality. Eliminate pressure points with a chaise pad design that fully supports the legs and choose from a variety of luxurious fabrics and colors that fit your home décor. Experience at-home comfort like never before!

FEATURES

- Two collections (Genesis LC-150 and Evolution LC-435), each with their own unique features and styles
- Genesis: Small, Medium sizes
- Evolution: Small, Medium, Large, XL, XXL sizes
- Easy-to-use toggle hand control
- Magnetic hand control dock on both arms (featured on all sizes except XL and XXL)
- Three-position lift and recline
- Evolution Small, Medium, and Large sizes offer an enhanced three-position experience with additional resting positions
- Chaise pad offers complete support of the legs to eliminate pressure points
- Battery back-up in case of power failure

OPTIONS & ACCESSORIES

- Lithium battery backup (standard on XL & XXL)
- Head and arm cover kit
- Non-skid leg levelers
- Cleaning kit