

Genesis Collection¹



Evolution Collection¹



MODEL	LC-150S¹	LC-150M¹	LC-435S¹,⁴	LC-435M¹,⁴	LC-435L¹,⁴	LC-435XL¹,⁴	LC-435XXL¹,⁴	
Position	Traditional 3-Position	Traditional 3-Position	Enhanced 3-Position	Enhanced 3-Position	Enhanced 3-Position	Traditional 3-Position	Traditional 3-Position	
Size	Small	Medium	Small	Medium	Large	Heavy-Duty	Very Heavy-Duty	
Suggested Height Range³	5'3" & Below	5'4" - 5'9"	5'3" & Below	5'4" - 5'9"	5'9" - 6'	5'8" - 6'1"	5'8" - 6'1"	
Weight Capacity	325 lbs.	325 lbs.	400 lbs.	400 lbs.	400 lbs.	500 lbs.	600 lbs.	
Back Style	Split Back	Split Back	Split Back	Split Back	Split Back	Split Back	Split Back	
Seating Type	Sinuous Spring	Sinuous Spring	Coil Pocketed Spring	Coil Pocketed Spring	Coil Pocketed Spring	Total Comfort²	Total Comfort²	
Footrest Extension	N/A	N/A	Standard	Standard	Standard	N/A	N/A	
Fully Padded Chaise	Yes	Yes	Yes	Yes	Yes	Split Ottoman	Split Ottoman	
Additional Features	No	No	Cup Holder, Storage Bin, USB C&A Charger Connections			Lithium Battery Backup	Lithium Battery Backup	
Seat to Floor Height	18"	20"	18"	20"	22"	20"	20"	
Seat Width	20"	22"	20"	22"	24"	25.5"	29.5"	
Seat Depth	18"	20"	18"	20"	22"	22"	22"	
SeatDepth(withfootrestextended)	37"	38"	39"	41"	43.5"	38"	38"	
Top of Back to Seat	24"	26"	24"	26"	28"	27"	27"	
Overall Width	32"	33"	38.5"	40.5"	42"	41.5"	46.5"	
Overall Depth	35.5"	36"	36"	39"	43.5"	40"	40"	
Overall Height	41.5"	44.5"	41"	46"	48.5"	46.5"	46"	
Distance from Wall	17"	18.5"	16"	17.5"	19"	19"	19"	
Chair Weight	111 lbs.	116 lbs.	136 lbs.	141 lbs.	146 lbs.	187 lbs.	196 lbs.	
Warranty	2-Yr. Limited	2-Yr. Limited	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	
FABRIC SELECTION GUIDE	Standard	Sophie <sup>QS</sup> LATTE SAPPHIRE SHALE CABERNET		Cosi <sup>QS</sup> MOCHA GRAVEL BLUE DIAMOND SYRAH			Cosi <sup>QS</sup> MOCHA GRAVEL	
	Optional	Lexis Vinyl <sup>QS</sup> CHESTNUT		Sta-Kleen <sup>QS</sup> CHESTNUT			Sta-Kleen <sup>QS</sup> CHESTNUT	

PARTNER TRADEMARKS

ICON KEY

<sup>QS</sup> = Indicates shortest available lead times

FOOTNOTES:

- Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
- Unique spring design and a one-piece soft foam seat for outstanding relaxation.
- Height ranges should be used only as a guide due to varying factors of each individual user.
- Evolution Collection chairs have a Lifetime Warranty on lift mechanism, shell and recline mechanism (except LC-435XXL).

Like all furniture, Pride® Power Lift Recliners are individually built and there may be differences in

ADDITIONAL NOTES:

- All chairs feature premium Furniture-Grade Laminate/Hardwood frames.
- All chairs have integrated emergency battery backup (batteries not included).
- All chairs can have fiber fill added or removed for custom comfort adjustment.
- All chairs are shipped in an impact-resistant, reinforced, corrugated box.

POSITION GUIDE

PRODUCT VIEW	LIFTED	SEATED	READING	NAPPING	RELAXING
Traditional 3-Position	The Traditional 3-Position full recline is available on most Pride® Power Lift Recliner models and puts you in a position ideal for napping or relaxing.				
Enhanced 3-Position	The Enhanced 3-Position experience allows the leg rest to extend before the back reclines, offering additional resting positions.				

The lifted, seated and reading positions are achievable on all Pride® Power Lift Recliners.

NOTE: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.



GENESIS  
LC-150

EVOLUTION  
LC-435

Sitting Pretty™ power lift recliners are FDA Class II Medical Devices designed to aid individuals with mobility impairments.



Storage  
Compartment  
Only available on  
LC-435 S/M/L



Cup  
Holder  
Only available on  
LC-435 S/M/L



New Sitting Pretty power lift recliners from Pride Mobility® exude style and sophistication. These three-position recliners offer superior craftsmanship and functionality. Eliminate pressure points with a chaise pad design that fully supports the legs and choose from a variety of luxurious fabrics and colors that fit your home décor. Experience at-home comfort like never before!

#### FEATURES

- Two collections (Genesis LC-150 and Evolution LC-435), each with their own unique features and styles
- Genesis: Small, Medium sizes
- Evolution: Small, Medium, Large, XL, XXL sizes
- Easy-to-use toggle hand control
- Magnetic hand control dock on both arms (featured on all sizes except XL and XXL)
- Three-position lift and recline
- Evolution Small, Medium, and Large sizes offer an enhanced three-position experience with additional resting positions
- Chaise pad offers complete support of the legs to eliminate pressure points
- Battery back-up in case of power failure

#### OPTIONS & ACCESSORIES

- Lithium battery backup (standard on XL & XXL)
- Head and arm cover kit
- Non-skid leg levelers
- Cleaning kit