

SUPERB COMFORT

# EXCEPTIONAL VALUE



GENESIS  
LC-150



EVOLUTION  
LC-435

Sitting Pretty™ power lift recliners are FDA Class II Medical Devices designed to aid individuals with mobility impairments.



Storage  
Compartment  
Only available on  
LC-435 S/M/L



Cup  
Holder  
Only available on  
LC-435 S/M/L

## Sitting Pretty™

New Sitting Pretty power lift recliners from Pride Mobility® exude style and sophistication. These three-position recliners offer superior craftsmanship and functionality. Eliminate pressure points with a chaise pad design that fully supports the legs and choose from a variety of luxurious fabrics and colors that fit your home décor. Experience at-home comfort like never before!

PRIDEMOBILITY.CA

### FEATURES

- Two collections (Genesis LC-150 and Evolution LC-435), each with their own unique features and styles
- Genesis: Small, Medium sizes
- Evolution: Small, Medium, Large, XL, XXL sizes
- Easy-to-use toggle hand control
- Magnetic hand control dock on both arms (featured on all sizes except XL and XXL)
- Three-position lift and recline
- Evolution Small, Medium, and Large sizes offer an enhanced three-position experience with additional resting positions
- Chaise pad offers complete support of the legs to eliminate pressure points
- Battery back-up in case of power failure

### OPTIONS & ACCESSORIES

- Lithium battery backup (standard on XL & XXL)
- Head and arm cover kit
- Non-skid leg levelers
- Cleaning kit

Sitting Pretty™ power lift recliners are FDA Class II Medical Devices designed to aid individuals with mobility impairments.



	GENESIS		EVOLUTION				
Size	Small	Medium	Small	Medium	Large	X-Large	XX-Large
Model	LC-150S	LC-150M	LC-435S	LC-435M	LC-435L	LC-435XL	LC-435XXL
Position Type	Traditional 3-position	Traditional 3-position	Enhanced 3-position	Enhanced 3-position	Enhanced 3-position	Traditional 3-position	Traditional 3-position
Hand Control	Toggle	Toggle	Toggle	Toggle	Toggle	Toggle	Toggle
Back Type	Blown Fiber	Blown Fiber	Blown Fiber	Blown Fiber	Blown Fiber	Blown Fiber	Blown Fiber
Seating Type	Sinous Spring	Sinous Spring	Coil Pocketed Spring	Coil Pocketed Spring	Coil Pocketed Spring	Total Comfort	Total Comfort
Footrest Extension	Standard	Standard	Standard	Standard	Standard	Standard	Standard
Additional Features	N/A	N/A	Cup Holder, Storage Bin, USB C&A Charger Connections			Lithium Battery Back-up	
User Height Range	5'3" & under	5'4" - 5'9"	5'3" & under	5'4" - 5'9"	5'9" - 6'	5'8" - 6'1"	5'8" - 6'1"
Seat-to-Floor	18"	20"	18"	20"	22"	20"	20"
Seat Width	20"	22"	20"	22"	24"	25.5"	29.5"
Seat Depth	18"	20"	18"	20"	22"	22"	22"
Seat Depth with footrest extended	37"	38"	39"	41"	43.5"	38"	38"
Top of Back to Seat	24"	26"	24"	26"	28"	27"	27"
Chair Width	32"	33"	38.5"	40.5"	42"	41.5"	46.5"
Chair Depth	35.5"	36"	36"	39"	43.5"	40"	40"
Chair Height	41.5"	44.5"	41"	46"	48.5"	46.5"	46"
Distance from Wall	17"	18.5"	18"	17.5"	19"	19"	19"
Chair Weight	111 lbs.	116 lbs.	136 lbs.	141 lbs.	146 lbs.	187 lbs.	196 lbs.
Weight Capacity	325 lbs.	325 lbs.	400 lbs.	400 lbs.	400 lbs.	500 lbs.	600 lbs.
Warranty	2-year Limited	2-year Limited	7-year Pro-rated	7-year Pro-rated	7-year Pro-rated	7-year Pro-rated	7-year Pro-rated


Fabric Options

LC-150S / LC-150M

Sophie Fabrics



Latte



Sapphire




Shale



Cabernet


Lexis Vinyl<sup>1</sup>



Chestnut

LC-435S / LC-435M / LC-435L


Cosi Fabrics




Mocha



Gravel




Blue Diamond



Syrah


Sta-Kleen<sup>1</sup>



Chestnut

LC-435XL / LC-435XXL

Cosi Fabrics




Mocha



Gravel

Sta-Kleen<sup>1</sup>



Chestnut



pridemobility.ca



NOTE: Measurements may vary +/- 1". The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

1) Additional charges will apply for this fabric.