

End-to-End

**Genesis Collection¹**



**Evolution Collection¹**



MODEL	LC-150S¹	LC-150M¹	NEW LC-435PW¹,⁴	LC-435S¹,⁴	LC-435M¹,⁴	LC-435L¹,⁴	LC-435XL¹,⁴	LC-435XXL¹,⁴
Position	Traditional 3-Position		Enhanced 3-Position				Traditional 3-Position	
Size	Small	Medium	Petite Wide	Small	Medium	Large	Heavy-Duty	Very Heavy-Duty
Suggested Height Range³	5'3" & Below	5'4" - 5'9"	5'3" & Below	5'3" & Below	5'4" - 5'9"	5'9" - 6'	5'8" - 6'1"	5'8" - 6'1"
Weight Capacity	325 lb	325 lb	400 lb	400 lb	400 lb	400 lb	500 lb	600 lb
Back Style	Split Back	Split Back	Split Back	Split Back	Split Back	Split Back	Split Back	Split Back
Seating Type	Sinuous Spring	Sinuous Spring	Coil Pocketed Spring	Coil Pocketed Spring	Coil Pocketed Spring	Coil Pocketed Spring	Total Comfort²	Total Comfort²
Footrest Extension	N/A	N/A	Standard	Standard	Standard	Standard	N/A	N/A
Fully Padded Chaise	Yes	Yes	Yes	Yes	Yes	Yes	Split Ottoman	Split Ottoman
Additional Features	None		Cup Holder, Storage Bin, USB C&A Charger Connections				Lithium Battery Backup	
Seat to Floor Height	18"	20"	20"	18"	20"	22"	20"	20"
Seat Width	20"	22"	24"	20"	22"	24"	25.5"	29.5"
Seat Depth	18"	20"	19"	18"	20"	22"	22"	22"
Seat Depth (with footrest extended)	37"	38"	38"	39"	41"	43.5"	38"	38"
Top of Back to Seat	24"	26"	23"	24"	26"	28"	27"	27"
Overall Width	32"	33"	41"	38.5"	40.5"	42"	41.5"	46.5"
Overall Depth	35.5"	36"	36"	36"	39"	43.5"	40"	40"
Overall Height	41.5"	44.5"	40"	41"	46"	48.5"	46.5"	46"
Distance from Wall	17"	18.5"	17"	16"	17.5"	19"	19"	19"
Chair Weight	111 lb	116 lb	130 lb	136 lb	141 lb	146 lb	187 lb	196 lb
Warranty	2-Yr. Limited	2-Yr. Limited	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated

FABRIC SELECTION GUIDE	STANDARD FABRICS	<p>Sophie </p> <p>LATTE SAPPHIRE SHALE CABERNET</p>	<p>Cosi </p> <p>MOCHA GRAVEL BLUE DIAMOND SYRAH</p>	<p>Cosi </p> <p>MOCHA GRAVEL</p>
	OPTIONAL FABRIC	<p>Lexis Vinyl® </p> <p>CHESTNUT</p>	<p>Sta-Kleen® </p> <p>CHESTNUT</p>	<p>Sta-Kleen® </p> <p>CHESTNUT</p>

PARTNER TRADEMARKS



ICON KEY

= Indicates shortest available lead times

FOOTNOTES:

1. Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
2. Unique spring design and a one-piece soft foam seat for outstanding relaxation.
3. Height ranges should be used only as a guide due to varying factors of each individual user.
4. Evolution Collection chairs have a Lifetime Warranty on lift mechanism, shell and recline mechanism (except LC-435XXL).

Like all furniture, Pride® Power Lift Recliners are individually built and there may be differences in actual measurements due to variations in fabric, upholstery, and filling, in addition to the method of measurement.

ADDITIONAL NOTES:

- All chairs feature premium Furniture-Grade Laminate/Hardwood frames.
- All chairs have integrated emergency battery backup (batteries not included).
- All chairs can have fiber fill added or removed for custom comfort adjustment.
- All chairs are shipped in an impact-resistant, reinforced, corrugated box.

POSITION GUIDE

PRODUCT VIEW	LIFTED	SEATED	READING	NAPPING	RELAXING
Traditional 3-Position	The Traditional 3-Position full recline is available on most Pride® Power Lift Recliner models and puts you in a position ideal for napping or relaxing.				
Enhanced 3-Position	The Enhanced 3-Position experience allows the leg rest to extend before the back reclines, offering additional resting positions.				
The lifted, seated and reading positions are achievable on all Pride® Power Lift Recliners.					