

# Genesis Collection<sup>1</sup>









Evolution Collection<sup>1</sup>





							i .	
	MODEL	LC-150S1	LC-150M <sup>1</sup>	LC-435S <sup>1,4</sup>	LC-435M <sup>1,4</sup>	LC-435L <sup>1,4</sup>	LC-435XL <sup>1,4</sup>	LC-435XXL <sup>1,,4</sup>
	Position	Traditional 3-Position	Traditional 3-Position	Enhanced 3-Position	Enhanced 3-Position	Enhanced 3-Position	Traditional 3-Position	Traditional 3-Position
	Size	Small	Medium	Small	Medium	Large	Heavy-Duty	Very Heavy-Duty
Suggested Height Range <sup>3</sup>		5'3" & Below	5′4″ - 5′9″	5′3″ & Below	5′4″ - 5′9″	5′9″ - 6′	5′8″ - 6′1″	5′8″ - 6′1″
Weight Capacity  Back Style  Seating Type		325 lbs.	325 lbs.	400 lbs.	400 lbs.	400 lbs.	500 lbs.	600 lbs.
		Split Back	Split Back	Split Back	Split Back	Split Back	Split Back	Split Back
		Sinuous Spring	Sinuous Spring	Coil Pocketed Spring	Coil Pocketed Spring	Coil Pocketed Spring	Total Comfort <sup>2</sup>	Total Comfort <sup>2</sup>
Footre	est Extension	N/A	N/A	Standard	Standard	Standard	N/A	N/A
Fully Pac	dded Chaise	Yes	Yes	Yes	Yes	Yes	Split Ottoman	Split Ottoman
Additional Features		No	No	Cup Holder, Storage Bin, USB C&A Charger Connections			Lithium Battery Backup	Lithium Battery Backup
Seat to F	Floor Height	18"	20"	18"	20"	22"	20"	20"
	Seat Width	20"	22"	20"	22"	24"	25.5"	29.5″
	Seat Depth	18"	20"	18"	20"	22"	22"	22"
SeatDepth(withfoo	otrestextended)	37"	38"	39"	41"	43.5"	38"	38"
Top of	f Back to Seat	24"	26"	24"	26"	28"	27"	27"
Ov	verall Width	32"	33"	38.5"	40.5"	42"	41.5"	46.5"
Ov	verall Depth	35.5"	36"	36"	39"	43.5"	40"	40"
Ov	verall Height	41.5"	44.5"	41"	46"	48.5"	46.5"	46"
Distan	nce from Wall	17"	18.5″	16"	17.5″	19"	19"	19"
C	Chair Weight	111 lbs.	116 lbs.	136 lbs.	141 lbs.	146 lbs.	187 lbs.	196 lbs.
	Warranty	2-Yr. Limited	2-Yr. Limited	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated
出	Standard	Sophie 😎		Cosi as			Cosi OS	
SELECTION GUIDE		LATTE SAPPHIRE SHALE CABERNET		MOCHA GRAVEL BLUE SYRAH DIAMOND			MOCHA GRAVEL	
FABRIC SELEC	Optional	Lexis Vinyl* CHESTNUT		Sta-Kleen Vinyl  CHESTNUT			Sta-Kleen Vinyl  CHESTNUT	

PARTNER TRADEMARKS

sta\*kleen.

ICON KEY = Indicates shortest available lead times

 $NOTE: Measurements\ may\ vary\ +/-1.''Specifications\ subject to\ change\ without\ prior\ notice. The\ information\ contained\ herein\ is\ correct\ at\ the\ time\ of\ publication;\ we\ reserve\ the\ right\ to\ alter\ specifications\ without\ prior\ notice.$ 

- $1.\ Pride^*FDA\ Class\ II\ Medical\ Devices\ are\ designed\ to\ aid\ individuals\ with\ mobility\ impairments.$
- 2. Unique spring design and a one-piece soft foam seat for outstanding relaxation. 3. Height ranges should be used only as a guide due to varying factors of each individual user.
- 4. Evolution Collection chairs have a Lifetime Warranty on lift mechanism, shell and recline mechanism (except LC-435XXL).

Like all furniture, Pride® Power Lift Recliners are individually built and there may be differences in

- All chairs feature premium Furniture-Grade Laminate/Hardwood frames.
- $\bullet \ \ \text{All chairs have integrated emergency battery backup (batteries not included)}.$
- All chairs can have fiber fill added or removed for custom comfort adjustment. • All chairs are shipped in an impact-resistant, reinforced, corrugated box.

POSITION GUIDE										
PRODUCTVIEW	LIFTED	SEATED	READING	NAPPING	RELAXING					
Traditional 3-Position	and the second control of the second control									

The lifted, seated and reading positions are achievable on all Pride® Power Lift Recliners.





New Sitting Pretty power lift recliners from Pride Mobility® exude style and sophistication. These three-position recliners offer superior craftsmanship and functionality. Eliminate pressure points with a chaise pad design that fully supports the legs and choose from a variety of luxurious fabrics and colors that fit your home décor. Experience at-home comfort like never before!

### **FEATURES**

- Two collections (Genesis LC-150 and Evolution LC-435), each with their own unique features and styles
- Genesis: Small, Medium sizes
- Evolution: Small, Medium, Large, XL, XXL sizes
- Easy-to-use toggle hand control
- Magnetic hand control dock on both arms (featured on all sizes except XL and XXL)
- Three-position lift and recline
- Evolution Small, Medium, and Large sizes offer an enhanced three-position experience with additional resting positions
- Chaise pad offers complete support of the legs to eliminate pressure points
- Battery back-up in case of power failure

## **OPTIONS & ACCESSORIES**

- Lithium battery backup (standard on XL & XXL)
- Head and arm cover kit
- Non-skid leg levelers
- Cleaning kit