

VivaLIFT!®

POWER RECLINER COLLECTION¹

THE MOST COMFORTABLE RECLINER IN THE WORLD!



NEW Exclusive Weightless Tilt



Easily accessible heat and massage buttons on the inside of the left armrest



VIVALIFT!® ULTRA¹

When you need to soothe away the stress of your day, the **VivaLift! Ultra** Collection by VivaLift!® Power Recliners is there for you! Exclusive Weightless Tilt technology reclines so your knees are elevated above your heart. The **advanced massage** system delivers gentle pressure on your lower back. Choose from three intensity settings for the perfect massage!



Scan here for more information on **VivaLift!®**



PRIDE®
MOBILITY
LIVE YOUR BEST®

VIVALIFT!® ULTRA¹

VivaLIFT!
POWER RECLINER COLLECTION¹

COLLECTION	ULTRA ¹	ULTRA ¹	ULTRA ¹
SIZE	Small	Medium	Large
MODEL	PLR-4955S	PLR-4955M	PLR-4955L
POSITION	Infinite-Position, Weightless Tilt	Infinite-Position, Weightless Tilt	Infinite-Position, Weightless Tilt
HAND CONTROL	Toggle Remote	Toggle Remote	Toggle Remote
BACK TYPE	Deluxe Biscuit	Deluxe Biscuit	Deluxe Biscuit
SEAT TYPE	Total Comfort	Total Comfort	Total Comfort
STANDARD FEATURES	Yes	Yes	Yes
ADDITIONAL FEATURES	Heat, Massage, Cup Holder, Wireless Charger ³		
USER HEIGHT RANGE ²	5'4" and below	5'4"-5'10"	5'10" and above
SEAT-TO-FLOOR	18"	19.5"	22.5"
SEAT WIDTH	20"	22.5"	23.5"
SEAT DEPTH	18"	21"	22"
SEAT DEPTH WITH FOOTREST ELEVATED	40.5"	38"	39"
TOP OF BACK TO SEAT	25"	28"	30"
CHAIR WIDTH	37"	39"	40"
CHAIR DEPTH	41"	43"	49"
CHAIR HEIGHT	41"	40"	42"
DISTANCE FROM WALL	15"	15"	21.5"
CHAIR WEIGHT	158LB/71.7KG	158LB/71.7KG	176LB/79.8KG
WEIGHT CAPACITY	400LB/158.1KG	400LB/158.1KG	400LB/158.1KG
WARRANTY	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated

STANDARD FEATURES

- Infinite positioning enables the backrest and footrest to move independently from one another so that a variety of positions can be reached
- The power headrest lifts your head, neck, and shoulders for maximum support
- The upward lift offers convenient stand assist for getting in and out of the recliner
- The power lumbar enables you to personalize your comfort
- The footrest extension allows you to relax and stretch out completely
- Lithium battery backup ensures your recliner is operational when the power goes out
- **Newly** designed Toggle Remote:
 - Modern toggle design for effortless control
 - Easy-to-use toggle in the middle of the remote performs the up / down and recline positions
 - Memory button to program your favorite position
 - Tilt buttons to enhance seat angle
 - USB charger

ADDITIONAL FEATURES

- **Weightless Tilt**, tilts the lift chair into a position where one's knees are elevated above your heart, providing the feeling of weightlessness
- **Advanced massage** system with a touch of a button provides gentle pressure to your lower back for a soothing massage with three intensity settings
- Three heating pads located in the shoulder, lumbar and seat allow you to adjust your heat levels for your comfort
- Easily accessible heat and massage buttons on the inside of the left armrest
- Hideaway wireless charging³ station conveniently charges your phone
- Hideaway cup holder keeps your beverage close at hand

1. VivaLift!® Power Recliners are FDA Class II Medical Devices designed to aid individuals with mobility impairments.
2. User height ranges should be used as a guide only due to varying factors of each individual user.
3. The wireless charger is only compatible with smartphones which have wireless charging capabilities. Additionally, no other smart devices, (i.e., smartwatches) should be used with the wireless charger. The charging pad may not be suitable for all phone types. If the charging pad is not maintaining a suitable charge level for your device, the standard phone charger should be used.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

 = Indicates shortest available lead time

STANDARD FABRIC: **CAPRICCIO** 



SLATE CAPPUCCINO SMOKE DOVE



PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

pridemobility.ca