

VivaLIFT!®

POWER RECLINER COLLECTION¹

THE MOST COMFORTABLE RECLINER IN THE WORLD!

Power headrest

Power backrest

Power lumbar

NEW Advanced massage system

Ergonomically designed toggle remote with one programmable memory position

Lithium battery backup

Shown in Sorrento Coffee 100% Leather

Heating pads in the shoulders, lumbar and seat

Heat and massage buttons

Cup holder

Wireless charging station

Footrest extension

Exclusive VivaLift!® Tilt

Ergonomically designed toggle remote with programmable memory position button and USB port for charging devices



Exclusive VivaLift!® Tilt



Easily accessible heat and massage buttons on the inside of the left armrest

FDA Class II Medical Device¹



VIVALIFT!® ULTRA¹

The **Ultra** collection by VivaLift!® Power Recliners is loaded with features to help you experience pure relaxation. VivaLift! Tilt technology helps you recline with your knees elevated about your heart. The **advanced massage** system delivers gentle pressure on your lower back, and three heating pads soothe away the stress of the day.



Scan here for more information on **VivaLift!®**



PRIDE®
MOBILITY
LIVE YOUR BEST®

VIVALIFT!® ULTRA¹

VivaLIFT!
POWER RECLINER COLLECTION¹


COLLECTION	ULTRA ¹	ULTRA ¹	ULTRA ¹
MODEL	PLR-4955S	PLR-4955M	PLR-4955L
POSITION	True-Infinite	True-Infinite	True-Infinite
SIZE	Small 5'4" and below	Medium 5'5"-6'	Large/Tall 6' and above
WEIGHT CAPACITY	400 lbs.	400 lbs.	400 lbs.
SEAT-TO-FLOOR	18"	19.5"	21.5"
SEAT DEPTH	18"	21"	22"
SEAT WIDTH	20"	22.5"	23"
TOP OF BACK TO SEAT	25"	28"	30"
FULLY PADDED CHAISE	Yes	Yes	Yes
FOOTREST EXTENSION	Standard	Standard	Standard
SEATING TYPE	Total Comfort	Total Comfort	Total Comfort
CHAIR HEIGHT	41"	43"	47"
DISTANCE FROM WALL	15"	15"	21.5"
OVERALL WIDTH	37"	39"	39"
CHAIR WEIGHT	158 lbs.	158 lbs.	176 lbs.
WARRANTY	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated

FEATURES AND BENEFITS

- New!** VivaLift!® Tilt reclines the lift chair into a recline position where one's knees are elevated above your heart, providing the feeling of weightlessness
- New!** Advanced massage system with a touch of a button provides gentle pressure to your lower back for a soothing massage with three intensity settings
- Newly** designed toggle remote:
 - Easy-to-use toggle in the middle of the remote performs the up /down and recline positions
 - One memory button to program your customer's favorite position
 - Tilt button
 - Modern toggle design for effortless control
 - USB charger
- Easily accessible heat and massage buttons on the inside of the left armrest
- Three heating pads located in the shoulder, lumbar and seat allow you to adjust your heat levels for your comfort
- True-infinite positioning enables the backrest and footrest to move independently from one another so that a variety of positions can be reached
- Hideaway wireless charging station conveniently charges your phone
- Hideaway cup holder keeps your beverage close at hand
- The power headrest lifts your head, neck, and shoulders for maximum support
- The upward lift offers convenient stand assist for getting in and out of the recliner
- The power lumbar enables you to personalize your comfort
- The footrest extension allows you to relax and stretch out completely
- Lithium battery backup ensures your recliner is operational when the power goes out

STANDARD FABRICS:

CAPRICCIO



SLATE

CAPPUCCINO

SMOKE

DOVE

NEW OPTIONAL FABRIC:

SORRENTO*



COFFEE

1. VivaLift!® Power Recliners are FDA Class II Medical Devices designed to aid individuals with mobility impairments.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

*Upgrade charges apply.

PRIDE
MOBILITY
LIVE YOUR BEST®

pridemobility.ca