

A SMOOTH RIDE for Your Active Lifestyle



Victory® 10 LX
shown with solid tires

VICTORY® 10
LX WITH CTS SUSPENSION

CTS SUSPENSION



Front Suspension



Rear Suspension



The **Victory® 10 LX with CTS Suspension** offers a smooth ride with performance to match. Comfort-Trac Suspension (CTS) is the next generation in mobility scooter ride comfort. This industry first technology incorporates advanced suspension components for improved absorption over rough surfaces and a more comfortable ride experience over various terrains.

PRIDE®
MOBILITY
SCOOTERS

| pridemobility.ca |



Specifications 4W: S710LX

Weight Capacity	400 lb (181.4 kg)
Maximum Speed ¹	Up to 5.6 mph (9 km/h)
Ground Clearance ²	2.25" at motor
Turning Radius ²	61.25"
Overall Length ²	47"
Overall Width ²	22.25"
Seat-to-Deck Height ²	16.25"-18.25"
Seat-to-Ground Height ²	23.5"-25.5"
Front Tires ²	3.5" x 10" solid or pneumatic
Rear Tires ²	3.5" x 10" solid or pneumatic
Suspension	Front and rear
Range Per Charge ^{1,4}	Up to 16.2 mi (26 km)
Total Weight ⁶ Without Batteries	156.25 lb (71 kg)
Heaviest Piece When Disassembled ³	78 lb (35.4 kg) (front section)
Standard Seating ²	Type: Reclining high-back Weight: 45.5 lb (20.6 kg) Material: Black Vinyl Dimensions: 18"w x 18"d
Drive System	Rear-wheel drive sealed transaxle with 24-volt DC motor
Dual Braking System	Electronic, regenerative and electromechanical
Battery Requirements ^{3,5}	Size: (2) 40 AH standard Weight: 32 lb (14.5 kg) each
Battery Charger	5-amp, off-board (40 AH)
Colors	  Candy Apple Red Viper Blue

¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).

² Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or 3%).

³ AGM or gel cell type required.

⁴ Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

⁵ Battery weight may vary based on manufacturer.

⁶ Includes standard seat

Features

- Comfort-Trac front and rear suspension
- Exclusive one-hand feather touch disassembly provides the ultimate quick transport system
- Adjustable Delta tiller with wraparound handles' easy-to-read LED battery meter, and a courtesy light under tiller
- Long-lasting, high-intensity LED lighting offers turn signal, 4-way flasher, and headlight for optimal pathway illumination
- Depth-adjustable limited recline, black vinyl, swivel, high back with height adjustable headrest, flip-up height and width adjustable armrests
- Storage compartment below tiller
- Easy-grip tiller adjustment knob
- Easy access tiller-mounted charger port
- Removeable front mounted basket
- Unoccupied transit tie-down points

Accessories

- Cane/crutch holder (*single, double*)
- Cell phone holder
- Cup holder
- Additional rear-view mirror
- Forearm crutch holder
- Oxygen tank holder
- Quad cane holder
- Rear basket
- Saddlebag
- Safety flag
- Walker holder
- Weather cover
- Wishbone crutch holder
- XLR USB charger
- Lap belt

NOTE: All specifications subject to change without notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Dimensions listed are for scooter only. Overall measurements will vary based on seating and accessory selections.