

Model	Victory® 10.2 <sup>1</sup>
Product Number	SC6102 / SC7102
Weight Capacity	400 lb / 181 kg
Maximum Speed <sup>2</sup>	Up to 5.2 mph / 8.4 km/h
Ground Clearance <sup>4</sup>	2.5" at transaxle
Turning Radius <sup>4</sup>	3-wheel 42" / 4-wheel 45.5"
Overall Length <sup>4</sup>	3-wheel 45" / 4-wheel 47.25"
Overall Width <sup>4</sup>	22.25"
Seat-to-Ground Height	22"-24"
Seat-to-Deck Height	15.25"-17.25"
Front Tires	3.5" x 10" solid
Rear Tires	3.5" x 10" solid
Range Per Charge <sup>2,3</sup>	Up to 11 mi / 18 km
Total Weight (without batteries) <sup>6</sup>	3-wheel 130.25 lb / 59 kg 4-wheel 139 lb / 63 kg
Heaviest Piece When Disassembled <sup>4</sup>	56.5 lb (front section)
Standard Seating	Type: Foldable Weight: 28.5 lb / 13 kg Material: Black Vinyl Dimensions: 18" width, 17" depth
Drive System	Rear-wheel drive sealed transaxle with 24-volt DC motor
Dual Braking System	Regenerative and electromechanical
Battery Requirements <sup>5,7</sup>	Size: U-1 (2) standard, Weight: 22.5 lb / 10.25 kg each
Battery Charger	3.5-amp, off-board (U1)
Warranty	Frame: lifetime limited Drivetrain: 2-year limited Electronics: 2-year limited Battery: 13-months

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

**PRIDE**  
MOBILITY



- <sup>1</sup> Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
- <sup>2</sup> Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- <sup>3</sup> Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- <sup>4</sup> AGM or gel cell type required.
- <sup>5</sup> Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- <sup>6</sup> Battery weight may vary based on manufacturer.
- <sup>7</sup> Includes standard seat weight.

CA\_Victory102\_Trifold\_052825

**CONTEMPORARY FEATURES  
ENHANCE YOUR ACTIVE LIFESTYLE**

FDA Class II Medical Device<sup>1</sup>

**LIVE YOUR BEST®**

# CONVENIENCE AND VERSATILITY



18" x 17"  
Foldable molded,  
black vinyl seat  
with 8-point  
swivel and  
sliders for  
increased  
comfort

Integrated  
storage/cup  
holder and USB  
charging port

Standard  
front basket

Easy-to-grip  
tiller adjustment  
knob

Long-lasting,  
maximum-  
intensity LED  
headlight  
offers optimal  
pathway  
illumination

Exclusive  
stylish,  
lightweight,  
non-scuffing,  
black, low-  
profile tires

FDA Class II Medical Device<sup>1</sup>



Frosted LED headlights

Under-tiller puddle light

Ambient voltmeter

## CONVENIENCE AND VERSATILITY

The Victory® 10.2 offers great convenience and versatility with exclusive feather-touch disassembly and a delta tiller with wrap-around handles. Contemporary features such as the integrated storage/cup holder on the tiller ensure ample storage options, while a new under-tiller puddle light offers improved visibility to both the storage area and USB charging port.

### Features

- Exclusive one-hand feather touch disassembly
- Integrated tiller charger port lets you conveniently charge your scooter and mobile devices
- Delta tiller with ergonomic wraparound handles lets you operate the scooter with one hand and rest your wrist
- EZ Turn technology on the 4-wheel scooter allow for the tight turning radius of the 3 wheel
- Integrated storage/cup holder in tiller

### Accessories

- Rearview mirror
- Cell phone holder
- Cup holder
- Saddlebag
- 50", 60", 70" Lap belt
- Safety flag
- Rear basket
- Pneumatic tires
- Weather cover
- Cane clip attachment
- Cane/crutch holder (single, double)
- Forearm crutch holder
- Quad cane holder
- Walker holder
- Wishbone crutch holder
- Oxygen tank holder

### Optional Seating

Foldable molded, black vinyl, flip-up height and width adjustable armrests:

- 20" x 18" (w/sliders)

Limited recline, black vinyl, high back with height adjustable headrest, flip-up height and width adjustable armrests:

- 18" x 18" (w/sliders)
- 22" x 18"

### Couloers:

Candy  
Apple  
Red



Viper  
Blue

