

VivaLIFT!®

POWER RECLINER COLLECTION¹

THE MOST
COMFORTABLE
RECLINER
IN THE WORLD!

THE 4-MOTOR RECLINER



Shown in Astro Brown Fabric

COUNTLESS
POSITION
COMBINATIONS

TRANQUIL 2¹

There is no greater joy than reading a good book in a comfortable, cozy recliner. The **Tranquil 2 Collection** by VivaLift!® Power Recliners combines quiet comfort with gorgeous fabrics, for a rich and upscale look. Get the support you need for your neck, shoulders and lower back as you dive into your favorite story. The Tranquil Collection exudes style and class with contrast stitching. It's the perfect addition to your home.



Scan here for
more information
on **VivaLift!®**











PRIDE®
MOBILITY
LIVE YOUR BEST®

TRANQUIL 2¹

VivaLIFT![®]
POWER RECLINER COLLECTION¹

COLLECTION	TRANQUIL 2 ¹	TRANQUIL 2 ¹	TRANQUIL 2 ¹	TRANQUIL 2 ¹
SIZE	Petite Wide	Small	Medium	Large/Tall
MODEL	PLR-935PW	PLR-935S	PLR-935M	PLR-935LT
POSITION	Infinite-Position	Infinite-Position	Infinite-Position	Infinite-Position
HAND CONTROL	Toggle Remote	Toggle Remote	Toggle Remote	Toggle Remote
BACK TYPE	Waterfall	Waterfall	Waterfall	Waterfall
SEAT TYPE	Total Comfort	Total Comfort	Total Comfort	Total Comfort
STANDARD FEATURES	Yes	Yes	Yes	Yes
USER HEIGHT RANGE ²	5'4" and below	5'4" and below	5'4"-5'10"	5'10" and above
SEAT-TO-FLOOR	18"	18"	19.5"	21.5"
SEAT WIDTH	24"	20"	22"	23.5"
SEAT DEPTH	18.5"	18.5"	20.5"	22.5"
SEAT DEPTH WITH FOOTREST ELEVATED	36"	36"	38"	40"
TOP OF BACK TO SEAT	24.5"	25"	27"	30"
CHAIR WIDTH	36"	33.8"	37.5"	38"
CHAIR DEPTH	39.5"	36.8"	42"	43"
CHAIR HEIGHT	38"	40"	42"	47"
DISTANCE FROM WALL	12.5"	13"	16"	16"
CHAIR WEIGHT	143LB/64.9KG	137LB/62.1KG	145LB/65.8KG	155LB/70.3KG
WEIGHT CAPACITY	400LB/158.1KG	400LB/158.1KG	400LB/158.1KG	400LB/158.1KG
WARRANTY	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated

STANDARD FEATURES


-  The power headrest lifts your head, neck and shoulders for maximum support.
-  The upward lift offers convenient stand assist for getting in and out of the recliner.
-  The power lumbar enables you to personalize your comfort.
-  The footrest extension allows you to relax and stretch out completely.
-  The backrest moves up and down independently to help you find the perfect position.
-  Infinite positioning enables the backrest and footrest to move independently so that a variety of positions can be reached including raising your feet above your heart.
-  The ergonomically designed Toggle Remote has two programmable memory position buttons, performs the lift, up/down recline and footrest functions, and has a convenient USB charger port that allows you to charge your smart device.
-  Lithium battery backup ensures your recliner is operational when the power goes out.

1. VivaLift![®] Power Recliners are FDA Class II Medical Devices designed to aid individuals with mobility impairments.
2. User height ranges should be used as a guide only due to varying factors of each individual user.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

 = Indicates shortest available lead time

STANDARD FABRIC: ASTRO



GREY

BROWN

MUSHROOM

OPTIONAL FABRICS:

CRYPTON[®] ARIA

CAPRICCIO



COOL GREY

ESPRESSO

SMOKE

DOVE



PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

pridemobility.ca