

THE MOST COMFORTABLE RECLINER IN THE WORLD!



COLLECTION	ULTRA ^{1,2,3,7}			RADIANCE ^{1,2,3}				ATLAS PLUS 2 ^{1,2}	TRANQUIL 2				ELEGANCE 2		LEGACY 2	METRO 2	URBANA 2	PERFECTA 2		
SIZE	Large	Medium	Small	Large/Tall	Medium	Petite/Wide	Small	Medium	Large/Tall	Medium	Petite/Wide	Small	Large	Medium	Medium	Medium	Medium	Medium		
MODEL NUMBER	PLR4955L	PLR4955M	PLR4955S	PLR3955LT	PLR3955M	PLR3955PW	PLR3955S	PLR2985M	PLR935LT	PLR935M	PLR935PW	PLR935S	PLR975L	PLR975M	PLR958M	PLR925M	PLR965M	PLR945M		
POSITION*	Infinite Position Weightless Tilt	Infinite Position Weightless Tilt	Infinite Position Weightless Tilt	Infinite Position	Infinite Position	Infinite Position	Infinite Position	Infinite Position	Infinite Position	Infinite Position	Infinite Position	Infinite Position	Infinite Position	Infinite Position	Infinite Position	Infinite Position	Infinite Position	Infinite Position	Infinite Position	
HAND CONTROL	Toggle Remote ⁶	Toggle Remote ⁶	Toggle Remote ⁶	Toggle Remote ⁵	Toggle Remote ⁵	Toggle Remote ⁵	Toggle Remote ⁵	Toggle Remote ⁴	Toggle Remote ⁴	Toggle Remote ⁴	Toggle Remote ⁴	Toggle Remote ⁴	Toggle Remote ⁴	Toggle Remote ⁴	Toggle Remote ⁴	Toggle Remote ⁴	Toggle Remote ⁴	Toggle Remote ⁴	Toggle Remote ⁴	
BACK TYPE	Deluxe Biscuit	Deluxe Biscuit	Deluxe Biscuit	Deluxe Biscuit	Deluxe Biscuit	Deluxe Biscuit	Deluxe Biscuit	Club	Waterfall	Waterfall	Waterfall	Waterfall	Tailored	Tailored	Button	Split	Sewn Pillow	Bergere		
SEAT TYPE	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	
STANDARD FEATURES	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	
ADDITIONAL FEATURES	Heat, Massage, Cup Holder, Wireless Charger			Heat, Cup Holder, Wireless Charger				Cup Holder, Wireless Charger	NA	NA	NA	NA	NA	NA	NA	NA	NA	NA	NA	NA
USER HEIGHT RANGES**	5'10" and above	5'4"-5'10"	5'4" and below	5'10" and above	5'4"-5'10"	5'4"-5'10"	5'4" and below	5'4"-5'10"	5'10" and above	5'4"-5'10"	5'4" and below	5'4" and below	5'10" and above	5'4"-5'10"	5'4"-5'10"	5'4"-5'10"	5'4"-5'10"	5'4"-5'10"	5'4"-5'10"	
SEAT-TO-FLOOR	22.5"	19.5"	18"	22.5"	19.5"	19.5"	18"	19.5"	21.5"	19.5"	18"	18"	20.5"	20"	19.5"	19"	20"	19"	19"	
SEAT WIDTH	23.5"	22.5"	20"	23.5"	22.5"	24"	20"	22"	23.5"	22"	24"	20"	22"	21"	22"	21"	21"	21"	22"	
SEAT DEPTH	22"	21"	18"	22"	21"	18"	18"	21"	22.5"	20.5"	18.5"	18.5"	22"	21"	21.5"	21"	20"	20"	20"	
SEAT DEPTH WITH FOOTREST ELEVATED	39"	38"	40.5"	39"	38"	34"	40.5"	39"	40"	38"	36"	36"	39"	39"	39"	39.5"	39"	39"	43"	
TOP OF BACK TO SEAT	30"	28"	25"	31"	28"	23"	25"	26"	30"	27"	24.5"	25"	30"	27"	25"	26.5"	25"	25"	25"	
CHAIR WIDTH	40"	39"	37"	40"	39"	41"	37"	39"	38"	37.5"	36"	33.8"	37.5"	35.5"	35"	33"	33"	33"	31"	
CHAIR DEPTH	40"	42"	41"	40"	42"	40"	41"	42.5"	43"	42"	39.5"	36.8"	42"	41.5"	42.5"	40.5"	41"	41"	39"	
CHAIR HEIGHT	49"	43"	41"	49"	43"	40"	41"	42"	47"	42"	38"	40"	42"	42"	42"	41"	42"	42"	43"	
DISTANCE FROM WALL	21.5"	15"	15"	21.5"	15"	14"	15"	16"	16"	16"	12.5"	13"	16"	14"	16"	14"	14"	14"	15"	
CHAIR WEIGHT	176 lb/79.8 kg	158 lb/71.7 kg	158 lb/71.7 kg	178 lb/80.7 kg	158 lb/71.7 kg	161 lb/73 kg	156 lb/70.8 kg	164 lb/74.4 kg	155 lb/70.3 kg	145 lb/65.8 kg	143 lb/64.9 kg	137 lb/62.1 kg	157 lb/71.2 kg	152 lb/68.9 kg	153 lb/69.4 kg	142 lb/64.4 kg	145 lb/65.8 kg	142 lb/64.4 kg		
WEIGHT CAPACITY	400 lb/181.5 kg	400 lb/181.5 kg	400 lb/181.5 kg	400 lb/181.5 kg	400 lb/181.5 kg	400 lb/181.5 kg	400 lb/181.5 kg	400 lb/181.5 kg	400 lb/181.5 kg	400 lb/181.5 kg	400 lb/181.5 kg	400 lb/181.5 kg	400 lb/181.5 kg	400 lb/181.5 kg	400 lb/181.5 kg	400 lb/181.5 kg	400 lb/181.5 kg	400 lb/181.5 kg	400 lb/181.5 kg	
WARRANTY	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	

STANDARD FABRIC SELECTION GUIDE	CAPRICCIO				CANYON				BADLANDS	ASTRO			BADLANDS	SAVILLE	SAVILLE	STONEWASH	MERINO		
OPTIONAL FABRIC SELECTION GUIDE	SORRENTO								CRYPTON ARIA	CRYPTON ARIA		CAPRICCIO							

STANDARD FEATURES:



*Not all units achieve all positions. Please refer to specific product.

**Height ranges should be used as a guide only due to varying factors of each individual user.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice.

1. Unit has a stowaway cup holder and cell phone charging pad.
2. The wireless charger is only compatible with smart phones that have wireless charging capabilities. Additionally, no other smart devices, (i.e. smart watches) should be used with the wireless charger.
3. Unit has three different heating pads in the shoulders, lumbar and seat.
4. Toggle remote has two programmable memory buttons.
5. Toggle remote has two programmable memory buttons and one heat button. (Three heat levels standard)
6. Toggle remote has one programmable memory button.
7. Unit has Heat and Massage controls built into the armrest of the chair. (Three Heat and three Massage levels standard)

POSITION GUIDE: *





POWER HEADREST



3.

POWER LUMBAR



FOOTREST EXTENSION

4.



BACKREST

5.



6.

INFINITE POSITION



2.

UPWARD LIFT



TOGGLE REMOTE

7.



8.

THE MOST COMFORTABLE RECLINER IN THE WORLD!

INDEPENDENT FUNCTIONS



(STANDARD ON ALL MODELS)

1. The power headrest lifts your head, neck and shoulders for maximum support and improved sight-line for TV or conversing with others.
2. The upward lift offers convenient stand assist for getting in and out of the recliner.
3. The power lumbar enables you to personalize your comfort and back support.
4. The footrest extension allows you to relax and stretch out completely.
5. The backrest moves up and down independently to help you find the perfect position.
6. Infinite position allows the backrest and footrest to move independently from one another so that a variety of positions can be reached.
7. Toggle remote with memory button feature offers a USB port for charging your mobile device.
8. Power supply with lithium-ion battery backup ensures your recliner is operational when the power goes out.

ADDITIONAL FEATURES

(CHECK SPECIFIC MODELS FOR AVAILABILITY)

- A. Retractable cup holder neatly tucked in the armrest holds your beverage secure when required.
- B. Retractable mobile device wireless charging pad help keep your device ready to go.
- C. Three independent heating pads located at the shoulder, lumbar, and front portion of the seat, help keep you cosy.
- D. Advanced massage system provides gentle pressure to your lower back for a soothing massage.



PATENT PENDING HEADREST/LUMBAR MECHANISM WITH MOTORS

Independently operates the headrest and lumbar functions for enhanced head, neck and back support. Anti-entrapment on headrest for safety.

STRONG WOOD FRAME

Premium furniture-grade hardwood is robust, durable and long-lasting.



STEEL SEAT BOX

Springs are attached to heavy-duty steel frame.

TOTAL COMFORT SEATING

Sinuous springs with plastic-lined clips are attached to morley hinge.

LIFT AND RECLINE MOTORS

Independently adjusts the recliner to achieve relaxed, reclined, lay flat or upward lift positioning.

STEEL FRAME

Heavy-duty metal framework forged to last.



VivaLIFT!
POWER RECLINER COLLECTION

PRIDE
MOBILITY
LIVE YOUR BEST®

pridemobility.ca |

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. VivaLift!® Power Recliners are FDA Class II Medical Devices designed to aid individuals with mobility impairments.