

Zero Turn 10, 4-Wheel

Weight Capacity	400 lbs.
Maximum Speed ¹	Up to 7.3 miles at 400 lbs. Up to 7.5 miles at 200 lbs.
Ground Clearance ³	2" at motor
Turning Radius ³	43"
Overall Length ³	48"
Overall Width ³	24.25"
Seat-to-Ground Height	23"-25"
Seat-to-Deck Height	18"-20"
Front Tires	9" solid
Rear Tires	10.75" solid
Suspension	Front and rear (CTS Suspension)
Range Per Charge ^{1,2}	Up to 15.39 miles at 400 lbs. Up to 21.2 miles at 200 lbs.
Total Weight ⁵	Without Batteries: 185.8 lbs. With Batteries: 246.3 lbs.
Heaviest Piece When Disassembled ³	75 lbs. (front section)
Standard Seating ³	Type: Foldable, deluxe stadium, compact Weight: 34 lbs. Material: Black vinyl Dimensions: Width: 18", Depth: 17"
Drive System	24-volt DC motor, rear-wheel drive, dual in-line
Dual Braking System	Regenerative and electromechanical
Battery Requirements ^{4,6}	Size: (2) 40 Ah included Weight: 30.25 lbs. each
Battery Charger	Off-board, 3.5 amp
Warranty	Frame: lifetime limited Drivetrain: 2-year limited Electronics: 2-year limited Battery: 6-months

1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.

2 Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.





4 Battery weight may vary based on battery manufacturer.

5 Includes standard seat weight.

6 AGM or gel cell type required.

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

Beamsville, ON • 888-570-1113

pridemobility.ca |    

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

ZERO
TURN10

INDEPENDENCE
REIMAGINED



LIVE YOUR BEST[®]

pridemobility.ca |



CUTTING-EDGE VERSATILITY

THE FIRST 4-WHEEL INDOOR/OUTDOOR SCOOTER!



18"x17"
Black vinyl,
deluxe stadium
seat with swivel
and sliders,
height, width,
depth and angle
adjustable armrests

Tiller
adjustment
lever

Ambient voltmeter

Standard large
front basket

Full LED light
package includes
LED headlights
and rear lights
with sequential
turn signals and
puddle light
under the tiller



Front LED lights



Tiller adjustment lever



Pride's patent-pending, intelligent-turning, iTurn Technology™ provides 4-wheel stability with 3-wheel maneuverability

With contemporary styling and cutting-edge technology, the **Zero Turn 10** is the perfect combination of aggressive performance and maneuverability. Powerful dual motors and two-wheel drive offers aggressive traction, easily tackling rugged terrain on trails. When you're done spending time outdoors, the Zero Turn 10 transitions seamlessly indoors. Maneuver smoothly through the halls and doorways with EZ Turn Technology™.

Additional Features

- Three-speed control lets you adjust your speed when moving from indoor to outdoor driving
- Delta tiller with wraparound handles
- Easy access tiller-mounted XLR charger port
- USB charger integrated in the tiller for convenient charging of smart phones and portable electronic devices
- Integrated storage/cup holder in the tiller

Available Colours



Accessories*

- Cane/crutch holder (*single, double*)
- Cell phone holder
- Cup holder
- Forearm crutch holder
- High-back seat
- 50", 60", 70" lap belt
- Oxygen tank holder
- Quad cane holder
- Rear basket
- Rearview mirror
- Saddlebag
- Walker holder
- Weather cover
- Wishbone crutch holder

*Only one rear accessory can be used at a time.

**Matte finish