

# 4-Wheel Stability with 3-Wheel Maneuverability

## ZERO TURN8



Patented, intelligent-turning

**iTURN**<sup>®</sup>  
TECHNOLOGY



Get a 4-wheel scooter with a 3-wheel turning radius with the Zero Turn 8 and Pride's patented, intelligent-turning, iTurn Technology™. Navigate tight corners and small spaces effortlessly with a 38.25" turning radius. Illuminate your ride with bright LED lighting. Dual motors and CTS suspension ensure a smooth and comfortable drive.

PRIDE<sup>®</sup> MOBILITY SCOOTERS  
LIVE YOUR BEST<sup>®</sup>


| [pridemobility.ca](http://pridemobility.ca) |





## Specifications

Weight Capacity	325 lbs.
Maximum Speed <sup>1</sup>	5.9 mph
Ground Clearance <sup>2</sup>	1.3" at motor
Turning Radius <sup>2</sup>	37.25"
Overall Length <sup>2</sup>	45.75"
Overall Width <sup>2</sup>	20"
Seat-to-Deck Height	17"-18.5"
Seat-to-Ground Height	21.5"-23.5"
Front Tires	2" x 7" solid
Rear Tires	2.5" x 9" solid
Range Per Charge <sup>2,4</sup>	Up to 11.6 miles at 325lbs/13.1 miles at 200lbs
Total Weight Without Batteries (standard seat)	120 lbs.
Heaviest Piece When Disassembled	48 lbs. (rear section)
Standard Seating	Type: Foldable, compact, molded plastic Weight: 27 lbs. Material: Black vinyl Dimensions: Width: 17" Depth: 17"
Motor Output	24-volt DC motor, rear-wheel drive, two-motor
Brakes	Electronic, regenerative and electromechanical
Battery Requirements <sup>3,5</sup>	Size: (2) 22they AH Weight: Battery box 31.5 lbs. 13.5 lbs. each
Battery Charger	Off-board, 3.5 amp
Delta Tiller	Yes
Colors	*Note color may appear more silver with hues of green in some lighting 



## Features

- The first 4-wheel indoor scooter
- Patented iTurn Technology™ provides 4-wheel stability with a 3-wheel turning radius
- Charging port built into the tiller for easy charging
- Contemporary new rim design makes the Zero Turn 8 look as good as it drives
- CTS suspension for comfortable long range use
- Delta tiller, ergonomic wraparound handles enables easy-to-use left or right hand operation
- Dual motors for added traction in rough terrain
- Feather touch disassembly permits simple frame separation with only one hand
- Standard front basket
- Full LED lighting package for safe operation in low light settings
- USB charger built into the tiller for convenient charging of smart phones and portable devices

## Accessories\*

- Rear basket
- Oxygen holder
- Cell phone holder
- Cup holder
- Forearm crutch holder
- 50" lap belt
- Micro, black, weather cover
- Rear view mirror
- Saddle bag
- Safety flag
- Single rear mount, 3" cane holder
- Double rear mount, 3" cane holder
- Walker holder
- Wishbone crutch holde

\*Only one rear accessory can be used at a time.

## Option

- 22 AH battery box (additional spare)
- Dual under seat storage bins

<sup>1</sup> Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).

<sup>2</sup> Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

<sup>3</sup> AGM or gel cell type required.

<sup>4</sup> Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

<sup>5</sup> Battery weight may vary based on manufacturer.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for scooter only. Overall measurements will vary based on seating and accessory selections.