

Part Number	SC54ALX
Weight Capacity	325 lbs.
Maximum Speed ²	Up to 4.8 mph
Ground Clearance ⁴	1.39" at motor
Turning Radius ⁴	38.39"
Overall Length ⁴	42.48"
Overall Width ⁴	20.5"
Seat-to-Deck Height	16.18" – 17.72"
Seat-to-Ground Height	21.85" – 23.54"
Suspension	Full front and rear suspension
Tires	Front: 2"x 8" Rear: 2.5"x 8"
Range Per Charge ^{2,3}	16 Ah: Up to 18 miles at 200 lbs. 16 Ah: Up to 13.6 miles at 325 lbs. 8 Ah: Up to 9 miles at 200 lbs. 8 Ah: Up to 6.8 miles at 325 lbs.
Product Weight	83.2 lbs. without batteries
Heaviest Piece when Disassembled ⁴	35.2 lbs.
Standard Seating	Type: Mesh Back Material: Black Vinyl Dimensions: 17"x17" Weight: 23 lbs.
Motor Output	24-volt DC motor, sealed transaxle
Brakes	Regenerative and electromechanical
Battery Requirements ⁵	Type: Lithium battery box Size: 8 Ah or 16 Ah
Battery Weight	8 Ah: 8 lbs. (3.5 kg) 16 Ah: 12.2 lbs. (5.5 kg)
Battery Charger	Off-board 3.5 amp lithium battery charger
Warranty	Lifetime limited on frame; 2-year limited on electronics; 2-year limited on drivetrain 2-year limited on battery

HCPCS Code: K0899

- ¹ Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
- ² Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- ³ Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- ⁴ Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- ⁵ Battery weight may vary based on manufacturer.



GOGO[®]
ENDURANCE AL+

Lightweight mobility **plus**
so much more



PRIDE
MOBILITY[®]

800-800-8586 • pridemobility.com



Lightweight
aluminum
frame



Unrivaled
tight turning



Charge360
magnetic
charging



Front
and rear CTS
suspension

*K0899 Power Mobility Device not coded by DME PDAC. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage, and bundling guidelines may apply.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

081425

PRIDE
MOBILITY[®]



pridemobility.com

The Go Go® Endurance AL+ is an FDA Class II Medical Device designed to aid individuals with mobility impairments.