

Go Go® Super Portable

Model Number	SC15M
Weight Capacity	300 lbs.
Turning Radius	51.6"
Overall Width ³	21.5"
Overall Length ³	37.2"
Maximum Speed ²	Up to 4.5 mph
Ground Clearance ³	3.7" at mid-deck 1.85" at motor
Front Wheels	7" solid
Rear Wheels	8" solid
Braking System	Regenerative and electromechanical
Seat Size	15.7"W x 14.7"D
Seat-to-Deck	16.6"
Per Charge Range ^{3,4} with 11.7 Ah battery	Up to 10 miles at 220 lbs. Up to 8.8 miles at 300 lbs.
Total Weight	41.6 lbs. (without battery) 49.9 lbs. with 11.7 Ah battery
Battery Requirements	11.7Ah Lithium-ion
Battery Charger	Off-board
Warranty	Frame: 3-year limited Electronics: 2-year limited Drive Motors: 2-year limited Batteries: 2-year limited

- ¹ FDA Class II Medical Devices are designed to aid individuals with mobility impairments
- ² Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- ³ Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- ⁴ Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- ⁵ This accessory must be removed prior to using the travel bag.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.



UNFOLD



YOUR



WORLD



PRIDE MOBILITY®

800.800.8586

| pridemobility.com |



HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

The information contained herein is correct at time of publication; we reserve the right to alter specifications without prior notice.

012926

GO GO™ SUPER PORTABLE



UNFOLD YOUR WORLD

pridemobility.com |



FDA Class II Medical Device¹

UNFOLD YOUR WORLD



Our exclusive Charge360 magnetic charging port

Optional Armrests

Carbon fiber seat support

11.7 Ah lithium-ion battery is airline compliant**

Magnesium rims

USB port for charging your mobile devices

Front storage bag (standard)

Adjustable and telescoping carbon fiber tiller

Convenient handle for towing while unit is folded

Storage compartment under seat

Front and rear lighting



KEY FEATURES

- Offers a simple folding process without needing to remove the seat.
- Overall depth is reduced to 12", making transport and storage fast and convenient.
- Weighs only **41.6 lbs.** without the battery
- Our **exclusive Charge360 technology** requires minimal dexterity
- 11.7 Ah **airline-compliant lithium battery** makes this an ideal choice for air travel**
- Front and rear lighting for both seeing and being seen
- Integrated USB port for charging mobile devices
- Stands upright when folded and features a handle for easy towing
- Storage compartment located under seat
- Accepts optional armrests
- Front storage bag



COLOR OPTIONS



Black

Mint

True Blue

Orange

ACCESSORIES

- Arm rests
- Cup holder⁵
- Cell phone holder⁵
- Rear view mirror(s)⁵



GO GO® TRAVEL BAG BUNDLE*

Conveniently tow your Go Go® Super Portable scooter through airports and parking lots with this durable, protective travel bag bundle.

*Currently only compatible with the Go Go® Super Portable

FDA Class II Medical Device¹

** Transport of this item on commercial aircraft is subject to federal regulations and airline policy. Please consult your airline prior to flight to ensure your Go Go® Super Portable and its 11.7 Ah lithium-ion battery are permitted.



Folds to a depth of only 12" without having to remove the seat.

Our exclusive Charge360 technology features easy-to-use magnetic charging.

Front and rear lighting for increased visibility.