

Model Number	SC493
Weight Capacity	300 lbs.
Maximum Speed²	Up to 4.3 mph
Ground Clearance⁴	1.29" at motor 3.5" under deck
Turning Radius⁴	48.75"
Overall Length⁴	41.50"
Overall Width⁴	19.25"
Seat-to-Deck Height	16.5" - 18"
Seat-to-Ground Height	20.75" - 22.25"
Tires	Front: 2"x 7" Solid Rear: 2"x 8" Solid
Range^{2,3}	Up to 7.1 miles at 300 lbs Up to 8.9 miles at 220 lbs.
Product Weight	83.4 lbs. without batteries 103.4 lbs. with batteries
Heaviest Piece⁴	36.2 lbs. (front section)
Standard Seating	Type: Compact foldable Material: Black Vinyl Dimensions: 17"x17" Weight: 23 lbs.
Motor Output	Rear-wheel drive, single motor, 24-volt DC motor
Brakes	Regenerative & electromechanical
Batteries⁵	(Two) 12-Ah, sealed lead-acid, non-spillable deep-cycle
Battery Weight	12Ah battery pack: 20 lbs. each
Battery Charger	2 amp off-board charger
Warranty	Lifetime limited on frame; 2-year limited on Electronics; 2-year limited on Drivetrain; 6-month limited on Batteries

HCPCS Code: K0800 Group 1 Standard

- 1 Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
- 2 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- 3 Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 4 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- 5 Battery weight may vary based on manufacturer.

Ask your Authorized Pride® Dealer for a test drive today!



**PRIDE®
MOBILITY**

pridemobility.com | 800.800.8586

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

GO GO®
ULTRA X GEN 3



pridemobility.com

ROOM FOR THE RIDE



NEW deck design that offers more leg room and increased seating comfort.



Feather-touch disassembly easily separates into 5 lightweight pieces



Soft-touch wraparound delta tiller

The **Go Go[®] Ultra X Gen 3** is equipped with an ergonomically designed deck for increased leg room and more comfortable seating. Its wraparound delta tiller provides additional comfort, while feather-touch disassembly makes it easy to store and transport.

FEATURES

- **NEW** ergonomic deck allows for more pass-through space, leg room, and comfort
- **Soft-touch, easy-to-grip** wraparound delta tiller
- 17" x 17" compact foldable seat
- Standard front basket
- **Feather-touch disassembly** separates the unit into five lightweight pieces
- Charging ports are located on the tiller and battery box

ACCESSORIES*

- | | |
|---------------------|----------------------------|
| • Cell phone holder | • Saddle bag |
| • Cup holder | • Safety flag |
| • Lap belt | • Single/Double rear mount |
| • Military patches | 3" cane/crutch holder |
| • Oxygen holder | • Under seat storage |
| • Rear basket | • Walker holder |
| • Rear view mirror | • Weather cover |

* Only one rear accessory can be used at a time.

WARRANTY

FRAME: Lifetime limited
ELECTRONICS: 2-year limited
DRIVETRAIN: 2-year limited
BATTERY: 6-month limited