

MODEL

Jazzy Air MED ¹	JAZZYAIRMED
Weight Capacity	300 lbs.
Turning Radius	24.41"
Overall Width ³	23.25"
Overall Length ³	33.86" / 41.54" (with footplate)
Maximum Speed ²	Up to 3.71 mph (elevated and non-elevated position)
Ground Clearance ³	1" at motor bracket
Front Wheels	6" solid casters
Drive Wheels	10" solid
Rear Wheels	6" solid casters
Drivetrain	Two motor, Mid-Wheel 6 [®] drive
Braking System	Regenerative and electromechanical
Suspension Type	Active-Trac [®] (ATX) with Mid-Wheel 6 [®] Technology
High-back Seat Sizes	16"x16"-18", 18"x18'-20", 20"x 18"-20"
Seat-to-Floor Range	21.25" – 33"
Per Charge Range ^{3,5}	Up to 13 miles at 300 lbs. Up to 14.1 miles at 220 lbs.
Standard Electronics	40A, PG VR2
Component Weights ⁷	197.2 lbs. without batteries With batteries: 241.6 lbs.
Battery Requirements ^{4,6}	(2) 12 volt, U-1
Battery Charger	5 amp, off-board
Warranty	Frame: 5-year limited Drivetrain: 13-month Electronics: 13-month Battery: 13-month

HCPCS Code: K0830/K0831

- Jazzy[®] Power Chairs are FDA Class II Medical Devices designed to aid individuals with mobility impairments.
- Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (10%-5%).
- Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- AGM or gel cell type required.
- Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- Battery weight may vary based on manufacturer.
- Weight includes base, seat and batteries.
- WC-19 Occupied Transit.



FDA Class II Medical Device¹

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

(US) 800-800-8586 • (Canada) 888-570-1113

FOLLOW US      

pridemobility.com

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

The information contained herein is correct at time of publication; we reserve the right to alter specifications without prior notice.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

020725

Jazzy¹
air[®]
MED

**ELEVATION
MADE SIMPLE**



12"
**OF
LIFT**

FDA Class II
Medical Device¹

AMERICA'S FAVORITE POWER CHAIR

FOLLOW US      

pridemobility.com

REACHING NEW HEIGHTS

High-back or solid seat pan, limited recline, memory foam captain seat. Seat features 2" of depth adjustment. Choice of 16"x16"-18", 18"x18"-20", 20"x18"-20" seat sizes

Height and angle adjustable headrest



VR2 controller

12" of power adjustable seat height

23.25" base width

Durable, color-through shroud with matte finish

airMED

FDA Class II Medical Device¹



Mid-Wheel 6[®] Drive and Active-Trac[®] ATX Suspension



Non-scuffing, flat free tires



12" seat elevation allows access to hard-to-reach areas

Raise your Expectations!

Elevating up to 12", the Jazzy Air[®] MED offers users greater functionality in completing activities of daily living. This medical device drives safely up to 3.71 mph while elevated, allowing users to enjoy face-to-face social engagement. With Mid-Wheel 6[®] Drive and Active-Trac[®] ATX Suspension, the Jazzy Air MED offers enhanced performance and stability both indoors and out.

- 300 lb. weight capacity
- Up to 3.71 mph maximum speed while elevated
- Angle, width, height and depth adjustable armrests
- Foot platform with 2" of adjustment for comfort in both the elevated and non-elevated positions
- Dual-sided caster forks for enhanced maneuverability

Available Color

Matte Black

Accessories

- Cane holder
- Cell phone holder
- Cup holder
- Flag kit
- 50", 70" Lap belt
- Oxygen tank holder
- Rear basket
- Saddlebag
- Walker holder
- Weather cover
- XLR USB charger

Options

- Swing-away joystick
- Stealth[®] cushions
- Elevating leg rests
- Swing-away leg rests
- Heel loops

Military Patches

