

Part Number	Jazzy® Carbon 27X
Weight Capacity	275 lbs.
Maximum Speed ¹	Up to 3.7 mph
Ground Clearance ³	3"
Turning Radius ³	35.4"
Overall Length ³	41"
Overall Width ³	23"
Seat-to-Ground Height	19.7"
Tires	Front: 7" Rear: 10"
Range Per Charge ^{1,2}	Up to 8 miles per charge
Product Weight ^{3,5}	27.5 lbs. without batteries
Standard Seating ³	Cushioned back with removable seat cushion
Folded Dimensions	34" L x 12" W x 30" H
Motor Output	24V 150W
Brakes	Regenerative and electromechanical
Battery Requirements ⁶	24V 10Ah
Battery Weight ⁴	3.5 lbs
Battery Charger	3.5A
Warranty	5-year limited on frame; 13-month limited on electronics 13-month limited on drivetrain 13-month on battery

HCPCS Code: K0899

POWER MOBILITY DEVICE, NOT CODED BY DME PDAC

1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.

2 Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.

4 Battery weight may vary based on battery manufacturer.

5 Includes standard seat weight.

6 Transport of this item on commercial aircraft is subject to federal regulations and airline policy. Please consult your airline prior to your flight to ensure your Jazzy® and its lithium-ion battery are permitted.

* Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.

Jazzy®

carbon 27X

Weighs just
27lbs
without the
battery

PRIDE MOBILITY®

800-800-8586 • pridemobility.com



HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.



FDA Class II
Medical Device*



FDA Class II
Medical Device*

Jazzy[®] carbon 27X

Weighs just
27lbs
without the
battery



Folds & unfolds
for easy transport



Advanced durable
carbon fiber frame



X-frame design for
vertical folding



Latch for locking unit
open or closed

YOUR MOBILITY X-FACTOR.

The Jazzy[®] Carbon 27X folds vertically in just a few seconds thanks to its X-frame design, and it weighs a mere 27.5 lbs. without the battery. Comprised of super-strong carbon fiber, this power wheelchair is durable enough to withstand all your adventures and light enough to keep you moving.

FEATURES

- Unit weighs 27.5 lbs. without the battery
- 275-lb. weight capacity
- Top speed up to 3.7 mph
- X-frame design allows the unit to fold nearly flat in vertical position
- Advanced carbon fiber frame offers durability and light weight
- Folded dimensions: 34" L x 12" W x 30" H
- Flip up lightweight carbon fiber arms
- Adjustable Foot Rests
- Airline-compliant lithium battery** (10 Ah only)
- Storage pocket in seat back
- Cushion and adjustable lap belt (standard)
- Programmable quick disconnect joystick (standard)
- Carbon fiber Storage with Cup Holder (standard)
- Collapsible handles for easy storage

COLOR OPTIONS

Black

OPTIONAL ACCESSORIES

Additional 10 Ah
airline-compliant battery**
Battery Travel Bag
Jazzy Travel Kit (Coming Soon)
Military Patches
Simplicity Cushion 18" x 18"
Mesh Cup Holder