

Jazzy Elite 14-2S-SS¹	Group 2: Standard (K0822)
Jazzy Elite 14-2S-C¹	Group 2: Standard (K0823)
Weight Capacity	300 lbs.
Turning Radius	24"
Overall Width³	24.25"
Overall Length³	35.5"/42.5" (with footplate)
Maximum Speed²	Up to 4 mph
Ground Clearance³	3" at battery tray
Front Wheels	3" solid anti-tips
Drive Wheels	14" foam filled
Rear Wheels	6" solid casters
Drivetrain	Two motor, front-wheel drive
Braking System	Regenerative and electro-mechanical
Suspension Type	Limited
High-back Seat Sizes	18" x 18-20", 20" x 18-20" and 22" x 20-22"
Seat-to-Floor Range	20"-22"
Per Charge Range^{3,5}	Up to 18 miles @ 300 lbs.; Up to 21 miles @ 200 lbs.
Standard Electronics	40A, Dynamic LiNX controller
Component Weights	Base: 99 lbs. Seat: 18" wide: 51.8 lbs. 20" wide: 54 lbs. 22" wide: 57.6 lbs. Batteries: 37 lbs. each
Battery Requirements^{4,6}	(2) 12 volt, NF-22 batteries
Battery Charger	5 amp, off-board
Battery HCPCS Code	E2361
Warranty	Frame: 5-year limited Drivetrain: 13-month Electronics: 13-month Battery: 13-month

- Jazzy® Power Chairs are FDA Class II Medical Devices designed to aid individuals with mobility impairments.
- Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- AGM or gel cell type required.
- Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- Battery weight may vary based on manufacturer.



**PRIDE[®]
MOBILITY**

800-800-8586 • pridemobility.com



NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

The information contained herein is correct at time of publication; we reserve the right to alter specifications without prior notice.

060526

Jazzy[®]
ELITE 14

**PEAK PERFORMANCE
INDOORS AND OUT**



FDA Class II
Medical Device

800-800-8586 • pridemobility.com



Jazzy®

ELITE 14

Depth-adjustable, limited-recline, fold-flat, high-back seat

Height, width, and angle-adjustable armrests

Jazzy® Red color-through shroud

Convenient rear access to batteries

FDA Class II Medical Device

Intuitive 40A Dynamic LiNX controller

Large 14" drive wheels

Flip-up foot platform

FEATURES

- Large 14" drive wheels provide enhanced outdoor performance
- Color-through shroud
- Depth-adjustable, limited-recline, fold-flat, high-back seat
- Universal, fixed, in-line joystick mount
- Flip-up foot platform
- Black, flat-free, non-marking tires
- Easy-access free-wheel lever
- 60" lap belt

OPTIONAL SEATING

Deluxe contour, limited recline, high back with solid seat pan, available sizes:

- 16-20" (W) x 16-20" (D)

Limited recline high back, available sizes:

- 16" (W) x 16" (D)

OPTIONS

- Bell shaped foot platform
- Elevating leg rests
- Heel loops
- Height adjustable foot platform bracket
- Residual limb support
- Swing-away legrests
- Swing-away joystick

ACCESSORIES

- Attendant control bracket
- Cane/crutch holder
- Cup holder
- Flag holder
- Lap belt (50/70/80/90")
- Military patches
- Oxygen tank holder
- Rear basket
- Saddle bag
- Walker holder
- Weather cover
- XLR USB charger



Easily accessible freewheel levers



Intuitive 40A Dynamic LiNX controller



Convenient rear access to batteries