

Jazzy® EVO 613¹	Group 2: Standard (K0822/K0823)
Jazzy® EVO 613 Li¹	Group 2: Standard (K0822/K0823)
Weight Capacity	300 lbs.
Turning Radius³	22"
Overall Width³	22"
Overall Length³	41" with footrest (35" without)
Maximum Speed²	Up to 5.1 mph
Ground Clearance³	3"
Front Casters³	5"
Drive Wheels³	13" solid
Rear Casters³	5"
Drivetrain	Two motor, in-line, mid-wheel drive
Braking System	Regenerative and electromechanical
Suspension Type	Active-Trac® (ATX)
High-back Seat Sizes	16"x16-18" 18"x18-20" 20"x18-20"
Seat-to-Floor Range	19.75" - 21.25"
Per Charge Range^{3,5}	EVO 613: 14 miles at 300 lbs. (17.1 miles at 200 lbs.) EVO 613 Li: 17.6 miles at 300 lbs. (21 miles at 200 lbs.)
Electronics	40A, PG GC3
Component Weights⁶	Base: 94 lbs. Seat: 48 lbs. (18" seat) U1 Batteries: (2) 23.5 lbs. each; Lithium battery: (1) 18.5 lbs.
Total Weight⁷	EVO 613: 189 lbs. EVO 613 Li: 160.5 lbs.
Batteries	EVO 613: (2) 12 volt, U-1 EVO 613 Li: (1) Lithium iron phosphate
Battery Charger	3.5 amp, off-board
Warranty	Frame: 5-year limited Drivetrain: 13-month limited Electronics: 13-month limited SLA Battery: 13-month limited Lithium Battery: 3-year limited

- ¹ Jazzy® Power Chairs are FDA Class II Medical Devices designed to aid individuals with mobility impairments.
- ² Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- ³ Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- ⁴ AGM or gel cell type required.
- ⁵ Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- ⁶ Battery weight may vary based on manufacturer.
- ⁷ Total weight includes base, seat and batteries.



**PRIDE®
MOBILITY**

800-800-8586 • pridemobility.com



HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

The information contained herein is correct at time of publication; we reserve the right to alter specifications without prior notice.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

060926

**Jazzy®
EVO 613 SERIES**

Always insist on **lithium!**



800-800-8586 • pridemobility.com



The Jazzy® EVO 613 is an FDA Class II Medical Device designed to aid individuals with mobility impairments.

THIS IS JAZZY EVOLVED

Height, width,
depth and angle
adjustable armrests

GC3 controller with
quick disconnect

Depth-adjustable,
limited recline,
memory foam
high-back seat
with swivel

Narrow base
width of 22"

Convenient access
to the controller and
batteries while user
remains seated

13" drive wheels

Removable
color-through
shrouds

Mid-Wheel 6® Technology
provides maximum stability

FEATURES

- Maximum speed up to 5.1 mph
- 13" mid-wheel drive tires
- Narrow base width of only 22"
- Removable matte-finish color-through shrouds
- Access to the controller and batteries without the user needing to exit the chair
- Depth-adjustable, limited recline, memory foam single seat post allows the user to swivel
- Pride's patented Active-Trac® (ATX) Suspension
- Mid-Wheel 6® Technology

AVAILABLE COLORS (MATTE)



Robin's
Egg Blue



Iceberg
Blue



Sugar
Plum



Black



White



Red

ACCESSORIES

- Attendant control bracket
- Cane/crutch holder
- Cell phone holder
- Cup holder
- Flag kit
- Lap belt
- Military patches
- Oxygen tank holder
- Rear basket
- Saddlebag
- Walker holder
- Weather cover
- 2.1A XLR USB charger
- Stealth® cushions: Zen/Essence/Glacial

OPTIONS

- Depth-adjustable solid seat pan, limited recline, high-back with swivel
- Elevated leg rests
- Swing-away leg rests
- Unoccupied transit
- Heel loops
- Residual limb support



Lithium Iron Phosphate Battery*
(Available on Jazzy® EVO 613 Li)

- Faster charging time
- Sustained power performance
- Lighter weight
- More useable energy
- Longer life

* NOT AIRLINE APPROVED.