

<b>Model</b>	Victory® Platinum (4-Wheel)
<b>Part Number</b>	SC710PLA
<b>Weight Capacity</b>	400 lbs.
<b>Maximum Speed</b> <sup>2</sup>	Up to 6 mph
<b>Ground Clearance</b> <sup>4</sup>	2.75" at motor & transaxle
<b>Turning Radius</b> <sup>4</sup>	45.5"
<b>Overall Length</b> <sup>4</sup>	47"
<b>Overall Width</b> <sup>4</sup>	24"
<b>Seat-to-Deck Height</b>	17.5" – 19"
<b>Seat-to-Ground Height</b>	23.5" – 25"
<b>Suspension</b>	Full front and rear suspension
<b>Tires</b>	Front: 3.5"x 10.4" Rear: 3.5"x 11"
<b>Range Per Charge</b> <sup>2,3</sup>	Up to 15.5 miles at 400 lbs.
<b>Product Weight</b>	145 lbs. without batteries 188.2 lbs. with batteries
<b>Heaviest Piece when Disassembled</b> <sup>4</sup>	67 lbs. (front section)
<b>Standard Seating</b>	Type: Compact, foldable Material: Black Vinyl Dimensions: 18"x17" Weight: 30.4 lbs.
<b>Motor Output</b>	Rear-wheel drive, sealed transaxle, 24-volt DC motor
<b>Brakes</b>	Regenerative and electromechanical
<b>Battery Requirements</b> <sup>5</sup>	(Two) 12-volt, sealed lead-acid, non-spillable deep-cycle
<b>Battery Weight</b>	21.6 lbs. each
<b>Battery Charger</b>	3.5 amp magnetic charger
<b>Warranty</b>	Lifetime limited on frame; 2-year limited on Electronics; 2-year limited on Drivetrain

## HCPCS Code: K0899

- 1 Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
- 2 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- 3 Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 4 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- 5 Battery weight may vary based on manufacturer.



FDA Class II Medical Device<sup>1</sup>

**PRIDE  
MOBILITY®**

800-800-8586 • [pridemobility.com](http://pridemobility.com)



\*K0899 Power Mobility Device not coded by DME PDAC. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage, and bundling guidelines may apply.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

**VICTORY®**  
PLATINUM

PROUDLY MADE IN  
**AMERICA**  
OF FOREIGN AND DOMESTIC PARTS

**EXPERIENCE THE  
PLATINUM DIFFERENCE**



FDA Class II  
Medical Device<sup>1</sup>

# MADE IN AMERICA

- BUILT IN PONTOTOC, MISSISSIPPI WITH FOREIGN AND DOMESTIC PARTS •



FDA Class II  
Medical Device<sup>1</sup>



Easy one-handed  
tiller adjustment



Featuring Pride's patented **EZ TURN**  
technology for tight turning



Magnetic charging with  
**CHARGE360** technology

Proudly assembled in the U.S., the Victory® Platinum 4-wheel scooter offers a premium mobility experience. Patented EZ Turn technology helps you navigate tight spaces with 3-wheel precision. Exclusive Charge360 technology keeps you powered up with magnetic charging. And with full front and rear suspension, take on uneven terrain with confidence.

## FEATURES

- **MADE IN AMERICA**, in Pontotoc, Mississippi from foreign and domestic parts
- **EXCLUSIVE EZ TURN** technology for an impressive turning radius of only 45.5"
- **CHARGE360** magnetic charging port requires minimal dexterity
- One-handed adjustable tiller makes it easier to adjust the angle of the tiller
- **FULL FRONT AND REAR SUSPENSION** with optional suspension seat post
- Feather-touch disassembly for quick and convenient transport or storage
- Water resistant console
- Front basket and pair of rearview mirrors (standard)

## COLOR OPTIONS

MATTE

Cardinal Red

MATTE

Iceberg Blue

## ACCESSORIES\*

- |                         |                             |
|-------------------------|-----------------------------|
| • Cane Clip Attachment  | • Saddlebag                 |
| • Cane Crutch Holder    | • Safety Flag               |
| • Cell Phone Holder     | • Single/Double Cane Holder |
| • Cup Holder            | • Suspension seat post      |
| • Forearm Crutch Holder | • Walker Holder             |
| • Lap Belt              | • Weather Cover             |
| • Military Patches      | • Wishbone Crutch Holder    |
| • Oxygen Tank Holder    | • XLR USB Charger           |
| • Rear Basket           |                             |

- |   |   |
|---|---|
| • Cup holder, rear basket, and saddlebag are available as a bundle. | • High back seat, with limited recline, is available as an upgrade. |
|---|---|

\* Only one rear accessory can be used at a time.