

Accu-Trac

ADVANCED TRACKING TECHNOLOGY


Accu-Trac Advanced Tracking Technology creates strikingly intuitive power chair drive performance by utilizing motor encoders that are in constant communication with Q-Logic 3 and Q-Logic 2 Drive Controls. It's advanced technology with a simple goal: a power chair that tracks straighter, therefore, requiring fewer steering corrections for users with specialty controls.

PRODUCT FEATURES

- Automatically distributes torque for enhanced obstacle climbing
- Synchronizes power for straighter tracking
- Requires fewer steering corrections for specialty controls

CLIENT BENEFITS

- More efficient operation
- Reduced fatigue
- Enhanced independence



Edge® 3

Available on most
Quantum® bases.


**PROUDLY
BUILT IN
CANADA**
WITH IMPORTED PARTS


QUANTUM®
#1 FOR REHAB POWER

 | quantumrehab.ca

Beamsville, ON • 888-570-1113