

ADI 2.0 | STYLES

Aluminum  
**AFT**  
*Flat Top*

Aluminum  
**AL**  
*Curved Top*



**Solid Contour Support**  
Available in: **2in, 4in & 6in** Depths

*Widths* **10 to 24 in.**

*Heights*



*Contour Options*

**Active**

2"



**Deep**

4"



**Extra Deep**

6"



*Heights*



*Contour Options*

**Active**  
2"



**Deep**  
4"




**Extra Deep**  
6"




ADI 2.0 | CONTOUR



2"



4"



6"

*Widths*  
**10" to 14"**

*Heights*  
**Niño** ● 8"  
**Low** ● 10"  
**Medium** ● 13"  
**Standard** ● 16"




*Mounting Option*  
**BR1129**  
*Quick Release Pro™*



*Alternative Mounting*  
**Adapter Plate**

**SSAPGEN-PED**

Fits Merit®, Permobil®, Rovi® and other PWCs  
with similar mounting configurations.



**SSAPT84-PED**


Fits Quantum Rehab® TB3 & TB4



# PWC MOUNTING

ALTERNATIVE MOUNTING

**SIMPLE SOLUTIONS ADAPTER PLATE (INTERFACE)**



**SSAPT B4**

Fits Quantum Rehab® TB4 and TB3 with a weight limit of 350lbs.



**SSAP GEN**

Fits Merit®, Permobil®, Rovi® and other PWCs with similar mounting configurations. Available also for pediatric sizes with item no. **SSAP GEN-PED**



**SSAPT B4-HD**

Fits Quantum Rehab® TB3-HD with an increased weight limit of 450lbs.

# CANE MOUNTING

Quick Release

**BR1102 BR1129 PRO**



Removable  
Width | Depth | Angle  
Adjustable

2-Point Pro

**BR1103**



Fixed  
Width | Depth | Angle  
Adjustable

4-Point Pro

Width | Depth | Angle - Adjustable

**BR1114**



**BR1114T**



## ADI™ 2.0 BACKS DIMENSIONS AVAILABILITY

		Width:	10	11	12	13	14	15	16	17	18	19	20	22	24	
OVERALL HEIGHT	Niño <i>(Pediatric)</i> 8in	Active	AL	AL	AL	AL	AL	-	-	-	-	-	-	-	-	
		Deep	AL	AL	AL	AL	AL	-	-	-	-	-	-	-	-	
		Extra Deep	-	-	-	-	-	-	-	-	-	-	-	-	-	
	Low 10in	Active	AL	AL	AL	AL	AL	AL	AL	AL	AL	AL	AL	AL	AL	AL
		Deep	AL	AL	AL	AL	AL	AL	AL	AL	AL	AL	AL	AL	AL	AL
		Extra Deep	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	Medium 13in	Active	-	-	AL	AL	AL	AL	AL	AL	AL	AL	AL	AL	AL	AL
		Deep	-	-	AL	AL	AL	AL	AL	AL	AL	AL	AL	AL	AL	AL
		Extra Deep	-	-	-	-	-	-	AL	AL	AL	AL	AL	AL	AL	AL
	Standard 16in	Active	-	-	-	AL • AFT	AL • AFT	AL • AFT	AL • AFT	AL • AFT	AL • AFT	AL • AFT	AL • AFT	AL • AFT	AL • AFT	AL • AFT
		Deep	-	-	-	AL • AFT	AL • AFT	AL • AFT	AL • AFT	AL • AFT	AL • AFT	AL • AFT	AL • AFT	AL • AFT	AL • AFT	AL • AFT
		Extra Deep	-	-	-	-	-	-	AL • AFT	AL • AFT	AL • AFT	AL • AFT	AL • AFT	AL • AFT	AL • AFT	AL • AFT
	Tall 18in	Active	-	-	-	-	-	-	AL • AFT	AL • AFT	AL • AFT	AL • AFT	AL • AFT	AL • AFT	AL • AFT	AL • AFT
		Deep	-	-	-	-	-	-	AL • AFT	AL • AFT	AL • AFT	AL • AFT	AL • AFT	AL • AFT	AL • AFT	AL • AFT
		Extra Deep	-	-	-	-	-	-	AL • AFT	AL • AFT	AL • AFT	AL • AFT	AL • AFT	AL • AFT	AL • AFT	AL • AFT
Extra Tall 20in	Active	-	-	-	-	-	-	AL • AFT	AL • AFT	AL • AFT	AL • AFT	AL • AFT	AL • AFT	AL • AFT	AL • AFT	
	Deep	-	-	-	-	-	-	AL • AFT	AL • AFT	AL • AFT	AL • AFT	AL • AFT	AL • AFT	AL • AFT	AL • AFT	
	Extra Deep	-	-	-	-	-	-	AL • AFT	AL • AFT	AL • AFT	AL • AFT	AL • AFT	AL • AFT	AL • AFT	AL • AFT	

us Funding Pediatric

● Backs ≤ 14in wide are classified as pediatric | AFT available in 3 heights (Standard, Tall, Extra Tall)