



# CF BR 12 L A R

Shell type

Back rest

Width

Height

Contour

Mounting

## CF - CF Series

### Note:

Shell Assembly Weight

CF Series - .820lb / 370gm

- above weights are of a 16" STD

- all weights approx.

12" / 30.5cm

13" / 33cm

14" / 35.5cm

15" / 38cm

16" / 40.5cm

17" / 43cm

18" / 46cm

19" / 48cm

**L**  
low  
(10"/25cm)

**M**  
medium  
(13"/33cm)

**S**  
standard  
(16"/40.5cm)

**A**  
active  
(2.5"/63.5mm)  
**D**  
deep  
(4"/100mm)

**F**  
Fixed Elite  
(.60lb / 272gm)

**2**  
2-Point Pro  
(1.25lb/566gm)

**4**  
4-Point Pro  
(1.25lb/566gm)

**R**  
Quick  
Release  
(1.30lb / 589gm)

**W**  
Width  
Adjustable Elite  
(.63lb / 285gm)

**B**  
Basic Release  
(1.50lb / 680gm)



# CF Series

## Back Supports

The CF Series Back Support utilizes a high quality carbon fiber shell that is not only ultra lightweight, but also highly durable and attractive.

### Fits Chair Widths

12" (30.5cm) - 19" (48cm)

### Support Heights

Low: 10" (25cm)

Medium: 13" (33cm)

Standard: 16" (41cm)

### Contour Depths

Active: 2.5" (63.5mm)

Deep 4" (100mm)

### Weight

Weight: .820lb (370 Grams)

(approx. weight of a 16" / 40.5cm standard height back)

### NOTE:

Custom heights available



# CF Series

## Low Height (10"/25cm), Positioning Back Supports

### Active Contour Depth

#### Performance Mounting Hardware Options

Width of Wheelchair	Fixed ELITE	Width Adjustable ELITE	Quick Release	2-Point PRO	4-Point PRO
12" / 30.5cm	CFBR12LAF	CFBR12LAW	CFBR12LAR	CFBR12LA2	CFBR12LA4
13" / 33cm	CFBR13LAF	CFBR13LAW	CFBR13LAR	CFBR13LA2	CFBR13LA4
14" / 35.5cm	CFBR14LAF	CFBR14LAW	CFBR14LAR	CFBR14LA2	CFBR14LA4
15" / 38cm	CFBR15LAF	CFBR15LAW	CFBR15LAR	CFBR15LA2	CFBR15LA4
16" / 40.6cm	CFBR16LAF	CFBR16LAW	CFBR16LAR	CFBR16LA2	CFBR16LA4
17" / 43cm	CFBR17LAF	CFBR17LAW	CFBR17LAR	CFBR17LA2	CFBR17LA4
18" / 45.7cm	CFBR18LAF	CFBR18LAW	CFBR18LAR	CFBR18LA2	CFBR18LA4
19" / 48.2cm	n/a	CFBR19LAW	CFBR19LAR	CFBR19LA2	CFBR19LA4

### Deep Contour Depth, (4"/100mm)

Width of Wheelchair	Fixed ELITE	Width Adjustable ELITE	Quick Release	2-Point PRO	4-Point PRO
12" / 30.5cm	CFBR12LDF	CFBR12LDW	CFBR12LDR	CFBR12LD2	CFBR12LD4
13" / 33cm	CFBR13LDF	CFBR13LDW	CFBR13LDR	CFBR13LD2	CFBR13LD4
14" / 35.5cm	CFBR14LDF	CFBR14LDW	CFBR14LDR	CFBR14LD2	CFBR14LD4
15" / 38cm	CFBR15LDF	CFBR15LDW	CFBR15LDR	CFBR15LD2	CFBR15LD4
16" / 40.6cm	CFBR16LDF	CFBR16LDW	CFBR16LDR	CFBR16LD2	CFBR16LD4
17" / 43cm	CFBR17LDF	CFBR17LDW	CFBR17LDR	CFBR17LD2	CFBR17LD4
18" / 45.7cm	CFBR18LDF	CFBR18LDW	CFBR18LDR	CFBR18LD2	CFBR18LD4
19" / 48.2cm	n/a	CFBR19LDW	CFBR19LDR	CFBR19LD2	CFBR19LD4



#### NOTE:

- 1) Each back support includes: support shell, mounting hardware, contoured foam cushion, cover, and operation manual.
- 2) 19" CF Series backs utilize the same width shell supports as 18" CF Series backs

# CF Series

## Medium Height(13" / 33cm), Positioning Back Support Active Contour Depth(2.5"/63.5mm)

### Performance Mounting Hardware Options

Width of Wheelchair	Fixed ELITE	Width Adjustable ELITE	Quick Release	2-point PRO	4-Point PRO
12" / 30.5cm	CFBR12MAF	CFBR12MAW	CFBR12MAR	CFBR12MA2	CFBR12MA4
13" / 33cm	CFBR13MAF	CFBR13MAW	CFBR13MAR	CFBR13MA2	CFBR13MA4
14" / 35.5cm	CFBR14MAF	CFBR14MAW	CFBR14MAR	CFBR14MA2	CFBR14MA4
15" / 38cm	CFBR15MAF	CFBR15MAW	CFBR15MAR	CFBR15MA2	CFBR15MA4
16" / 40.6cm	CFBR16MAF	CFBR16MAW	CFBR16MAR	CFBR16MA2	CFBR16MA4
17" / 43cm	CFBR17MAF	CFBR17MAW	CFBR17MAR	CFBR17MA2	CFBR17MA4
18" / 45.7cm	CFBR18MAF	CFBR18MAW	CFBR18MAR	CFBR18MA2	CFBR18MA4
19" / 48.2cm	n/a	CFBR19MAW	CFBR19MAR	CFBR19MA2	CFBR19MA4

## Deep Contour Depth(4"/100mm)

### Performance Mounting Hardware Options

Width of Wheelchair	Fixed ELITE	Width Adjustable ELITE	Quick Release	2-point PRO	4-Point PRO
12" / 30.5cm	CFBR12MDF	CFBR12MDW	CFBR12MDR	CFBR12MD2	CFBR12MD4
13" / 33cm	CFBR13MDF	CFBR13MDW	CFBR13MDR	CFBR13MD2	CFBR13MD4
14" / 35.5cm	CFBR14MDF	CFBR14MDW	CFBR14MDR	CFBR14MD2	CFBR14MD4
15" / 38cm	CFBR15MDF	CFBR15MDW	CFBR15MDR	CFBR15MD2	CFBR15MD4
16" / 40.6cm	CFBR16MDF	CFBR16MDW	CFBR16MDR	CFBR16MD2	CFBR16MD4
17" / 43cm	CFBR17MDF	CFBR17MDW	CFBR17MDR	CFBR17MD2	CFBR17MD4
18" / 45.7cm	CFBR18MDF	CFBR18MDW	CFBR18MDR	CFBR18MD2	CFBR18MD4
19" / 48.2cm	n/a	CFBR19MDW	CFBR19MDR	CFBR19MD2	CFBR19MD4



#### NOTE:

1) Each back support includes: support shell, mounting hardware, contoured foam cushion, cover, and operation manual.

2) 19" CF Series backs utilize the same width shell supports as 18" CF Series backs

# CF Series

## Standard Height(16" / 41cm), Positioning Back Support

### Active Contour Depth(2.5"/63.5mm)

#### Mounting Hardware Options

Width of Wheelchair	Fixed ELITE	Width Adjustable ELITE	Quick Release	2-point PRO	4-Point PRO
12" / 30.5cm	CFBR12AF	CFBR12SAW	CFBR12SAR	CFBR12SA2	CFBR12SA4
13" / 33cm	CFBR13AF	CFBR13SAW	CFBR13SAR	CFBR13SA2	CFBR13SA4
14" / 35.5cm	CFBR14AF	CFBR14SAW	CFBR14SAR	CFBR14SA2	CFBR14SA4
15" / 38cm	CFBR15AF	CFBR15SAW	CFBR15SAR	CFBR15SA2	CFBR15SA4
16" / 40.6cm	CFBR16AF	CFBR16SAW	CFBR16SAR	CFBR16SA2	CFBR16SA4
17" / 43cm	CFBR17AF	CFBR17SAW	CFBR17SAR	CFBR17SA2	CFBR17SA4
18" / 45.7cm	CFBR18AF	CFBR18SAW	CFBR18SAR	CFBR18SA2	CFBR18SA4
19" / 48.2cm	n/a	CFBR19SAW	CFBR19SAR	CFBR19SA2	CFBR19SA4

### Deep Contour Depth(4"/100mm)

#### Performance Mounting Hardware Options

Width of Wheelchair	Fixed ELITE	Width Adjustable ELITE	Quick Release	2-point PRO	4-Point PRO
12" / 30.5cm	CFBR12SDF	CFBR12SDW	CFBR12SDR	CFBR12SD2	CFBR12SD4
13" / 33cm	CFBR13SDF	CFBR13SDW	CFBR13SDR	CFBR13SD2	CFBR13SD4
14" / 35.5cm	CFBR14SDF	CFBR14SDW	CFBR14SDR	CFBR14SD2	CFBR14SD4
15" / 38cm	CFBR15SDF	CFBR15SDW	CFBR15SDR	CFBR15SD2	CFBR15SD4
16" / 40.6cm	CFBR16SDF	CFBR16SDW	CFBR16SDR	CFBR16SD2	CFBR16SD4
17" / 43cm	CFBR17SDF	CFBR17SDW	CFBR17SDR	CFBR17SD2	CFBR17SD4
18" / 45.7cm	CFBR18SDF	CFBR18SDW	CFBR18SDR	CFBR18SD2	CFBR18SD4
19" / 48.2cm	n/a	CFBR19SDW	CFBR19SDR	CFBR19SD2	CFBR19SD4



#### NOTE:

- 1) Each back support includes: support shell, mounting hardware, contoured foam cushion, cover, and operation manual.
- 2) 19" CF Series backs utilize the same width shell supports as 18" CF Series backs



# Lateral Supports

Available on Active Contour AL Series and CF Series Backs\*

Integrates with mounting hardware for less weight

Custom offsets available upon request

## Mounting Options

Stationary

Swing away

\*Lateral Supports are not compatible with Fixed ELITE hardware



## Adjustable Lateral Bracket



## Pad Sizes



Small 3" H x 3" D



Medium 3" H x 4" D



Large 4.5" H x 4.5" D



X Large 4" H x 6" D



# PSIS Lumbar Support Insert

Posterior, Superior, Iliac Spine (PSIS) When positioning an ADI back, the interior portion of the back support should be positioned to make contact at the PSIS to provide optimal pelvic support and trunk control. If additional support is required, ADI offers the support insert shown below as an optional feature.

**To install, unzip the fabric cover and insert the lumbar support between the yellow foam surface and the fabric cover. Position vertically as desired.**

**Note:**

For a secure hold, double sided Velcro is provided.



Lumbar Support Insert



# Boards & Gloves



## Anti-slip Transfer Boards

Sliding boards are used to promote independence and safety during transfers.

For those with impaired balance and upper body strength, we highly recommend ADI's **Anti-Slip Transfer Boards** as a transfer-assist tool, often providing the needed assistance to enable therapist/caregiver performance of a safer and biomechanical friendly transfer.

### Anti-Slip Material

The bottom side of each transfer board has an anti-slip surface to keep the board in place with little to no movement during each transfer.

### Hand Hole

Three options of board hand hole cutouts for optimal user manage-ability. Your choice of either ver-tical, horizontal or no hand holes.

### Hardwood Materials

Durable, layered hardwood with multiple coats of clear sealer make each board smooth and aesthetically pleasing.

### Tapered Ends

Allows the user to easily slide the board underneath his or her leg/bottom for proper placement during transfer and to aid in reduction of skin sheering



21"(53.cm) 24"(61cm) 29"(73.6cm)

### Length Options



Vertical

Horizontal

None

### Hand Hole

Note: maximum weight capacity is 200lbs. (91kg)

## Push & Transfer Gloves

ADI's **Push & Transfer Gloves** have an open finger, half thumb, wrap-around design for easy don and doff. Made of leather with padded palms and wrists, the gloves save wear and tear on users' hands throughout a long day of active chair use. Three inches of Velcro closure provide a secure fit as well as allows the gloves to fit most hand sizes.



## Grip Options

All Leather / Rubber Grip

Open finger, half thumb design enables independent donning and doffing by persons with limited hand functions

Made of 100% leather

Non-slip rubber or leather surface enhances grip for transfer, push, and braking

Note: One size fits most

Padded palms and wrists

3" of dual Velcro closure for a secure fit



# Gloves

## Push & Transfer Gloves

**NOTE:** Color: Black Size: One size fits most.

Item #:	Description
ADG	Wheelchair Push & Transfer Glove w/ rubber grip palm (pair) back, all leather, open finger, 2" wrist protection
ADGL	Wheelchair Push & Transfer Glove w/ leather grip palm (pair) black, all leather, open finger, 2" wrist protection





# Transfer Boards

## Anti-Slip Transfer Boards


Item #:	Description
ADV21	21" Anti-Slip transfer board w/ two vertical hand holes
ADV24	24" Anti-Slip transfer board w/ two vertical hand holes
ADV29	29" Anti-Slip transfer board w/ two vertical hand holes
ADH21	21" Anti-Slip transfer board w/ two horizontal hand holes
ADH24	24" Anti-Slip transfer board w/ two horizontal hand holes
ADH29	29" Anti-Slip transfer board w/ two horizontal hand holes
NH21	21" Anti-Slip transfer board w/ no hand holes
NH24	24" Anti-Slip transfer board w/ no hand holes
NH29	29" Anti-Slip transfer board w/ no hand holes

### Three Options In Length




### Three Options in Hand Holes

- A. Vertical
- B. Horizontal
- C. No Hand Holes



### Anatomy of an ADI Transfer Board





# Mounting Hardware

ADI offers multiple options in back support mounting hardware, each designed to fit specific needs.

Whether you want the capability to remove the backrest from the chair, need the flexibility of adjustment, or just want to "set it and forget it", we've got your mounting hardware!

## Fixed Elite™



Stationary mounting to back post, and fixed-in-place positioning (no moving parts).

### Features:

- Exclusive for Carbon Fiber (CF) Series backs
- Compatible with Active and Deep contours
- Fits on backs with 14in (35.56cm) through 18in (45.72cm) widths
- Fits on backs with 10in (25.4cm), 13in (33.02cm) and 16in (40.64cm) heights

## Width Adjustable Elite™



Stationary mounting to back post, and width adjustable fixed positioning hardware.

### Features:

- Compatible with Aluminum and Carbon Fiber backs
- Fits All Contours, All Widths, All Height

## 2-Point Pro™



Fixed mounting width adjustable positioning (angle, depth and width). Two point mounting with four point stability.

### Features:

- Performance Stationary Mounting
- Adjustable Positioning
- Angle, Width and Depth
- Compatible with all ADI Backs and Contour Depths
- 10in (25.40cm) through 23in (58.42cm) Widths and All Heights

# Mounting Hardware

## 4-Point Pro™



Fixed mounting width adjustable positioning (angle, depth and width).

Basic Mounting offers removable mounting with angle, depth and width adjustment at an economical price.

### Features:

- Performance Stationary Mounting
- Adjustable Positioning
- Angle, Width and Depth
- Compatible with All ADI Backs, Contour Depths, Widths and Heights

## Quick Release



Removable mounting width adjustable positioning (angle, depth and width).

### Features:

- Performance Mounting
- Angle: 20° of adjustment without the loss of seat depth
- Depth: Range of adjustment is 1 1/4in / 3.18cm
- Width: Approximately 1in / 2.5cm

## Quick Release PRO



Removable mounting width adjustable positioning (angle, depth and width).

### Features:

- Performance Mounting
- Angle: 20° of adjustment without the loss of seat depth
- Depth: Range of adjustment is 1 1/4in / 3.18cm
- Width: Approximately 1in / 2.5cm



®© 1999-2017  
1.800.965.9229 | fax 1.800.806.1225  
info@stealthproducts.com  
www.stealthproducts.com  
104 John Kelly Drive, Burnet ,TX 78611

Distributed by:  
Pride Mobility Products Canada  
Beamsville Ontario  
(888) 570-1113 - pridemobility.ca