



AL (contoured top) and AFT (flat top) options; the AFT allowing for easier mounting of over the shoulder supports or chest harnesses.

The Q-Back 2.0 offers unmatched comfort. The 4-way stretch outer cover uses patented Coolcore® material for optimal breathability and moisture-wicking keeping users cooler in the summer and warmer in the winter while the 1 ¼" Medium Soft Sunmate® FRG insert, formulated with 33-36% natural plant-based ingredients and patented All Natural Fire-Retardant materials, offers a stable and consistent comfort.

Q-Back 2.0 back supports seamlessly integrates with the Quantum TRU-Balance® seating systems; Installation can be done in just a few minutes and uses less than 1/2" of seat depth, and is available in a wide range of widths, heights, and depths to fit most user profiles. Optional adjustable contour extensions are available for those that need customized support.

AL



AFT



POWERED BY **COOLCORE**



Adjustable contour extensions

AVAILABLE CONTOURS	Fixed contours: Active - 2.5" deep Deep - 4" deep Extra Deep - 6" deep Adjustable contour extensions: (for use with Active contour back) Deep - 4" deep Extra Deep - 6" deep
AVAILABLE WIDTHS	10" to 22" in 1" increments and 24"
AVAILABLE HEIGHTS	13", 16", 18", 20", 22"
WEIGHT CAPACITY	Up to 450 lb / 204 kg

