

QUANTUM® **EDGE 4**

THE LEADING EDGE IN
MID-WHEEL DRIVE JUST
TOOK A QUANTUM LEAP.



FDA Class II
Medical Device*

FEATURING 9 NEW
COLOR OPTIONS!



“ THE EDGE 4 ISN'T JUST MOBILITY,
IT'S MOMENTUM. IT CARRIES YOU
FORWARD WITH POWER, GRACE, AND THE
KIND OF CONFIDENCE THAT MAKES
THE WORLD FEEL WIDE OPEN. ”

STEWART LUNDY



QUANTUM®



quantumrehab.ca

Specifications

DRIVE WHEELS	14" pneumatic or flat free tires Optional: 14" knobby flat-free
CASTER WHEELS	7" flat free
MAXIMUM SPEED ¹	Up to 11.26 km/h (7 mph) (7.25 km/h (4.5 mph) with iLevel®)
GROUND CLEARANCE	2.87"
TURNING RADIUS	22.8"
LENGTH	38.27"
BASE WIDTH	24.5"

SEATING SIZES ⁵

TRU Balance® 4 Anterior Tilt	Width: 16-22" Depth: 14-22"
TRU Balance® 3 Power Tilt TRU Balance® 4 Power Tilt	Width: 12-22" Depth: 12-22"
TRU Balance® 3 Tilt and iLevel®	Width: 12-22" Depth: 12-22"
TRU Balance® Flex Seating	Width: 10-22" Depth: 10-20"

SEAT-TO-FLOOR HEIGHTS ⁵

TRU Balance® 4 Anterior Tilt	18.25-21.25"
TRU Balance® 3 Power Tilt TRU Balance® 4 Power Tilt	17.75-19.75"
TRU Balance® 3 Tilt & iLevel®	17.5-19.5"
TRU Balance® Flex Seating	16.25-18.25"

BATTERY SIZE ³	C5 = 65 Ah C20 = 75 Ah
BATTERY WEIGHT ⁷	47.8 lb (21.68 kg) each
AVAILABLE ELECTRONICS ⁶	Q-Logic NES (standard) Q-Logic 3 (optional) Q-Logic 3e (optional)
BATTERY CHARGER	8A Off-board
MOTOR PACKAGES	4-Pole
WEIGHT CAPACITY ⁴	300 lb (136 kg)
BASE WEIGHT	154 lb (70 kg) without batteries

- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
 2) Without front riggings
 3) Battery capacity may change when adding power positioning systems.
 4) Weight capacity is dictated by the selection of seating.
 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details.
 Listed measurements are in the lowest position; does not include elevating system unless specified.
 6) All Quantum Rehab® units can be ordered with specialty controls.
 7) Battery weight may vary +/- 3 lbs. based on manufacturer.

Features


- 4-pole motors with up to 11.26 km/h (7 mph) driving speeds standard
- Independent SRS (Smooth Ride Suspension)
- iLevel® option available with up to 12" of seat elevation at walking speed up to 7.25 km/h (4.5 mph)
- Standard accent lights and USB charger included
- Wide array of color choices
- 7" Caster and 14" drive wheels
- Group 24 batteries
- Easy battery access



Options

- iLevel® option allows 12" of seat elevation while at walking speed 7.25 km/h (4.5 mph)
- Q-Logic 3 Drive Controls
- TRU Balance® 3 & 4 Seating
- Depth-adjustable vent tray
- Alternative drive controls
- WC19 crash-test approved occupied transit kit (pending test results)
- Accu-Trac Advanced Tracking Technology (coming soon)

Colors



QUANTUM®

quantumrehab.ca • Follow us on social media for our latest innovations:

