

# EDGE 3

4.5 MPH iLEVEL® READY

**NOW AVAILABLE**

Max. speed up to 7 mph (11.3 km/h)  
with new optional 4-pole motors



Upgraded SRS (Smooth Ride Suspension) provides a smooth, comfortable ride and advanced stability. Optional iLevel® is available at 4.5 mph (7.2 km/h) with up to 12" of power adjustable seat height. Standard features include front and rear LED fender lights, a USB charger, and drive wheel color accents. The Edge 3 is available in a total of 19 color options, including seven bold colors.



RASPBERRY BERET

**Q**  
**QUANTUM®**  
#1 FOR REHAB POWER



quantumrehab.ca

# Specifications

<b>DRIVE WHEELS</b>	Standard: 14" Solid Optional: 14" Pneumatic; 14" Knobby
<b>CASTER WHEELS</b>	Front: 6" Solid Rear: 6" Solid
<b>ANTI-TIP WHEELS</b>	N/A
<b>MAXIMUM SPEED</b> <sup>1</sup>	Standard: Up to 6.25 mph (10 km/h) <sup>11</sup> Optional: Up to 7 mph (11.3 km/h) <sup>12</sup> Up to 4.5 mph (7.2 km/h) w/iLevel® activated
<b>GROUND CLEARANCE</b>	2.7" (battery tray)
<b>TURNING RADIUS</b> <sup>2</sup>	20.5"
<b>OVERALL LENGTH</b> <sup>2</sup>	32.23" without front riggings
<b>BASE WIDTH</b>	24.4"
<b>SEATING SIZES</b>	
<b>Synergy®/TRU-Balance® 3 Static Seat</b>	Width: 10-22" Depth: 10-22"
<b>TRU-Balance® 3 Power Tilt</b>	Width: 12-22" Depth: 12-22"
<b>TRU-Balance® 3 Tilt and iLevel®</b>	Width: 12-22" Depth: 12-22"
<b>TRU-Balance® Flex Seating</b>	Width: 15-22" Depth: 15-20"
<b>SEAT-TO-FLOOR HEIGHTS</b> <sup>10</sup>	
<b>Synergy®/TRU-Balance® 3 Static Seat</b> <sup>5</sup>	Static: 16.5 - 18.5" Synergy: 15.875 - 20.875"
<b>TRU-Balance® 3 Power Tilt</b> <sup>5</sup>	17.5-19.5"
<b>TRU-Balance® 3 Tilt &amp; iLevel®</b> <sup>5</sup>	17.5-19.5"
<b>TRU-Balance® Flex Seating</b>	16.5 - 18.25"
<b>MANUAL RECLINE</b> <sup>6</sup>	Yes
<b>BATTERY SIZE</b> <sup>3</sup>	NF-22 (2)
<b>BATTERY WEIGHT</b> <sup>7</sup>	38 lb each
<b>AVAILABLE ELECTRONICS</b> <sup>8</sup>	Q-Logic 3 EX Q-Logic 3e QLNES
<b>BATTERY CHARGER</b>	8A Off-board
<b>MOTOR PACKAGES</b>	2-pole, Mid-Wheel 6® 4-pole, Mid-Wheel 6®
<b>WEIGHT CAPACITY</b> <sup>4</sup>	300 lb
<b>BASE WEIGHT</b>	138.6 lb

1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.  
 2) Without front riggings  
 3) Battery capacity may change when adding power positioning systems.  
 4) Weight capacity is dictated by the selection of seating.  
 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.  
 6) Maximum dimensions of the Manual Recline vary by power base and seating system. Call for details.  
 7) Battery weight may vary +/- 3 lbs. based on manufacturer.  
 8) All Quantum Rehab® units can be ordered with specialty controls.  
 9) Call for availability  
 10) 16.5" seat-to-floor height is not available with power tilt or iLevel®. 19.5" seat-to-floor height is not available with static seat or power recline.  
 11) 6.25 mph (10 km/h) speed available only with Q-Logic 3.  
 12) Up to 7 mph (11.3 km/h) speed available only with 4-pole motor option and Q-Logic 3 or Q-Logic 3e.


# Features

- Upgraded SRS (Smooth Ride Suspension)
- Fender lights and USB charger included
- Bold, contemporary colors
- Proudly built in Canada

# Options

- TB4 seating system with up to 30 degrees of anterior tilt
- iLevel® with speeds up to 4.5 mph (7.2 km/h) while fully elevated
- **Up to 7 mph (11.3 km/h)** with new optional 4-pole motors<sup>12</sup>

# Colours



# Hydrographics

