

# Q6 EDGE 3 *Stretto*<sup>™</sup>

*Every other chair is  
at least 2.5" wider!*



The Q6 Edge 3 Stretto<sup>™</sup> features a narrow width for outstanding maneuverability and is equipped with independent SRS (Smooth Ride Suspension), making it ideal for children, teenagers and small adults. The Stretto is available with 12.5" or 14" drive wheels and optional 12" of iLevel<sup>®</sup> power adjustable seat height.



Narrowest  
Base in  
North America

  
**QUANTUM<sup>®</sup>**  
#1 FOR REHAB POWER



quantumrehab.ca

# Specifications

<b>DRIVE WHEELS</b>	12.5", 14" or 14" knobby tire pneumatic or flat free
<b>CASTER WHEELS</b>	<b>Front:</b> 5" flat free <b>Rear:</b> 6" flat free
<b>MAXIMUM SPEED</b> <sup>1</sup>	Up to 6.25 mph Up to 3.5 mph w/iLevel activated
<b>GROUND CLEARANCE</b>	2.9"
<b>TURNING RADIUS</b> <sup>2</sup>	21"
<b>OVERALL LENGTH</b> <sup>2</sup>	32.5" without front riggings
<b>BASE WIDTH</b>	21.75" with 14" drive wheels 20.75" with 12.5" drive wheels
<b>SEATING SIZES</b>	
<b>TRU-Balance® 3 Static Seat</b>	<b>Width:</b> 12-20" <b>Depth:</b> 12-20"
<b>TRU-Balance® 3 Power Tilt</b>	<b>Width:</b> 12-20" <b>Depth:</b> 12-20"
<b>TRU-Balance® 3 Tilt and iLevel®</b>	<b>Width:</b> 12-20" <b>Depth:</b> 12-20"
<b>TRU-Balance® Flex Seating</b>	<b>Width:</b> 15-20" <b>Depth:</b> 15-20"
<b>SEAT-TO-FLOOR HEIGHTS</b>	
<b>TRU-Balance® 3 Static Seat</b> <sup>5</sup>	16 -18"
<b>TRU-Balance® 3 Power Tilt</b> <sup>5</sup>	17.5-19.5"
<b>TRU-Balance® 3 Tilt &amp; iLevel®</b> <sup>5</sup>	17.5-19.5"
<b>TRU-Balance® Flex Seating</b>	16.5-18.25"
<b>MANUAL RECLINE</b> <sup>6</sup>	Yes (TRU-Balance® Flex only)
<b>BATTERY SIZE</b> <sup>3</sup>	40 AH (2) or 55 AH (2)
<b>BATTERY WEIGHT</b> <sup>7</sup>	40 AH: 31 lb each 55 AH: 40 lb each
<b>AVAILABLE ELECTRONICS</b> <sup>8</sup>	Q-Logic 3 EX Q-Logic 3e QLNES
<b>BATTERY CHARGER</b>	8A Off-board
<b>MOTOR PACKAGES</b>	2-Pole, Mid-Wheel 6 <sup>®</sup>
<b>WEIGHT CAPACITY</b> <sup>4</sup>	300 lb
<b>BASE WEIGHT</b>	126.4 lb
<b>MEDICARE CODE COVERAGE</b>	<b>Group 3:</b> Standard (K0848) Single Power (K0856) Multiple Power (K0861)

1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.  
 2) Without front riggings  
 3) Battery capacity may change when adding power positioning systems.  
 4) Weight capacity is dictated by the selection of seating.  
 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.  
 6) Maximum dimensions of the Manual Recline vary by power base and seating system. Call for details.  
 7) Battery weight may vary +/- 3 lbs. based on manufacturer.  
 8) All Quantum Rehab® units can be ordered with specialty controls.  
 9) Call for availability.


# Features

- Independent SRS (Smooth Ride Suspension)
- iLevel® option available with up to 12" of seat elevation at walking speed (up to 3.5 mph)
- Fender lights and USB charger included
- 17 bold, contemporary color choices



# Options

- TB4 Aterior tilt seating system
- iLevel® option allows 12" of seat elevation while at walking speed (3.5 mph)
- TRU-Balance® 3 Seating
- Q-Logic 3 EX Drive Controls
- Trilogy vent tray<sup>9</sup>
- Compact vent tray<sup>9</sup>
- 12" x 12" seat sizes with iLevel (SOF only)

# Colors



# Hydrographics



**QUANTUM**<sup>®</sup>  
 #1 FOR REHAB POWER

quantumrehab.ca |     