



OUTBACK

ADVANCED
OUTDOOR
REHAB POWER

BRYAN ANDERSON
Iraq War Veteran
and Quantum®
Spokesperson

QUANTUM® OUTBACK
shown with Q-Logic 3e
controller and optional
10" safe seat elevation

FDA Class II
Medical Device*



MAKE YOUR OWN PATH

Quantum Rehab® welcomes adventurers and nature enthusiasts to unleash their wild side with the release of our cutting-edge 4x4 all-terrain power chair, the Quantum® Outback! Designed for those who dare to defy conventional limits, this rugged mobility solution promises to elevate any outdoor adventure, whether navigating wooded trails or exploring uncharted terrain.



QUANTUM®
#1 FOR REHAB POWER



quantumrehab.ca

Specifications

DRIVE WHEELS	Standard: 15" pneumatic
MAXIMUM SPEED ¹	Up to 5 mph / 8 km/h (6 mph / 9.65 km/h available with optional upgrade)
GROUND CLEARANCE	4.92" (at motors)
TURNING RADIUS ²	49.2"
OVERALL LENGTH ²	42"
BASE WIDTH	28.8"
SUSPENSION	Independent SRS
RANGE ³	Up to 16 mi / 25.7 km at 300 lb / 136 kg Up to 17 mi / 27.3 km at 220 lb / 99.8 kg
SEATING SIZES	
Q-Captain Seat	Width: 16", 18", 20", 22" Depth: 16", 18", 20", 22"
TRU-Balance® 3 Power Tilt	Width: 16-22" Depth: 16-22"
TRU-Balance® 3 Power Lift & Tilt	Width: 16-22" Depth: 16-22"
SEAT-TO-FLOOR HEIGHTS	
Q-Captain Seat ⁵	19.5", 20.5", 21.5"
TRU-Balance® 3 Tilt ⁵	19", 20", 21"
TRU-Balance® 3 Lift & Tilt ⁵	19.5", 20.5", 21.5"
BATTERY SIZE ³	
	Group 27 (2)
BATTERY WEIGHT ⁶	
	58.6 lb / 26.6 kg (each)
AVAILABLE ELECTRONICS ⁷	
	120A Q-Logic 3
BATTERY CHARGER	
	8A Off-board
BRAKES	
	Regenerative and Electromechanical
DRIVETRAIN	
	4 motor, independent drive
MAXIMUM OBSTACLE CLIMBING	
	4"
WEIGHT CAPACITY ⁴	
	300 lb / 136 kg (TRU-Balance® 3) 350 lb / 158.8 kg (Q-Captain Seating)
BASE WEIGHT	
	168.6 lb / 76.5 kg

- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) Without front riggings
- 3) Battery capacity may change when adding power positioning systems.
- 4) Weight capacity is dictated by the selection of seating.
- 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and pneumatic tires. Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 6) Battery weight may vary +/- 3 lbs. based on manufacturer.
- 7) All Quantum Rehab® units can be ordered with specialty controls.
- 8) DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.
- 9) Extended lead times may apply

*The Quantum® Outback is an FDA Class II Medical Device designed to aid individuals with mobility impairments.

Standard Features

- 4-pole motors
- New camo color choices
- Full LED lighting package
- SRS (Smooth Ride Suspension)
- Unmatched driving performance and ride quality
- Advanced stability for smooth obstacle transitions

Options

- 10" of power adjustable safe seat elevation at speeds of up to 3.2 mph / 5.14 km/h
- Quantum® backup camera

Colors



FDA Class II
Medical Device*

QUANTUM®

#1 FOR REHAB POWER

(CAN) 888-570-1113 | quantumrehab.ca |